Masters Series™ Mechanical Heart Valve

Mechanical Heart Valve Sizer

REF

15 AHPJ-505, 15 MHPJ-505
905-15

EN: English	1
BG: Български	14
CS: Čeština	28
DA: Dansk	41
DE: Deutsch	54
EL: Ελληνικά	68
ES: Español	83
ET: Eesti keel	97
FI: Suomi	110
FR : Français	123
HR: Hrvatski	137
HU: Magyar	150
IT: Italiano	164
LT: Lietuvių	178
LV: Latviski	191
NL: Nederlands	204
NO: Norsk	218
PL: Polski	231
PT: Português	245
RO: Română	259
SK: Slovenčina	272
SL: Slovenščina	285
SV: Svenska	298
TR: Türkçe	311

Instructions for Use

This warranty is given by Abbott Medical. To make inquiries regarding this warranty, use the contact information for either the manufacturer or Australian Sponsor located on the back cover.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

The benefits to you given by this warranty are in addition to your other rights and remedies under the Australian Consumer Law.

Abbott Medical warrants to the buyer that the goods, for a period equal to the validated shelf life of the goods (the "Warranty Period"), shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Unless otherwise obligated by law or mandatory consumer guarantees, under this warranty Abbott Medical will, at its option, replace or repair the goods at its factory, if the buyer, at the buyer's expense, returns the goods within the Warranty Period to Abbott Medical to the address below and after Abbott Medical confirms that the goods are defective.

Abbott Medical
Roseville Distribution Center
Attn: Postmarket Surveillance
2305 Walnut Street
Roseville, MN 55113
Tel.: (651) 756-2000

To the maximum extent permitted by law, if a mandatory term is implied by law or a mandatory consumer guarantee applies to the goods, and the goods are not of a kind ordinarily acquired for personal, domestic or household use or consumption, Abbott Medical's liability for the breach of the term or guarantee is limited to, at Abbott Medical's option, either replacement or repair of the goods or payment of the costs of replacing or repairing the goods.

EXCEPT FOR MANDATORY TERMS IMPLIED BY LAW, MANDATORY CONSUMER GUARANTEES, OR AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED.

See the Terms and Conditions of Sale for further information.

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

Pat. <http://www.abbott.com/patents>

© 2024 Abbott. All Rights Reserved.

Device Description and Performance Characteristics

The Masters Series™ Mechanical Heart Valves are rotatable, bi-leaflet valves manufactured from pyrolytic carbon designed for implantation in the aortic or mitral position. Masters Series™ HP Mechanical Heart Valves are designed for supra-annular cuff placement. A sewing cuff of knitted double velour polyester is secured around the periphery of the valve to allow attachment of the valve to the native valve annulus. A metallic ring assembly within the sewing cuff allows the valve to be rotated *in-situ*.

Hemodynamic Plus (HP) models have been designed to maximize effective orifice area. This is made possible by removing the sewing cuff from the annulus.


Table 1. Reference Dimensions

Model	Tissue Annulus Diameter (mm)	Geometric Orifice Area (cm ²)	Implant Height Open (mm)	Overall Height Open (mm)	External Sewing Ring Diameter (mm)
15 MHPJ-505	15	1.3	2.2	8.5	20
15 AHPJ-505	15	1.3	4.6	8.5	20

Figure 1. Mitral valve measurements

Figure 2. Aortic valve measurements

Leaflets and Orifice Ring

Pyrolytic carbon coats the graphite substrates of the leaflets and orifice ring. Pyrolytic carbon was chosen for its biocompatibility and durability. The graphite substrate for the leaflets is tungsten-impregnated for radiopacity. X-ray or fluoroscopy may be used to visualize the valve. To obtain good x-ray visualization of the leaflets, orient the x-ray beam either parallel to the leaflet pivot mechanism axis, or perpendicular to the plane of the valve orifice. The orifice ring may be rotated to the surgeon's *in situ* preference after the valve is sutured in place.

Sewing Cuff

The double-velour knitted polyester fiber in the polyester cuff allows for rapid, controlled endothelial ingrowth over the entire sewing cuff. The cuff suture markers may be used to provide reference points for the valve's orientation or the placement of sutures. The aortic valve has three suture markers. The mitral valve has four suture markers.

Figure 3. Four suture markers on mitral sewing cuff


Figure 4. Three suture markers on aortic sewing cuff

Holder/Rotator

A holder/rotator is attached to each valve. Use this holder/rotator to rotate the valve *in situ*.

Intended Purpose

The Masters Series™ HP 15 mm Valve is intended to replace native aortic or mitral heart valves or previously implanted prosthetic valves.

The Valve Sizer is intended to aid in the selection of the appropriate mechanical heart valve size.

Indications

The Masters Series™ Mechanical Heart Valve is indicated for use as a replacement valve in patients with a diseased, damaged, or malfunctioning mitral or aortic heart valve or as a replacement for a mitral or aortic prosthetic heart valve.

The Sizer Set Model 905-15 is indicated to aid in the selection of the 15 AHPJ-505 and 15 MHPJ-505 valves.

Target Patient Population

Pediatric patients with a diseased, damaged, or malfunctioning aortic or mitral valve.

Intended User

Physicians trained in surgical mitral and/or aortic valve replacement in a pediatric population.

The device implantation should be performed in a standard operating room equipped for open-heart surgery.

Intended Clinical Benefits

Valve

Based on clinical evidence, the intended clinical benefits for the Masters Series™ Valves are:

- Low rate of reintervention
- Improved symptoms and functional capacity as assessed by New York Heart Association [NYHA] Classification
- Improved hemodynamic valve performance by decreasing valvular stenosis and/or valvular regurgitation
- Low rates of the following adverse events: valve related mortality, valve dysfunction, valve thrombosis, major paravalvular leak, hemolysis, major hemorrhage, thromboembolism, and endocarditis

Sizer

The Sizer Set Model 905-15 provides an indirect clinical benefit as an accessory to aid in the implantation of Masters Series™ 15 mm Mechanical Heart Valves.

Summary of Safety and Clinical Performance (SSCP)

A summary of the safety and clinical performance (SSCP) for this device is available at <https://ec.europa.eu/tools/eudamed>. Search for the device using the basic UDI-DI listed in the following table. This is the SSCP location after the launch of the European Database on Medical Devices/Eudamed.

Table 2. Masters Series™ Mechanical Heart Valves Basic UDI-DI

Basic UDI-DI	Basic UDI-DI Description
5415067SHV2001PL	Surgical Heart Valve Products, Mechanical Heart Valves

Limitations

No known limitations.

Contraindications

The Masters Series™ Mechanical Heart Valve is contraindicated for individuals unable to tolerate anticoagulation therapy.

Warnings

Valve

- For single use only. Attempts to reuse the valve may result in valve malfunction, inadequate sterilization, or patient harm.
- Use only mechanical heart valve sizers from Abbott Medical.
- Do not use if:
 - The valve has been dropped, damaged, or mishandled in any way.
 - The expiration date has elapsed.
 - The tamper-evident container seal or inner/outer tray seals are damaged, broken, or missing.
- Remove any residual tissue that may impair valve size selection, correct seating of the valve, rotation of the valve, or leaflet motion.
- Proper valve size selection is crucial. Do not oversize the valve. If the native annulus measurement falls between two Masters Series™ Mechanical Heart Valve sizes, use the smaller prosthetic valve size.
- The outer tray is not sterile, and should not be placed in the sterile field.
- To minimize direct handling of the valve during implantation, do not remove the holder/rotator until the valve has been seated in the annulus.
- Do not use hard or rigid instruments to test leaflet mobility, as this may result in structural damage to the valve or thromboembolic complications. Use a Model LT100 leaflet tester to gently test valve leaflet mobility.
- Place sutures in the outer half of the valve sewing cuff.
- Never apply force to the valve leaflets. Force may cause structural damage to the valve.
- Use only Valve Holder/Rotators from Abbott Medical to perform valve rotation. Use of other instruments could result in structural damage. The valve holder/rotator is intended for single use only and should be discarded after surgery.
- The two retention sutures on the valve holder/rotator must be cut and removed before the valve can be rotated.
- Do not pass catheters or other instruments through Masters Series™ mechanical heart valves. This could result in scratched or damaged valve components, leaflet fracture, or dislodgment.
- Cut suture ends short, especially in the vicinity of the pivot guards, to prevent leaflet impingement.

- Do not use cutting edge or cutting point needles in the sewing cuff. Use of standard round or taper point needles is recommended to prevent cutting the polyester fibers.

Precautions

Valve

- Do not touch the prosthetic valve unnecessarily, even with gloved hands. This may cause scratches or surface imperfections that may lead to thrombus formation.
- Be careful not to cut or tear the valve sewing cuff when removing the identification tag and the holder/rotator from the valve.
- Before placing sutures in the valve sewing cuff, verify that the valve is mounted correctly on the valve holder/rotator.
- To avoid structural damage, the mitral valve must be rotated in the fully open position.
- To avoid structural damage, the aortic valve must be rotated in the fully closed position.
- To minimize rotational torque, verify that the valve holder/rotator is properly seated in the valve, and that the valve holder handle is perpendicular to the valve.
- Remove any loose suture or thread, which may be a source of thrombus or thromboembolism.
- Implantation of a prosthetic valve too large for the annulus may result in increased risk of damage to the conductive system, obstruction of the left ventricular outflow tract, impairment of valve mobility, damage to the left circumflex artery, and damage to surrounding tissues or cardiac structures including obstruction and/or distortion of adjacent cardiac structures.

NOTE: Prospective data to support safety and effectiveness of the 15-mm hp valve implanted in the aortic position are not currently available.

- Safety and performance of the Masters Series™ Mechanical Heart Valve in the 15 mm size has not been established for the following specific populations:
 - Patients who are pregnant
 - Nursing mothers
 - Non-pediatric patients

Sizer

- Instruments must be cleaned and sterilized prior to use.
- Do not use cracked, deformed, discolored/rusted, or damaged instruments. If visual signs of deterioration are apparent, contact an Abbott customer care representative.
- Improper cleaning may result in an immunological or toxic reaction.
- Instrument sterilization temperature must not exceed 280°F (138°C).
- Do not bend flexible instrument handles beyond a 90° angle.
- Instruments must be sterilized in a tray or container that is permeable to steam.
- Do not expose instruments to cleaning or rinse agents that are not compatible with polysulfone or polyphenylsulfone.

Potential Adverse Events

Complications associated with replacement mechanical heart valves include, but are not limited to:

- adjacent cardiac structure interference
- death
- heart block requiring pacemaker implant
- heart failure
- hemolysis
- hemorrhagic complications secondary to anticoagulation therapy
- infections
- myocardial infarction
- prosthetic failure
- stroke
- thrombus or thromboembolism
- unacceptable hemodynamic performance
- valve dehiscence

Any of these complications may require reoperation or explantation of the device.

Device Lifetime

The product is designed and tested to support an expected lifetime of 15 years. Implant durations shorter than 15 years are expected in pediatric patients due to the need for reimplantation with a larger prosthesis to accommodate normal anatomical growth.

The device lifetime of the sizer and holder handles is for the duration of the implant procedure. The device can safely be re-used up to 100 cleaning cycles when cleaning using the specified reprocessing methods.

The reuse of the Sizers and Holder Handles is determined through visual inspections after each resterilization cycle and prior to use. Do not use cracked, deformed, discolored/rusted, or damaged instruments. If visual signs of deterioration are apparent, contact an Abbott customer care representative.

Magnetic Resonance (MR) Safety Information

A person with the Masters Series™ Mechanical Heart Valves and Valved Grafts may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.

Device Names	Masters Series™ Mechanical Heart Valves Masters Series™ Aortic Valved Graft Masters Series™ HP Valved Graft
Static Magnetic Field Strength (B0)	1.5T or 3.0T
Maximum Spatial Field Gradient	≤20 T/m (2000 gauss/cm)

RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Body coil
Operating Mode	Normal Operating Mode
Maximum Whole Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	N/A
Scan Duration	2 W/kg whole body average SAR for 15 minutes of continuous scanning
MR Image Artifact	The presence of this implant may produce an image artifact.

Packaging

The valve has been steam sterilized. The valve remains sterile until the “USE BY” date on the packaging, if package integrity is maintained.

The package includes:

- Sealed non-sterile outer tray
- Sealed sterile inner tray
- One (1) valve with identification tag, mounted on a plastic holder/rotator
- One (1) disposable support collar
- One (1) insert with Instructions for Use web address
- One (1) Medical Device Registration Form with attached Patient Identification Card and return envelope

The sizer is supplied non-sterile. Clean and sterilize the sizer before each use following the reprocessing guidelines described in these instructions.

Storage

Valve

There are no special storage requirements.

Sizer

Handle with care. Keep dry. Keep away from sunlight. Proper care and handling of the sizer set, sizer and handle is required to ensure that each reusable accessory maintains its appearance, corrosion resistance, and proper function over time.

Caution: Do not use cracked, deformed, discolored/rusted, or damaged instruments. If visual signs of deterioration are apparent, contact an Abbott customer care representative.

Valve Resterilization

If resterilization of the valve is necessary, use only the recommended steam cycles and follow the instructions below. **The valve must not be resterilized more than one (1) time.**

1. Remove the inner tray containing the valve and also remove the inner tray cover.
2. Place the inner tray inside a steam-permeable sterilization bag or wrap.

Table 3. Recommended Sterilization Cycle Parameters

Vacuum Cycle	Pre-Vacuum Steam	Pre-Vacuum Flash
Purge Time:	6 minutes	6 minutes
Pulses:	2	2
Pulse Pressure:	204.7 kPa (absolute)	204.7 kPa (absolute)
Pulse Vacuum:	23.4 kPa (absolute)	23.4 kPa (absolute)
Sterilize Time:	28 minutes	7 minutes
Sterilize Temperature:	122°C	132°C
Post-Vacuum:	13.3 kPa (absolute)	13.3 kPa (absolute)
Vacuum Dry Time:	10 minutes	10 minutes

Accessories

Sizer

The model 905-15 sizer may be used to facilitate proper valve size selection. The sizer handle may be bent in any direction up to a 90° angle without degradation. When the flexible handle is exposed to the heat produced during autoclave sterilization, it will revert to its original shape. The sizer is supplied non-sterile.

Alternatively, a Hegar dilator may be used to facilitate proper valve size selection.

Holder Handles

The following holder handles may be purchased separately. Consult the appropriate instrument Instructions for Use for full product description, and relevant cleaning and sterilization information. Instruments have been designed and tested for repeated use, however, if visual signs of deterioration become apparent, do not use the instrument, and contact Customer Service for a replacement.

Table 4. Optional Holder Handles

Description	Model
Flexible holder handle	905-HH
Rigid holder handle	905-RHH
Rigid mitral holder handle	905-MHH

Leaflet Tester

Use Leaflet Tester model LT100 to test leaflet mobility.

Figure 5. LT100 Leaflet Tester

Directions for Use

Removing the Native Valve

Excise the native valve and prepare the annulus for the valve replacement.


WARNING: Remove any residual tissue that may impair valve size selection, correct seating of the valve, rotation of the valve, or leaflet motion.

Sizing

If available, use sizer model 905-15 to confirm fit of the 15 AHPJ-505 or 15 MHPJ-505 valve. The 905-15 sizer is a double-ended tool with a cylindrical annulus sizing end and a flanged HP sewing cuff replica end. Use the cylindrical end of the tool to determine the appropriate valve size. The cylindrical end of the sizer must pass readily without resistance through the annulus. Clean and sterilize the sizer before use.

Alternatively, a Hegar dilator may be used to size the valve.


Figure 6. Cylindrical end of the 905-15 sizer in the annulus

Use the flanged end of the sizer to visualize supra-annular placement of the HP sewing cuff.

CAUTION: Do not pass the flanged end of the sizer through the annulus.

Figure 7. Flanged valve cuff replica end of the 905-15 sizer above the annulus

WARNING: Proper valve size selection is crucial. Do not oversize the valve. If the native annulus measurement falls between two Masters Series™ Mechanical Heart Valve sizes, use the smaller prosthetic valve size.

CAUTION: Implantation of a prosthetic valve too large for the annulus may result in increased risk of damage to the conductive system, obstruction of the left ventricular outflow tract, impairment of valve mobility, damage to the left circumflex artery, and damage to surrounding tissues or cardiac structures including obstruction and/or distortion of adjacent cardiac structures.


Circulating Nurse (Non-Sterile)

1. Remove the non-sterile outer tray from the product box.

WARNING: Do not use the device if the expiration date has elapsed.

2. Verify that the catalog number and serial number on the outer tray are identical to those on the box label. If the information is not identical, do not use the device, and contact Customer Service as soon as possible.
3. Open the outer tray as shown below. Do not touch the sterile inner tray cover.

Figure 8. Remove the cover from the outer tray


4. Holding the outer tray by the bottom, present the inner tray to the sterile scrub nurse or surgeon. Do not touch the sterile inner tray.
5. Complete the patient [medical device] registration information, as described in the "Patient Registration" section.

Scrub Nurse (Sterile)

1. Remove the sterile inner tray from the outer tray. Avoid contact with the exterior of the outer tray, which is not sterile.

Figure 9. Remove the inner tray

WARNING: Examine the packaging carefully, and ensure the tray is unopened and unbroken. If you note any damage, do not use the device, and contact customer service.

2. Holding the inner tray with the cover upward, grasp the tab and pull back to completely remove the inner tray cover.

Figure 10. Remove the inner tray cover

3. Press the sterilized mechanical valve holder handle into the valve holder/rotator. The valve holder handle should be straight when inserted into the valve holder. Verify that the valve holder handle is securely attached to the valve.

Figure 11. Insert the holder handle in the valve holder

4. To remove the valve from the inner tray, firmly lift the valve holder handle and the valve support collar.


Figure 12. Remove the valve and collar from the inner tray

NOTE: Always use the valve holder handle to remove the valve from the inner tray.

5. Prior to implanting the valve, remove the support collar from the valve by placing two fingers under the collar, pinching with the thumb and pulling back gently.

Figure 13. Remove the collar from the valve prior to implantation


6. An identification tag is sutured to the valve sewing cuff. Verify that the valve size and model number on the identification tag are identical to the valve size and model number listed on the packaging. If the size and model number do not match, do not use the valve.
7. Remove the identification tag and its suture prior to implantation. Retain the identification tag for the patient record.

CAUTION: Be careful not to cut or tear the valve sewing cuff when removing the identification tag from the valve.

Valve Implantation


1. Place sutures in the annulus. Suturing techniques may vary according to the preference of the implanting physician and the needs of the patient. Experience indicates that there are several satisfactory suturing methods.
2. Use the valve holder handle to align the valve in the annulus. **Orient the valve so that blood flow is always into the pivot guards.** Align the valve such that the pivot guards are in the desired orientation in the annulus.

Figure 14. Blood flow into pivot guards, mitral valve inflow mount

1. Flow
2. Pivot guards
3. Leaflets

Figure 15. Blood flow into pivot guards, aortic valve outflow mount


1. Pivot guards
2. Flow

NOTE: Do not use instruments other than the valve holder handle to align and seat the valve, and exercise extreme care to avoid putting stress on the valve orifice or leaflets.

3. Place sutures in the outer half of the valve sewing cuff. Incorporate adequate sewing cuff material into each suture pass to maintain the valve in the desired position. It may be helpful to tag the sutures adjacent to pivot guards when inserting them. To prevent cutting the fibers in the sewing cuff, use of standard round or taper point needles is recommended.
4. Parachute the valve into place.
5. Cut the two retention sutures from the holder/rotator, and gently withdraw the holder/rotator from the valve. Remove the valve retention sutures. Maintain the holder/rotator in the sterile field for further use as a rotation instrument.

Figure 16. Cut the valve retention suture


6. Using the Model LT100 Leaflet Tester, open the valve and inspect the area for any obstructive tissue. If visualization is inadequate, use the leaflet tester to confirm free leaflet motion.
WARNING: Do not use hard or rigid instruments to test leaflet mobility, as this may result in structural damage to the valve or thromboembolic complications.
7. Tie down the pivot guard sutures first, and then tie down the remaining sutures.
WARNING: Cut suture ends short, especially in the vicinity of the pivot guards, to prevent leaflet impingement.
8. Test leaflet motion again, and if desired, rotate the valve using the holder/rotator handle (see "Valve Rotation").
WARNING: The two retention sutures on the valve holder/rotator must be cut and removed before the valve can be rotated.

Valve Rotation

Using the valve holder/rotator and a valve holder handle, rotate the valve *in situ* to the desired position. The valve should rotate freely. If resistance is noted, the valve holder/rotator may not be properly seated in the valve, or the valve may be oversized. If the valve does not rotate freely, do not force valve rotation.

Figure 17. Use the Masters Series™ valve holder/rotator and the valve holder handle to rotate the valve

WARNING: Use only valve holder/rotators from Abbott to perform valve rotation. Use of other instruments could result in structural damage. The valve holder/rotator is intended for single use only and should be discarded after surgery.

CAUTION: To avoid structural damage, the mitral valve must be rotated in the fully open position.

CAUTION: To avoid structural damage, the aortic valve must be rotated in the fully closed position.

CAUTION: To minimize rotational torque, verify that the valve holder/rotator is properly seated in the valve, and that the valve holder handle is perpendicular to the valve.

Figure 18. Improperly seated valve holder/rotator



Figure 19. Properly seated valve holder/rotator

Postoperative Considerations

Echocardiography is recommended to assess valvular competency and performance. Fluoroscopy is particularly useful for determining leaflet motion in bileaflet valves.

WARNING: Do not pass catheters or other instruments through Masters Series™ Mechanical Heart Valves. This could result in scratched or damaged valve components, or leaflet fracture or dislodgment.

Patients who undergo dental or other procedures that are potentially bacteremic should receive endocarditis prophylactic antibiotic.

Anticoagulation/Antiplatelet Therapy

As there are insufficient data to indicate otherwise, Abbott recommends that patients implanted with the Masters Series™ 15 mm Mechanical Heart Valve be routinely maintained on anticoagulants unless, for other reasons, it is not medically indicated.

The recommended anticoagulation therapy for pediatric patients implanted with the Masters Series™ 15 mm Mechanical Heart Valve are based on the 2013 AHA Guidelines for the Prevention and Treatment of Thrombosis in Pediatric Patients ¹ in combination with experience from the HALO IDE Study (ClinicalTrials.gov Identifier: NCT02097420).

Abbott Medical recommends that patients implanted with the Masters Series™ 15 mm Mechanical Heart Valve be routinely maintained on anticoagulants to avoid the risk of thrombus formation and thromboembolic complications. Acceptable forms of anticoagulants include intravenous unfractionated heparin or oral warfarin. When using warfarin the recommended target INR is 2.5 to 3.5 for valves implanted in the mitral position and 2.0 to 3.0 for valves implanted in the aortic position in the absence of risk factors for thrombus formation. In the presence of risk factors consider using a higher target INR or adding aspirin to therapeutic warfarin. Risk factors include first 3 months post-implant, younger age (<6 months), small size (<6 kg), low flow state, unreliable oral intake, previous thromboembolism, and hypercoagulable condition. Home INR monitoring may be associated with a reduced rate of complications. Use of low molecular weight heparin may be associated with an increased rate of complications.

In patients receiving warfarin, the INR should be monitored daily until therapeutic levels are achieved, and then may be decreased in frequency when stable with a minimum of monthly testing. The INR should be tested whenever there is illness or with any changes in medication or diet.

Patient Identification Card

A Patient Identification Card is present within the packaging of the valve. Hospital personnel should complete the front side of the Patient Identification Card by printing the patient's name, implant date, along with the physician's name and healthcare facility information and phone number in the designated areas. A peel-away device label is provided in the sterile packaging which includes device information (model number, serial number, UDI number, and device name). Apply the peel-away device label to the back of the card in the designated area. If a peel-away label is not available, write the information in the designated area on the back of the card. In procedures where multiple implantable devices are implanted, include information for all implants in the designated areas.

NOTE: Write the valve serial number in the SN/LOT field.

After the card is completed, separate the top portion from the Patient Identification Card by tearing along the perforated line. Provide the completed Patient Identification Card to the patient. Instruct the patient to carry the Patient Identification Card at all times and to present the implant card to hospital staff in emergency situations. To obtain a replacement card if a patient loses or damages their card, contact Abbott Medical Technical Support. Once the replacement card is received, handwrite the required information in the respective areas and provide to the patient.

Product Materials

The following product materials are intended to come into contact with tissue.

Table 5. Valve Components and Material Breakdown

Component	Material Description	Substance	Substance Concentration within each component (% w/w)(maximum)	Max Weight of component (≤) in grams
Leaflets and Orifice	Graphite Substrate with Pyrolytic Carbon Coating	Graphite Substrate with Pyrolytic Carbon Coating	100	0.6
Cuff	Polyester Fabric	Polyester Blend	100	
Sutures	PTFE-impregnated braided polyester	PTFE	3-7	0.3
		PET	93-97	
Ring Retainers	Cobalt Chromium Alloy	Nickel	35	0.2
		Cobalt	35	
		Chromium	20	
		Molybdenum	10	
Spring	Cobalt Chromium Alloy	Cobalt	40	0.2
		Chromium	20	
		Nickel	15	
		Molybdenum	7	
		Manganese	2	
		Silicon	1.2 MAX	
		Carbon	0.15 MAX	
		Beryllium	0.10 MAX	
		Phosphorous	0.015 MAX	
		Iron	Balance	

CAUTION: One or more components of the Masters Series™ Mechanical Heart Valve contains the following substance defined as CMR 1B in a concentration above 0.1% weight by weight:

Cobalt; Chemical Abstracts Service CAS No. 7440-48-4; EC No. 231-158-0

Current scientific evidence supports that medical devices manufactured from metal alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

Incident Reporting

If, in the course of use of this device, you have reason to believe that a serious incident occurred, please report it to the manufacturer. For customers in the European Union, report the serious incident to your national authority as well as to the manufacturer.

Disposal

This instructions for use is recyclable. Dispose of all packaging materials as appropriate. Dispose of valves and accessories per standard solid biohazard waste procedures.

¹ Giglia, Therese M., M. Patricia Massicotte, James S. Tweddell, Robyn J. Barst, Mary Bauman, Christopher C. Erickson, Timothy F. Feltes et al. "Prevention and Treatment of Thrombosis in Pediatric and Congenital Heart Disease." Circulation 128, no. 24 (2013): 2622-2703.

Intra-Operative Assessment

The suggested method for assessing competence of the valve is with intra-operative Doppler echocardiography.

Sizer Reprocessing Requirements

These reprocessing instructions were validated for sizer model 905-15. If alternate methods are used, it is the user's responsibility to validate and subsequently ensure that their methods provide effective cleaning and sterilization to a sterility assurance level of at least 10^{-6} .

This instrument has been designed and tested for repeated use. However, if visual signs of deterioration become apparent, do not use the instrument and contact customer service for a replacement.

Testing indicates that this instrument may undergo 100 cleaning cycles when cleaned using the specified reprocessing method. Deviations from the specified method may result in reduced instrument life or inadequate cleaning.

Reprocess instruments as soon as reasonably possible following use. Handle the instruments as standard solid biohazard until point of reprocessing. After each reprocessing cycle inspect instruments for signs of cracking, crazing, or degradation that may affect function.

CAUTION: Instruments must be cleaned and sterilized prior to each use.

CAUTION: Do not use cracked, deformed, discolored/rusted, or damaged instruments.

CAUTION: Improper cleaning may result in an immunological or toxic reaction.

CAUTION: Do not expose instruments to cleaning or rinse agents that are not compatible with polyphenylsulfone.

Autoclavable sizer set trays do not provide a sterile barrier. Handle sizer sets in the same manner as other reusable instruments that require packaging or wrapping during autoclave sterilization.

All instruments must be sterilized in a tray or container that is permeable to steam.

Reprocessing Option 1

Manual and Enzymatic Cleaning Method

The following manual and enzymatic cleaning method was validated:

1. Wipe the instrument with a wet, disposable sponge to remove gross soil.
2. Rinse the instrument in hot, running tap water for a minimum of 20 seconds to remove visible blood soil.
3. Soak the instrument in a bath of Terg-A-Zyme[®] Enzymatic Cleaner (Alconox) for five minutes.
4. Items should be mechanically cleaned with appropriately sized non-abrasive brushes until visibly clean.
5. Rinse instrument with hot de-ionized (DI) water for a minimum of 20 seconds to remove remaining soil and cleaning agent.
6. Allow to air dry or dry the instrument with a clean, disposable, absorbent, non-shedding wipe.

Automated Cleaning, Washer-Disinfector Method

The automated washer-disinfector method described below was validated. Testing was conducted using HAMO[®] LS1000 washer-disinfector with Renuzyme[®] enzyme detergent.

Prior to placing the instruments in the washer-disinfector, remove visible blood soil using warm DI water and a brush. Secure instruments in the washer-disinfector and perform the following programmed cleaning cycle:

Table 6. Automated Washer-Disinfector Cycle

Step	Cycle	Duration	Tap Water Temperature	Comments
1	Pre-Wash	10 minutes	50°C / 122°F ±5°C / 9°F	Do not use detergents or additives. Drain when complete.
2	Wash	10 minutes	60°C / 140°F ±5°C / 9°F	Follow instructions from manufacturer of washer-disinfector for optimal mix ratio and concentrations of enzyme detergent. Drain when complete.
3	Rinse	10 minutes	90°C / 194°F ±5°C / 9°F	Drain when complete.
4	Cool	5 minutes	---	Remove parts from washer-disinfector and cool at room temperature.
5	Dry	---	---	Allow to air dry or dry the instrument with a clean, disposable, absorbent, non-shedding wipe.

Sterilization

Following cleaning and drying, package the device in a Tyvek[®] pouch or double-wrap with sterilization wrap.

Sterilize instruments with steam before each use. Do not sterilize using any method other than steam.

CAUTION: Sterilization temperatures must not exceed 138°C/280°F.

NOTE: Higher temperatures or longer sterilization times may be used; however, increasing the sterilization time or temperature (not to exceed 138°C/280°F) of the stated cycles may reduce instrument life.

Table 7. Sterilization Parameters

Pre-Vac Sterilization		
	Low Temperature Cycle	High Temperature Cycle
Temperature	121°C (250°F)	132°C (270°F)
Sterilization Time	30 minutes	3 minutes
Gravity Displacement Sterilization		

Table 7. Sterilization Parameters

	Low Temperature Cycle	High Temperature Cycle
Temperature	121°C (250°F)	132°C (270°F)
Sterilization Time	30 minutes	6 minutes

Reprocessing Option 2

Automated Cleaning, Washer-Disinfector Method

The automated washer-disinfector method described below was validated using STERIS[®] Synergy[®] washer-disinfector, model E3023-2.

Prior to placing instruments in the washer-disinfector, remove visible blood soil under cold tap water with a soft bristle brush. Secure the instruments in the washer-disinfector and perform the following programmed cleaning cycle.

After each cleaning cycle, inspect instruments for signs of cracking, crazing, or degradation that may affect function.

Table 8. Automated Washer-Disinfector Cycle

Step	Cycle	Duration	Cycle Requirements	Comments
1	Pre-Wash	2 minutes	Cold tap water	Do not use detergents or additives. Drain when complete.
2	Alkaline Detergent Wash	5 minutes	55-60°C / 131-140°F Neodisher [®] MediClean forte [®] diluted to 0.5% (v/v)	Drain when complete.
3	Neutralizing Rinse	5 minutes	~25°C / 77°F tap water	Drain when complete.
4	Intermediate Rinse	4 minutes	Cold tap water	Drain when complete.
5	Disinfection Rinse	5 minutes	93°C / 199°F tap water	Drain when complete.
6	Dry	--	--	Allow to air dry

Sterilization

Following cleaning and drying, package the device in a Tyvek[®] pouch or double-wrap with sterilization wrap.

Sterilize instruments with steam before each use. Do not sterilize using any method other than steam.

CAUTION: Sterilization temperatures must not exceed 138°C/280°F.

NOTE: Higher temperatures or longer sterilization times may be used; however, increasing the sterilization time or temperature (not to exceed 138°C/280°F) of the stated cycles may reduce instrument life.

Table 9. Sterilization Parameters

Pre-Vac Sterilization

Temperature	134°C/273°F
Time	5 minutes

Limited Warranty

Abbott Medical warrants that reasonable care has been used in the manufacturing of this device. THIS WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES NOT EXPRESSLY SET FORTH HEREIN, WHETHER EXPRESS OR IMPLIED BY OPERATION OF LAW OR OTHERWISE INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, since handling, storage, cleaning, and sterilization of this device as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Abbott's control directly affect this device and the results obtained from its use. ABBOTT SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE directly or indirectly arising from the use of this device other than the replacement of all or part of it. Abbott neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.



Some states in the United States do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you. This limited warranty gives you specific legal rights, and you may have other rights which vary from jurisdiction to jurisdiction.

Descriptions of specifications, appearing in Abbott literature, are meant solely to generally describe the device at the time of manufacture and do not constitute any express warranties.

Symbols

The symbols below and harmonized symbols may be found on the product or product label. For harmonized symbols, refer to the Universal Symbols Glossary at medical.abbott/manuals.

Symbol	Description
Poly	Standard Polyester Cuff
POLY HP	Standard Polyester HP Cuff
PTFE	Standard PTFE Cuff
EXP POLY	Expanded Polyester Cuff
EXP POLY HP	Expanded Polyester HP Cuff
MITRAL	Mitral
AORTIC	Aortic

Symbol	Description
	Rotatable
 medical.abbott/manuals	Follow instructions for use on this website
	Date
	Healthcare center or physician
	Patient identification
	Physician telephone
 UA.TR.101	Ukrainian mark of conformity to the technical regulations, where 101 is the assessment body.
	Medical Device
	Contents sterile: steam
 2797	European conformity, affixed in accordance with the relevant provisions of European Council Regulation 2017/745 (NB 2797). Hereby, Abbott Medical declares that this device is in compliance with the relevant provisions of this regulation.
 medical.abbott/manuals	Patient Information Website
	United Kingdom Conformity Assessed
	Responsible person in the United Kingdom
	Mechanical heart valve
	Importer
	Unique Device Identifier
	Effective Orifice Diameter
	Inflow Orifice Diameter
	Cuff Outer Diameter
 (7440-48-4)	Contains hazardous substances (Cobalt CAS 7440-48-4)
	Double-ended handle, sizer rings

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Abbott Medical UK Limited