**EN: FemoStop™ II Plus
Femoral Compression System**

**DE: FemoStop™ II Plus
Femoralis-Kompressionssystem**

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PPL600235697 B
2022-08

EN: English

This Instructions for Use pertains to the FemoStop™ II Plus Femoral Compression System.

Table 1

Product Name	Basic UDI-DI	Sterilization Method
FemoStop™ II Plus Femoral Compression System	8717648CCS003S5	Electron-beam Radiation (R)

This Instructions for Use (IFU) as per European Regulation (EU) 2017/745 can be found on the Abbott eLabeling website at vascular.eIFU.abbott.

1.0 CAUTION

The FemoStop™ II Plus Femoral Compression System should only be used for compression of the femoral artery or vein after vessel cannulation by or on the order of physicians (or allied healthcare professionals, authorized by, or under the direction of, such physicians) trained in femoral artery or vein compression procedures.


Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in complications.

2.0 DEVICE DESCRIPTION

The FemoStop™ II Plus Femoral Compression System consists of a reusable arch and pump with a manometer and a disposable set including a pneumatic dome, belt, and stopcock. The pressure dome is placed over the vessel puncture site in the groin. The belt is placed around the patient. The dome applies a mechanical pressure over the vessel puncture site to induce hemostasis. The pressure of the dome is controlled by the reusable pump and the manometer. The arch and the belt provide counter pressure for the dome. The temperature range on the label represents the temperature for long-term storage. This product is not made with natural rubber latex. (Figure 1)

The FemoStop™ II Plus Bilateral Adapter can be used with the FemoStop II Plus Femoral Compression System to allow for simultaneous compression of punctures on both the left and right femoral artery/vein. The FemoStop II Plus Pump has a pressure accuracy of ± 10 mmHg (over the full range of 0-300 mmHg).

Figure 1: FemoStop II Plus Femoral Compression System

3.0 INTENDED PURPOSE

The FemoStop II Plus Femoral Compression System is intended to restore normal blood vessel hemostasis.

4.0 INTENDED CLINICAL BENEFIT

The FemoStop II Plus Femoral Compression System is utilized in diagnostic and interventional cardiovascular procedures to provide the benefit of restoring femoral vessel hemostasis.

5.0 INDICATIONS

The FemoStop™ II Plus Femoral Compression System is indicated for use in the compression of the femoral artery or vein after vessel cannulation.

6.0 CONTRAINDICATIONS

- The FemoStop II Plus Bilateral Adapter may be used with vascular access sites with a distance of 100 mm in standard configuration and 125-170 mm apart in the extended configuration. Vascular access sites outside this range are contraindicated.
- Severe peripheral vascular disease due to the risk of arterial thrombosis.
- Critical limb ischemia.
- Overlying skin necrosis and/or infection.
- Arterial injuries above or near the inguinal ligament.
- The inability to adequately compress due to e.g., coexisting very large hematomas, excessive pain or discomfort (despite anesthetics/analgesics).
- Patients not suitable for compression of their femoral artery due to their leg edema, femoral nerve compression or arterial obstruction.
- Femoral artery graft or vein graft due to the risk of damage.

7.0 WARNINGS

- Do not use the FemoStop II Plus Femoral Compression System if the packaging or sterile barrier has been previously opened or damaged or if the components appear to be damaged or defective.

- DO NOT RESTERILIZE OR REUSE. The FemoStop II Plus Disposable Set (pressure dome, stopcock, belt) is for single use only.
- Reuse, reesterilization and repackaging of the FemoStop II Plus Disposable Set may result in patient/user infections or product deterioration leading to reduced concentration of dome pressure, causing bleeding.
- Adequate compression may not be obtained in morbidly obese patients (Body Mass Index $\geq 40 \text{ kg/m}^2$).
- Incorrectly connecting the pump may introduce air into the vasculature, which may result in patient injury or death. To avoid this risk, ensure that the pump is connected to only the FemoStop™ II Plus Femoral Compression System.
- Release of thrombotic material may result in embolization which could lead to patient injury or death.
- If compression period is greater than 3 hours, briefly release pressure to decrease possible superficial tissue damage and start the procedure again with a new dome.
- If arterial/venous hemostasis is not achieved, significant bleeding may occur which could result in patient injury or death.
- While removing the sheath, ensure that the pressure applied is kept low, to avoid damage to the vessel or a "milking" effect. Allowance for slight bleeding at the site is preferred to preclude introduction of thrombus to the vessel.
- To prevent limb ischemia, do not leave artery completely blocked for more than 3 minutes. Check pedal pulses (dorsalis

pedis and posterior tibial) periodically to confirm whether or not flow remains in the vessels.

- To minimize the risk for arterial/venous fistula formation, arterial hemostasis should be maintained prior to removal of the venous sheath. Ensure that limb reperfusion has returned prior to removal of venous sheath.
- Do not apply pressure to a femoral artery stent due to risk of damage.
- The FemoStop II Plus Bilateral Adapter should only be used together with the FemoStop II Plus Femoral Compression System for the compression of the femoral artery or vein after vessel cannulation.

8.0 PRECAUTIONS

- When removing the sheath for an Intra-Aortic Balloon Pump (IABP), the Instruction for Use for the IABP should be followed as appropriate.
- For successful compression, the system must be snug and secure around the patient's hips before pressure is applied. Do not over-tighten the belt. The dome should only rest on, not press into, the groin area.
- For successful compression, the system must be correctly positioned throughout the procedure so that pressure is applied to the optimal point of vessel control. Ensure that the dome is positioned over the actual vessel puncture site.

- On very obese patients (Body Mass Index ≥ 30 kg/m²), it may be necessary to tighten the belt slightly more to enhance downward compression.
- When using the system on obese patients, fatty tissue may be displaced giving a false impression of a developing hematoma.
- Placement of the system may not be suitable on large patients, or patients with very wide hips as the belt may be too short. An abdominal strap/tape may be used to pull excessive adipose tissue away from dome.
- The target inflation pressure should be 10-20 mmHg above the systolic pressure, or higher if necessary, to control arterial bleeding. Exceeding pressures of 200 mmHg may indicate the need to tighten the belt or reposition the dome.
- Careful monitoring of the dome pressure during the initial period of use is recommended, as the elastic material of the dome may stretch slightly during the first few minutes. You may notice a slight drop in pressure. If this occurs, reinflate to initial pressure.
- If using the stopcock, ensure that it is in an open position when increasing or decreasing the pressure using the pump.
- Ensure that the control knob on the pump is closed when increasing the pressure and open when decreasing the pressure.
- Use of the FemoStop™ II Plus Femoral Compression System is not intended to replace careful monitoring of the patient's puncture site. The patient should not be left unattended during the time of compression.

- After use, this product may be a potential biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.
- The temperature range on the label represents the recommended storage temperature prior to use.
- Persons with known history of allergy to any of the components of this device listed below may suffer an allergic reaction to this device. Prior to use on the patient, the patient should be counselled on the materials contained in the device and a thorough history of allergies must be discussed. This device contains an ethyl vinyl acetate (EVA) copolymer.

9.0 LIMITATIONS

As with any catheter-based procedures, vessel damage and / or device breakage is a possibility. Observe Warnings and Precautions at all times when using this device and be prepared for necessary intervention as per facility protocol.

10.0 ADVERSE EVENTS

Potential adverse events associated with the use of this device may include, but are not limited to the following:

- Allergic reaction, or hypersensitivity to device components
- Arteriovenous fistula
- Blistering of the skin / skin abrasion
- Bruising / ecchymoses
- Femoral artery and / or vein thrombosis

- Hematoma
- Hemorrhage / bleeding
- Hypotension
- Infection
- Inflammation
- Ischemia
- Nerve damage
- Pseudoaneurysm
- Tissue injury
- Tissue necrosis

11.0 REPORTING OF ADVERSE EVENTS

The user and / or patient should report any serious incident that has occurred in relation to the FemoStop II Plus Femoral Compression System to the manufacturer and the competent authority of the European Member State in which the user and / or patient is established.

12.0 PATIENT COUNSELING AND PATIENT INFORMATION

Physicians should consider the following in counseling patients about this device:


- Discuss the risks of femoral artery / vein complications in relation to vascular access and / or use of vascular access closure.

- Discuss the risks associated with the use of vascular access closure device in comparison to manual compression and / or surgical cutdown.
- Discuss the risks of allergic reaction or hypersensitivity to device components.

13.0 DIRECTIONS FOR USE

1. Prior to use, carefully examine the puncture site area for any pre-existing hematomas and ensure proper placement of the arch. This may be more difficult in larger patients due to excessive adipose tissue in this region.
2. Note patient's current blood pressure to determine proper initial inflation pressure.
3. Check pedal pulse(s) (dorsalis pedis and / or posterior tibial) to determine what is the baseline pedal pulse(s) for the patient.
4. Prior to sheath removal, if applicable, place belt under and around the patient's hips, so that it is pulled up equally on both sides and is directly in line with the puncture sites. Check that the belt has not become twisted or folded under the patient (see Figure 2). Thread the belt through the locks at one end of the arch by fully compressing the lever.

Figure 2: Place the Belt

5. If the FemoStop II Plus Bilateral Adapter is being used follow the steps below.

a) Place the FemoStop II Plus Bilateral Adapter on the FemoStop II Plus Compression Arch by unlocking the thumb lock and sliding on to the body of the arch. Do not lock the FemoStop II Plus Bilateral Adapter in place.

b) Palpate both sides where the introducer enters the vessel wall.



c) Hold the FemoStop II Plus Compression Arch over the patient and adjust the FemoStop II Plus Bilateral Adapter to align with both arteriotomy/venotomy sites.

NOTE: This distance between arteriotomy/venotomy sites may not be the same as the distance between the skin punctures.

NOTE: Take care not to contaminate the sterile puncture sites.

d) Lock the FemoStop II Plus Bilateral Adapter in place by pushing the thumb lock down onto the notch.

Figure 3: Bilateral Adapter Distance


6. Thread the belt through the locks at one end of the arch by fully compressing the levers.
7. Remove the pneumatic pressure dome from the packaging. Do not remove the protective peel-back lid from the sterile surface of the dome. Attach the dome to the compression

arch and lock by rotating the pneumatic pressure dome (counterclockwise). Ensure that the dome is locked in place (Figure 4).

a) If using the FemoStop II Plus Bilateral Adapter, attach a second dome and lock by rotating the dome. Do not remove the protective peel-back lid from the sterile surface of the dome. Confirm that the distance between the center of the two pneumatic pressure domes is the same as the distance between the arteriotomy/venotomy puncture sites.

Figure 4: Attach the Pneumatic Pressure Dome to the Compression Arch


8. Remove the protective lid from the sterile surface of the dome(s). Be careful to maintain the sterility of the dome(s).
9. Rest the sterile surface of the dome(s) on the puncture site(s). Tighten the straps of the device to the point where it does not shift while sheath is removed, if applicable. Withdraw

the introducer hubs about 2 cm or just enough to clear the rim of the dome. Palpate the entrance site of the sheath(s), or access site if no sheath is present, at the femoral artery/vein and place the center of the dome superior and medial to skin incisions, in order to ensure that the dome is positioned over the actual puncture site (Figure 5).


NOTE: The specific position of the arterial puncture site relative to the skin incision will vary depending on the angle of insertion of the introducer sheath and the distance from the skin to artery.

Figure 5: Place the Dome Over the Puncture Site


10. Thread the loose end of the belt through the open end of the FemoStop II Plus Compression Arch. Adjust the belt to a snug but comfortable fit around the patient. The arch should lie level and squarely across the groin area (Figure 6).

Figure 6: Adjust the Belt

11. Attach the FemoStop II Plus pump to the pneumatic pressure dome using the free end of the connecting tubing. Check connections to ensure a tight fit. The disposable stopcock should be connected between the pump tubing and the dome. Ensure that the stopcock is positioned so that its handle is accessible and that it is in an open position when increasing or decreasing the pressure using the pump. (Figure 7).

Figure 7: Connect the Stopcock

12. Close the control knob of the FemoStop II Pump by turning the knob clockwise and begin to inflate and remove the sheath, if present (Figure 8). When removing the arterial sheath, inflate the dome to 60 or 80 mmHg. When removing the venous sheath inflate the dome to 20 or 30 mmHg (Figure 9).

Figure 8: Increase Pressure



Figure 9: Remove the Sheath

13. As soon as the sheath has been removed (if one is present), continue to increase the pressure to 10-20 mmHg above systolic or higher as necessary to maintain initial hemostasis.
14. Maintain hemostasis pressure for about 1-3 minutes. Do not leave artery completely blocked for more than 3 minutes to prevent limb ischemia. Refer to Warnings.
15. Lower the FemoStop II Plus Pump pressure by opening the pressure gauge valve until pedal pulses are palpable, limb perfusion is restored, and hemostasis is maintained.
16. If the FemoStop II Plus Bilateral Adapter is being used, wait until the pressure has been reduced below systolic, patient pedal pulse is confirmed, and hemostasis is still maintained on the first side, before removing the sheath on the other side, if applicable. Repeat steps 12 – 15. Apply the following steps to both pneumatic pressure domes.
17. Continue the maintenance pressure according to facility guidelines specified for the type of patient, anticoagulant level, the interventional procedure performed, and the sheath French size used.
18. Once the above duration of time has elapsed, gradually lower the pressure in steps of 10–20 mmHg, as long as hemostasis is maintained. Observe the site for 2–3 minutes at each step of deflation. If hemostasis is not maintained, reinflate the pneumatic pressure dome to a pressure where hemostasis is confirmed, pedal pulses are palpable and limb perfusion is

restored. If hemostasis is maintained, continue until the pneumatic pressure dome is completely deflated. The FemoStop™ II Plus Femoral Compression System still permits visualization of the skin incision site, and adjustment of pressure and position.

19. Leave completely deflated pneumatic pressure dome in place for at least 2-3 minutes and confirm hemostasis. Once hemostasis is confirmed, proceed to remove the system as described below.

WARNING: If compression period is greater than 3 hours, briefly release pressure to decrease possible superficial tissue damage. Inappropriately long compression and/or immobilization may increase the risk for thrombosis or embolization which could lead to patient injury or death.

NOTE: If significant bleeding occurs, it is necessary to raise the pressure above the systolic pressure and start the procedure again with a new dome.


20. Carefully loosen the belt on one side of the patient (without totally removing it from the arch) and move the FemoStop II Plus Compression Arch so that the pneumatic pressure dome is clear of the puncture site. This will allow better visibility of the sites and still provide the control needed to quickly reposition and recompress the sites if needed.

21. Once complete hemostasis has been established remove the FemoStop™ II Plus Femoral Compression System from the patient by disconnecting the FemoStop II Plus Pump from the pneumatic pressure dome and removing the FemoStop II Plus Compression Arch completely off of the belt.

22. Detach the pneumatic pressure dome(s) from the FemoStop II Plus Compression Arch and FemoStop II Plus Bilateral Adapter. Discard the belt and the pneumatic pressure dome(s). If using the FemoStop II Plus Bilateral Adapter, detach the bilateral adapter from the FemoStop II Plus Compression Arch. Clean and store the FemoStop II Plus Compression Arch, FemoStop II Plus Pump and FemoStop II Plus Bilateral Adapter per Section 15.0.

14.0 COMPRESSION TIME EXAMPLE

The following chart provides an example of recommended times for hemostasis. Follow facility guidelines.

- A. Minutes
- B. Systolic BP
- C. Diastolic BP
- 1. 60-80 mmHg, Remove Arterial Sheath
- 2. Max 20 mmHg Suprasystolic
- 3. Palpable Pedal Pulse
- 4. -20 mmHg
- 5. -20 mmHg
- 6. ≈30 mmHg

7. 0 mmHg

* After three max minutes: lower to mean arterial pressure.
Check pedal pulse.

** The length of compression depends on factors which vary, such as sheath size and anti-coagulation status.

NOTE: If there is a venous sheath, inflate the dome to 20-30 mmHg and remove. To minimize the risk of AV-fistula formation, venous hemostasis should be achieved prior to removal of the arterial sheath

15.0 CLEANING AND MAINTENANCE

These instructions have been validated as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials, and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Limitations on Processing – The service life of the reusable components is 2 years or 30 uses, whichever comes first.

The following information applies to:

- FemoStop™ II Plus Compression Arch **Catalogue/REF 11168**
- FemoStop™ II Plus Bilateral Adapter (Optional)

Catalogue/REF 11169

CAUTION: When cleaning the FemoStop II Plus Compression Arch, avoid using cloths or solutions that include bleach-based disinfectants (sodium hypochlorite).

Cleaning / Disinfection – Using gloved hands, remove the FemoStop II Plus Compression Arch from the FemoStop II Plus Bilateral Adapter, if attached. Wipe down all surfaces to remove visible soil using a quaternary ammonium-based germicidal wipe. Keep components wet per manufacturer's instructions for contact duration. Allow the components to air dry.

Maintenance / Inspection – If visible/solid debris cannot be removed from the device or if there are visual signs of wear or damage that prevent the device from being used safely, do not use the device further and discard per the DISPOSAL INSTRUCTIONS section of this Instructions for Use.

The following information applies to:

- FemoStop™ II Plus Pump **Catalogue/REF 11170**

CAUTION: When cleaning the pump, avoid using cloths or solutions that include quaternary ammonium compounds (ammonium chlorides) or glutaraldehyde-based disinfectants.

Cleaning / Disinfection – Using gloved hands, wipe down all surfaces to remove visible soil using a bleach-based germicidal wipe. Keep component wet per manufacturer's instructions for contact duration. Allow the component to air dry.

Maintenance/Inspection – Routinely inspect the gauge for wear or other damage. Do not use if you see signs of damage. Periodically check calibration of the gauge by inspecting the pointer of the gauge at zero pressure and verifying that the pointer is within the oval surrounding the zero-pressure graduation on the dial. If the pointer is fully outside of the

tolerant zone, shown in the illustration below as the darkly shaded area, the device should not be used.

16.0 HOW SUPPLIED

The FemoStop™ II Plus Femoral Compression System consists of:

- FemoStop™ II Plus Disposable Set **Catalogue/REF 11166**
- FemoStop™ II Plus Pump **Catalogue/REF 11170**
- FemoStop™ II Plus Compression Arch **Catalogue/REF 11168**
- FemoStop™ II Plus Bilateral Adapter (Optional) **Catalogue/REF 11169**

Storage – Storage is recommended in a dry, dark, cool place. See the device label for recommended storage conditions.

Additional information for the FemoStop II Plus Disposable Set:

This single-use medical device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user. Do not resterilize.

Sterile – Sterilized using electron-beam radiation. Do not use if the package is open or damaged.

17.0 DISPOSAL INSTRUCTIONS

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

18.0 PRODUCT INFORMATION DISCLOSURE

Abbott has exercised reasonable care in the manufacturing of this device. Abbott excludes all warranties, whether expressed or implied, by operation of law or otherwise, including but not limited to, any implied warranties of merchantability of fitness, since handling and storage of this device, as well as factors

relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond the control of Abbott directly affect this device and the results obtained from its use. Abbott shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from the use of this device. Abbott neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

Graphical Symbols for Medical Device Labeling

Symbol	Definition
	Sterilized using irradiation
	Catalogue number
	Batch code
	Consult instructions for use or consult electronic instructions for use
	Do not reuse

Symbol**Definition**

Do not resterilize

Do not use if package is damaged and consult instructions for use

Use by date

Packaging unit

Keep away from sunlight

Keep dry

Nonsterile

Symbol**Definition**

Temperature Limitations

REMOVE PROTECTIVE LID BEFORE USE

Remove protective lid before use

Femoral Compression System

Femoral Compression System

Compression Arch

Compression Arch

Disposable Set

Disposable Set

Pump

Pump

Bilateral Adapter

Bilateral Adapter

Date of Manufacture

Manufacturer

EC REP

Authorized Representative in the European Community

Symbol

Definition

R_x
ONLY

CAUTION: Federal law (U.S.) restricts this device to sale by or on the order of a physician.

Importer in the European Community

Unique device identifier

Medical device

Single sterile barrier system

Responsible person in the United Kingdom

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Symbole	Définition
 The pictogram consists of a rectangular box divided into two sections. The left section contains the word "STERILE" in bold, uppercase letters. The right section contains the letter "R" in bold, uppercase letters.	Stérilisé par rayonnement
 The pictogram consists of a rectangular box containing the word "REF" in bold, uppercase letters.	Référence catalogue
 The pictogram consists of a rectangular box containing the word "LOT" in bold, uppercase letters.	Code du lot
 The pictogram consists of an open book icon with a lowercase letter "i" inside a circle on the right page.	Consulter le mode d'emploi ou le mode d'emploi électronique
 The pictogram consists of a circle containing a crossed-out infinity symbol (∞).	À usage unique
 The pictogram consists of a circle containing the number "2" above the word "STERILIZE" in uppercase letters.	Ne pas restériliser

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Symbool	Definitie
 The symbol consists of two adjacent rectangular boxes. The left box contains the word "STERILE" in bold, uppercase letters. The right box contains the letter "R" in bold, uppercase letters.	Gesteriliseerd met bestraling
 The symbol consists of a rectangular box containing the word "REF" in bold, uppercase letters.	Catalogusnummer
 The symbol consists of a rectangular box containing the word "LOT" in bold, uppercase letters.	Partijcode
 The symbol consists of an open book icon with a lowercase letter "i" inside the right page.	Raadpleeg de gebruiksaanwijzing of raadpleeg de digitale gebruiksaanwijzing
 The symbol consists of a circle with a diagonal slash over a symbol that represents reuse (two arrows forming a loop).	Niet hergebruiken
 The symbol consists of a circle with a diagonal slash over a symbol that represents re-sterilization (a number "2" above the word "STERILIZE").	Niet opnieuw steriliseren
 The symbol consists of a circle with a diagonal slash over a symbol that represents damaged packaging (a box with a lightning bolt striking it).	Niet gebruiken als de verpakking beschadigd is en raadpleeg de

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Argonstraat 1, 6422 PH
Heerlen
The Netherlands
TEL: +32-(0) 45-5435700

Abbott Medical UK Limited

Elder House
Central Boulevard
Blythe Valley Business Park
Solihull B90 8AJ
United Kingdom

**Abbott Vascular International
BVBA**

Park Lane, Culliganlaan 2B
1831 Diegem, BELGIUM
TEL: +32 2 714 14 11
FAX: +32 2 714 14 12

CUSTOMER SERVICE

TEL: (800) 227-9902
FAX: (800) 601-8874
Outside USA
TEL: +1 951-914-4669
Outside USA
FAX: +1 951-914-2531

Abbott

Reference Abbott website for patent markings:

www.abbott.com/patents

™ Indicates a trademark of the Abbott group of companies.

© 2022 Abbott. All Rights Reserved.