Masters Series™/Masters Series™ HP Mechanical Heart Valve

EN: English	1
BG: Български	15
CS: Čeština	30
DA: Dansk	44
DE: Deutsch	58
EL: Ελληνικά	73
ES: Español	88
ET: Eesti keel	102
FI: Suomi	116
FR : Français	130
HR: Hrvatski	144
HU: Magyar	158
IT: Italiano	172
LT: Lietuvių	187
LV: Latviski	201
NL: Nederlands	215
NO: Norsk	230
PL: Polski	244
PT: Português	258
RO: Română	272
SK: Slovenčina	286
SL: Slovenščina	300
SV: Svenska	314
TR: Türkçe	328

Instructions for Use

This warranty is given by Abbott Medical. To make inquiries regarding this warranty, use the contact information for either the manufacturer or the Australian Sponsor.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

The benefits to you given by this warranty are in addition to your other rights and remedies under the Australian Consumer Law.

Abbott Medical warrants to the buyer that the goods, for a period equal to the validated shelf life of the goods (the "Warranty Period"), shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Unless otherwise obligated by law or mandatory consumer guarantees, under this warranty Abbott Medical will, at its option, replace or repair the goods at its factory, if the buyer, at the buyer's expense, returns the goods within the Warranty Period to Abbott Medical to the address below and after Abbott confirms that the goods are defective.

Abbott Medical

Attn: Postmarket Surveillance

One St. Jude Medical Drive

St. Paul, MN 55117 USA

+1 651 756 2000

To the maximum extent permitted by law, if a mandatory term is implied by law or a mandatory consumer guarantee applies to the goods, and the goods are not of a kind ordinarily acquired for personal, domestic or household use or consumption, Abbott Medical's liability for the breach of the term or guarantee is limited to, at Abbott Medical's option, either replacement or repair of the goods or payment of the costs of replacing or repairing the goods.

EXCEPT FOR MANDATORY TERMS IMPLIED BY LAW, MANDATORY CONSUMER GUARANTEES, OR AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED.

See the Terms and Conditions of Sale for further information.

[™] Indicates a trademark of the Abbott group of companies.

[‡] Indicates a third-party trademark, which is property of its respective owner.

Pat. <http://www.abbott.com/patents>

© 2024 Abbott. All Rights Reserved.

EN: English

Masters Series™/Masters Series™ HP Mechanical Heart Valve

Instructions for Use


Device Description and Performance Characteristics

Masters Series™ Mechanical Heart Valves (Masters and Masters HP) are rotatable, bi-leaflet replacement heart valves manufactured from pyrolytic carbon designed for implantation in the aortic or mitral position. Masters Series™ Mechanical Heart Valves are designed for intra-annular cuff placement.

Masters Series™ HP Mechanical Heart Valves are designed for supra-annular cuff placement. A sewing cuff of knitted double velour polyester or knitted polytetrafluoroethylene (PTFE) is secured around the periphery of the valve to allow attachment of the valve to the native valve annulus. The Masters Series™ Mechanical Heart Valve includes a metallic ring assembly within the sewing cuff that allows the valve to be rotated *in-situ*.

Hemodynamic Plus (HP) models have been designed to maximize effective orifice area. This is made possible by removing the sewing cuff from the annulus.

Figure 1. Masters Series™ Standard and HP valves implanted in the aortic position

1. Tissue annulus diameter
2. Outer diameter


Table 1. Reference Dimensions

Model		Tissue Annulus Diameter (mm)	Geometric Orifice Area (cm ²)	Implant Height Open (mm)	Overall Height Open (mm)	External Sewing Ring Diameter (mm)
Standard Polyester Cuff	Expanded Polyester Cuff					
Masters Series™ Valves: HP Series Cuffed Valves						
Aortic						
17 AHPJ-505	17 AEHPJ-505	17	1.6	4.6	9.8	22.0/22.5
19 AHPJ-505	19 AEHPJ-505	19	2.1	5.2	10.6	24.0/24.0
21 AHPJ-505	21 AEHPJ-505	21	2.6	6.0	12.1	26.0/26.5
23 AHPJ-505	23 AEHPJ-505	23	3.1	6.2	13.3	28.5/29.0
25 AHPJ-505	25 AEHPJ-505	25	3.8	6.7	13.9	30.5/30.5
27 AHPJ-505	27 AEHPJ-505	27	4.4	7.5	14.6	32.5/33.0
Mitral						
17 MHPJ-505	--	17	1.6	2.4	9.8	22.0
19 MHPJ-505	--	19	2.1	2.7	10.6	24.0
21 MHPJ-505	--	21	2.6	3.0	12.1	26.0
23 MHPJ-505	--	23	3.1	3.2	13.3	28.5
25 MHPJ-505	--	25	3.8	3.8	13.9	30.5
27 MHPJ-505	--	27	4.4	4.6	14.6	32.0
Masters Series™ Valves: Polyester Cuff						
Aortic						
19 AJ-501	19 AECJ-502	19	1.6	3.5	9.8	22.5/23.0
21 AJ-501	21 AECJ-502	21	2.1	3.9	10.6	24.5/24.5

Table 1. Reference Dimensions

Model		Tissue Annulus Diameter (mm)	Geometric Orifice Area (cm ²)	Implant Height Open (mm)	Overall Height Open (mm)	External Sewing Ring Diameter (mm)
Standard Polyester Cuff	Expanded Polyester Cuff					
23 AJ-501	23 AECJ-502	23	2.6	5.2	12.1	27.0/28.0
25 AJ-501	25 AECJ-502	25	3.1	5.5	13.3	29.5/30.0
27 AJ-501	27 AECJ-502	27	3.8	5.8	13.9	31.0/32.0
29 AJ-501	29 AECJ-502	29	4.4	6.5	14.6	33.5/34.0
31 AJ-501	31 AECJ-502	31	5.1	7.3	16.0	35.5/37.0
Mitral						
19 MJ-501	19 MECJ-502	19	1.6	3.8	9.8	23.0/24.0
21 MJ-501	21 MECJ-502	21	2.1	3.7	10.6	25.0/25.5
23 MJ-501	23 MECJ-502	23	2.6	4.6	12.1	28.0/28.5
25 MJ-501	25 MECJ-502	25	3.1	5.2	13.3	30.0/30.5
27 MJ-501	27 MECJ-502	27	3.8	5.8	13.9	32.0/32.0
29 MJ-501	29 MECJ-502	29	4.4	6.7	14.6	34.5/35.0
31 MJ-501	31 MECJ-502	31	5.1	7.1	16.0	36.5/37.0
33 MJ-501	33 MECJ-502	33	5.1	8.3	16.0	38.0/38.5
35 MJ-501	--	35	5.1	8.2	16.0	40.5/NA
37 MJ-501	--	37	5.1	7.8	16.0	41.5/NA
Masters Series™ Valves: PTFE Cuff						
Aortic						
19 ATJ-503	--	19	1.6	3.5	9.8	24.0
21 ATJ-503	--	21	2.1	3.9	10.6	26.0
23 ATJ-503	--	23	2.6	5.2	12.1	28.5
25 ATJ-503	--	25	3.1	5.5	13.3	30.5
27 ATJ-503	--	27	3.8	5.8	13.9	32.5
29 ATJ-503	--	29	4.4	6.5	14.6	35.0
31 ATJ-503	--	31	5.1	7.3	16.0	37.0
Mitral						
19 MTJ-503	--	19	1.6	3.8	9.8	24.5
21 MTJ-503	--	21	2.1	3.7	10.6	26.5
23 MTJ-503	--	23	2.6	4.6	12.1	29.0
25 MTJ-503	--	25	3.1	5.2	13.3	31.0
27 MTJ-503	--	27	3.8	5.8	13.9	33.0
29 MTJ-503	--	29	4.4	6.7	14.6	35.5
31 MTJ-503	--	31	5.1	7.1	16.0	37.5
33 MTJ-503	--	33	5.1	8.3	16.0	39.5

Figure 2. Mitral valve measurements
Masters Valve

1. Tissue annulus diameter
2. Implant height, open
3. Overall height, open

Masters HP Valve




Figure 3. Aortic valve measurements
Masters Valve

1. Tissue annulus diameter
2. Overall height, open
3. Implant height, open

Masters HP Valve

Leaflets and Orifice Ring

Pyrolytic carbon coats the graphite substrates of the leaflets and orifice ring. Pyrolytic carbon was chosen for its biocompatibility and durability. The graphite substrate for the leaflets is tungsten-impregnated for radiopacity. X-ray or fluoroscopy may be used to visualize the valve. To obtain good x-ray visualization of the leaflets, orient the x-ray beam either parallel to the leaflet pivot mechanism axis, or perpendicular to the plane of the valve orifice. The orifice ring may be rotated to the surgeon's *in situ* preference after the valve is sutured in place.

Masters Series™ Valve Sewing Cuff

The following Masters Series™ Valve sewing cuff configurations are available:


- Standard Cuff, Polyester
- Expanded Cuff, Polyester
- Standard Cuff, PTFE
- HP Cuff, Polyester
- Expanded HP Cuff, Polyester

The double-velour knitted polyester fiber in the polyester cuff allows for rapid, controlled endothelial ingrowth over the entire sewing cuff. Masters Series™ valves with polyester sewing cuffs have suture markers that may be used to provide reference points for the valve's orientation or the placement of sutures. Aortic polyester cuff models have three suture markers; mitral polyester cuff models have four suture markers. Masters Series™ Valves with PTFE sewing cuffs do not have suture markers.

Figure 4. Three suture markers on aortic sewing cuff (polyester only)

Figure 5. Four suture markers on mitral sewing cuff (polyester only)

Expanded Cuff is constructed with more polyester cuff material and provides a bulkier configuration with additional suturing material. This translates into a larger or similar external sewing ring diameter for expanded cuff configuration in comparison to the standard cuff configuration.

Holder/Rotator

A holder/rotator is attached to each Masters Series™ Mechanical Heart Valve. Use this holder/rotator to rotate the valve *in situ*.

Intended Purpose

The Masters Series™ Mechanical Heart Valve is intended to replace native aortic or mitral heart valves or previously implanted aortic or mitral prosthetic valves.

Indications

The Masters Series™ Mechanical Heart Valve is indicated for use as a replacement valve in patients with a diseased, damaged, or malfunctioning aortic or mitral heart valve or as a replacement for a prosthetic heart valve.

Target Patient Population

Patients with a diseased, damaged, or malfunctioning aortic or mitral valve.

Intended User

Physicians trained in surgical mitral and/or aortic valve replacement.

The device implantation should be performed in a standard operating room equipped for open-heart surgery.

Intended Clinical Benefits

Based on clinical evidence, the intended clinical benefits for the Masters Series™ Valves are:

- Low rate of reintervention
- Improved symptoms and functional capacity as assessed by New York Heart Association [NYHA] Classification
- Improved hemodynamic valve performance by decreasing valvular stenosis and/or valvular regurgitation
- Low rates of the following adverse events: valve related mortality, valve dysfunction, valve thrombosis, major paravalvular leak, hemolysis, major hemorrhage, thromboembolism, and endocarditis

Summary of Safety and Clinical Performance (SSCP)

A summary of the safety and clinical performance (SSCP) for this device is available at <https://ec.europa.eu/tools/eudamed>. Search for the device using the basic UDI-DI listed in the following table. This is the SSCP location after the launch of the European Database on Medical Devices/Eudamed.

Table 2. Masters Series™ Mechanical Heart Valves Basic UDI-DI

Basic UDI-DI	Basic UDI-DI Description
5415067SHV2001PL	Surgical Heart Valve Products, Mechanical Heart Valves

Limitations

No known limitations.

Contraindications

The Masters Series™ Mechanical Heart Valve is contraindicated for individuals unable to tolerate anticoagulation therapy.

Warnings

- For single use only. Attempts to reuse the valve may result in valve malfunction, inadequate sterilization, or patient harm.
- Do not use if:
 - The valve has been dropped, damaged, or mishandled in any way.
 - The expiration date has elapsed.
 - The tamper-evident container seal or inner/outer tray seals are damaged, broken, or missing.
- Remove any residual tissue that may impair valve size selection, correct seating of the valve, rotation of the valve, or leaflet motion.
- Proper valve size selection is crucial. Do not oversize the valve. If the native annulus measurement falls between two Masters Series™ Mechanical Heart Valve sizes, use the smaller Masters Series™ Mechanical Heart Valve.
- Use only Sizer Set Model 905.
- The outer tray is not sterile, and should not be placed in the sterile field.
- To minimize direct handling of the valve during implantation, do not remove the holder/rotator until the valve has been seated in the annulus.
- Do not use hard or rigid instruments to test leaflet mobility, as this may result in structural damage to the valve or thromboembolic complications. Use a Model LT100 leaflet tester to gently test valve leaflet mobility.
- Place sutures in the outer half of the valve sewing cuff.

- Never apply force to the valve leaflets. Force may cause structural damage to the valve.
- Use only Valve Holder/Rotators from Abbott Medical to perform valve rotation. Use of other instruments could result in structural damage. The valve holder/rotator is intended for single use only and should be discarded after surgery.
- The two retention sutures on the valve holder/rotator must be cut and removed before the Masters Series™ Mechanical Heart Valve can be rotated.
- Do not pass catheters or other instruments through Masters Series™ mechanical heart valves. This could result in scratched or damaged valve components, leaflet fracture, or dislodgment.
- Cut suture ends short, especially in the vicinity of the pivot guards, to prevent leaflet impingement.
- Do not use cutting edge or cutting point needles in the sewing cuff. Use of standard round or taper point needles is recommended to prevent cutting the polyester fibers.

Precautions

- **Do not touch the prosthetic valve unnecessarily, even with gloved hands. This may cause scratches or surface imperfections that may lead to thrombus formation.**
- **Be careful not to cut or tear the valve sewing cuff when removing the identification tag and the holder/rotator from the Masters Series™ Mechanical Heart Valve.**
- **Before placing sutures in the valve sewing cuff, verify that the valve is mounted correctly on the valve holder/rotator.**
- **To avoid structural damage, the mitral valve must be rotated in the fully open position.**
- **To avoid structural damage, the aortic valve must be rotated in the fully closed position.**
- **To minimize rotational torque, verify that the valve holder/rotator is properly seated in the valve, and that the valve holder handle is perpendicular to the valve.**
- **Remove any loose suture or thread, which may be a source of thrombus or thromboembolism.**
- **Valve recipients should be maintained on anticoagulant therapy post procedure, per institutional standards and established guidelines, except when contraindicated, using individualized treatment as determined by their physician.**

Potential Adverse Events

Complications associated with replacement mechanical heart valves include, but are not limited to:

- adjacent cardiac structure interference
- death
- heart block requiring pacemaker implant
- heart failure
- hemolysis
- hemorrhagic complications secondary to anticoagulation therapy
- infections
- myocardial infarction
- prosthetic failure
- stroke
- thrombus or thromboembolism
- unacceptable hemodynamic performance
- valve dehiscence

Any of these complications may require reoperation or explantation of the device.

Device Lifetime

The product is designed and tested to support a minimum lifetime of 10 years. In addition, long-term clinical data published in the literature indicates Abbott's mechanical heart valves perform as intended for at least 20 years.

Magnetic Resonance (MR) Safety Information

A person with the Masters Series™ Mechanical Heart Valves and Valved Grafts may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.

Device Names	Masters Series™ Mechanical Heart Valves Masters Series™ Aortic Valved Graft Masters Series™ HP Valved Graft
Static Magnetic Field Strength (B0)	1.5T or 3.0T
Maximum Spatial Field Gradient	≤20 T/m (2000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Body coil
Operating Mode	Normal Operating Mode
Maximum Whole Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	N/A
Scan Duration	2 W/kg whole body average SAR for 15 minutes of continuous scanning
MR Image Artifact	The presence of this implant may produce an image artifact.

How Supplied

Packaging

The Masters Series™ Mechanical Heart Valve has been steam sterilized. The valve remains sterile until the "USE BY" date on the packaging, if package integrity is maintained.

The package includes:

- Sealed non-sterile outer tray
- Sealed sterile inner tray
- One (1) Masters Series™ Mechanical Heart Valve with identification tag, mounted on a plastic holder/rotator
- One (1) disposable support collar
- One (1) Instructions for Use
- One (1) Medical Device Registration Form with attached Patient Identification Card and return envelope

Storage

There are no special storage requirements.

Valve Resterilization

If resterilization of the Masters Series™ Mechanical Heart Valve is necessary, use only the recommended steam cycles and follow the instructions below. **The valve must not be resterilized more than one (1) time.**

1. Remove the inner tray containing the valve and also remove the inner tray cover.
2. Place the inner tray inside a steam-permeable sterilization bag or wrap.

Table 3. Recommended Sterilization Cycle Parameters

Vacuum Cycle	Pre-Vacuum Steam (Option 1)	Pre-Vacuum Flash (Option 2)
Purge Time:	6 minutes	6 minutes
Pulses:	2	2
Pulse Pressure:	204.7 kPa (absolute)	204.7 kPa (absolute)
Pulse Vacuum:	23.4 kPa (absolute)	23.4 kPa (absolute)
Sterilize Time:	28 minutes	7 minutes
Sterilize Temperature:	122°C	132°C
Post-Vacuum:	13.3 kPa (absolute)	13.3 kPa (absolute)
Vacuum Dry Time:	10 minutes	10 minutes

Accessories

Sizer Set

Use Sizer Set Model 905 to determine the appropriate Masters Series™ Mechanical Heart Valve size. The sizer set contains valve sizers, and a valve holder handle, Model 905-HH. Components of the sizer set have been designed and tested for repeated use. However, if visual signs of deterioration become apparent, do not use the sizer set, and contact Abbott Customer Service for a replacement. The sizer set is supplied non-sterile. Refer to the Model 905 Sizer Set Instructions for Use for cleaning and sterilization instructions.

Holder Handles

The following holder handles may be purchased separately. Consult the appropriate instrument Instructions for Use for full product description, and relevant cleaning and sterilization information. Instruments have been designed and tested for repeated use, however, if visual signs of deterioration become apparent, do not use the instrument, and contact Customer Service for a replacement.

Table 4. Optional Holder Handles

Description	Model
Flexible holder handle	905-HH
Rigid holder handle	905-RHH
Rigid mitral holder handle	905-MHH

Replacement Valve Holder/Rotators

Replacement aortic and mitral valve holder/rotators may be purchased separately. See the table below for model numbers and descriptions:

Table 5. Replacement Valve Holder/Rotators

Description	Model
Aortic (Standard)	A-RHR*
Aortic (HP)	AHP-RHR*
Mitral (Standard)	M-RHR*

Leaflet Tester

Use Leaflet Tester model LT100 to test leaflet mobility.

Figure 6. LT100 Leaflet Tester

Directions for Use

Removing the Native Valve

Excise the native valve and prepare the annulus for the valve replacement.

WARNING: Remove any residual tissue that may impair valve size selection, correct seating of the valve, rotation of the valve, or leaflet motion.

Sizing

Use Mechanical Heart Valve Sizer Set Model 905 to determine the appropriate Masters Series™ Mechanical Heart Valve size. The Model 905 sizer is a double-ended tool with a cylindrical annulus sizing end and a flanged HP sewing cuff replica end. Use the cylindrical end of the tool to determine the appropriate Masters Series™ Mechanical Heart Valve size. The cylindrical end of the sizer must pass readily without resistance through the annulus. When implanting HP valve models, use the flanged end of the sizer to visualize supra-annular placement of the HP sewing cuff.

CAUTION:For double-ended sizers only: Do not pass flanged sizers through the annulus.

Figure 7. Cylindrical end of the 905 sizer in the annulus



Figure 8. Flanged valve cuff replica end of the 905 sizer above the annulus

WARNING: Proper valve size selection is crucial. Do not oversize the valve. If the native annulus measurement falls between two Masters Series™ Mechanical Heart Valve sizes, use the smaller prosthetic valve size.

Pre-Implant Handling


Circulating Nurse (Non-Sterile)

1. Remove the non-sterile outer tray from the product box.

WARNING: Do not use the device if the expiration date has elapsed.

2. Verify that the catalog number and serial number on the outer tray are identical to those on the box label. If the information is not identical, do not use the device, and contact Customer Service as soon as possible.
3. Open the outer tray as shown in the following figure. Do not touch the sterile inner tray cover.

Figure 9. Remove the cover from the outer tray

4. Holding the outer tray by the bottom, present the inner tray to the sterile scrub nurse or surgeon. Do not touch the sterile inner tray.

5. Complete the patient [medical device] registration information, as described in the “Patient Registration” section.

Scrub Nurse (Sterile)

1. Remove the sterile inner tray from the outer tray. Avoid contact with the exterior of the outer tray, which is not sterile.

Figure 10. Remove the inner tray

WARNING: Examine the packaging carefully, and ensure the tray is unopened and unbroken. If you note any damage, do not use the device, and contact customer service.

2. Holding the inner tray with the cover upward, grasp the tab and pull back to completely remove the inner tray cover.

Figure 11. Remove the inner tray cover


3. Press the sterilized mechanical valve holder handle into the valve holder/rotator. The valve holder handle should be straight when inserted into the valve holder. Verify that the valve holder handle is securely attached to the valve.

Figure 12. Insert the holder handle in the valve holder

4. To remove the valve from the inner tray, firmly lift the valve holder handle and the valve support collar.


Figure 13. Remove the valve and collar from the inner tray

NOTE: Always use the valve holder handle to remove the valve from the inner tray.

5. Prior to implanting the valve, remove the support collar from the valve by placing two fingers under the collar, pinching with the thumb and pulling back gently.

Figure 14. Remove the collar from the valve prior to implantation


6. An identification tag is sutured to the valve sewing cuff. Verify that the valve size and model number on the identification tag are identical to the valve size and model number listed on the packaging. If the size and model number do not match, do not use the valve.
7. Remove the identification tag and its suture prior to implantation. Retain the identification tag for the patient record.

CAUTION: Be careful not to cut or tear the valve sewing cuff when removing the identification tag from the Masters Series™ Mechanical Heart Valve.

Valve Implantation

1. Place sutures in the annulus. Suturing techniques may vary according to the preference of the implanting physician and the needs of the patient. Experience indicates that there are several satisfactory suturing methods.
2. Size the annulus again to confirm valve size selection.
3. Use the valve holder handle to align the valve in the annulus. **Orient the valve so that blood flow is always into the pivot guards.** Align the valve such that the pivot guards are in the desired orientation in the annulus.

Figure 15. Blood flow into pivot guards, mitral valve inflow mount

1. Flow
2. Pivot guards
3. Leaflets


Figure 16. Blood flow into pivot guards, aortic valve outflow mount

NOTE: Do not use instruments other than the valve holder handle to align and seat the valve, and exercise extreme care to avoid putting stress on the valve orifice or leaflets.

4. Place sutures in the outer half of the valve sewing cuff. Incorporate adequate sewing cuff material into each suture pass to maintain the valve in the desired position. It may be helpful to tag the sutures adjacent to pivot guards when inserting them. To prevent cutting the fibers in the sewing cuff, use of standard round or taper point needles is recommended.
5. Parachute the valve into place.
6. Cut the two retention sutures from the holder/rotator, and gently withdraw the holder/rotator from the valve. Remove the valve retention sutures. Maintain the holder/rotator in the sterile field for further use as a rotation instrument.

Figure 17. Cut the valve retention suture

7. Using the Model LT100 Leaflet Tester, open the valve and inspect the area for any obstructive tissue. If visualization is inadequate, use the leaflet tester to confirm free leaflet motion.

WARNING: Do not use hard or rigid instruments to test leaflet mobility, as this may result in structural damage to the valve or thromboembolic complications.

8. Tie down the pivot guard sutures first, and then tie down the remaining sutures.

WARNING: Cut suture ends short, especially in the vicinity of the pivot guards, to prevent leaflet impingement.


9. Test leaflet motion again, and if desired, rotate the valve using the holder/rotator handle (see Valve Rotation (page)).

WARNING: The two retention sutures on the valve holder/rotator must be cut and removed before the valve can be rotated.

Valve Rotation

Using the valve holder/rotator and a valve holder handle, rotate the valve *in situ* to the desired position. The valve should rotate freely. If resistance is noted, the valve holder/rotator may not be properly seated in the valve, or the valve may be oversized. If the valve does not rotate freely, do not force valve rotation.

Figure 18. Use the Masters Series™ valve holder/rotator and the valve holder handle to rotate the valve

WARNING: Use only valve holder/rotators from Abbott to perform valve rotation. Use of other instruments could result in structural damage. The valve holder/rotator is intended for single use only and should be discarded after surgery.

CAUTION: To avoid structural damage, the mitral valve must be rotated in the fully open position.

CAUTION: To avoid structural damage, the aortic valve must be rotated in the fully closed position.

CAUTION: To minimize rotational torque, verify that the valve holder/rotator is properly seated in the valve, and that the valve holder handle is perpendicular to the valve.

Figure 19. Improperly seated valve holder/rotator



Figure 20. Properly seated valve holder/rotator

Postoperative Considerations

Echocardiography is recommended to assess valvular competency and performance. Fluoroscopy is particularly useful for determining leaflet motion in bileaflet valves.

WARNING: Do not pass catheters or other instruments through Masters Series™ Mechanical Heart Valves. This could result in scratched or damaged valve components, or leaflet fracture or dislodgment.

Patients who undergo dental or other procedures that are potentially bacteremic should receive endocarditis prophylactic antibiotic.

Anticoagulation/Antiplatelet Therapy

As there are insufficient data to indicate otherwise, Abbott recommends that patients implanted with the Masters Series™ Mechanical Heart Valve be routinely maintained on anticoagulants unless, for other reasons, it is not medically indicated.

Patient Identification Card

A Patient Identification Card is present within the packaging of the valve. Hospital personnel should complete the front side of the Patient Identification Card by printing the patient's name, implant date, along with the physician's name and healthcare facility information and phone number in the designated areas. A peel-away device label is provided in the sterile packaging which includes device information (model number, serial number, UDI number, and device name). Apply the peel-away device label to the back of the card in the designated area. If a peel-away label is not available, write the information in the designated area on the back of the card. In procedures where multiple implantable devices are implanted, include information for all implants in the designated areas.

NOTE: Write the valve serial number in the SN/LOT field.

After the card is completed, separate the top portion from the Patient Identification Card by tearing along the perforated line. Provide the completed Patient Identification Card to the patient. Instruct the patient to carry the Patient Identification Card at all times and to present the implant card to hospital staff in emergency situations. To obtain a replacement card if a patient loses or damages their card, contact Abbott Medical Technical Support. Once the replacement card is received, handwrite the required information in the respective areas and provide to the patient.

Product Materials

The following product materials are intended to come into contact with tissue.

Masters Series™ Valve with Polyester Cuff Components and Material Breakdown

Component	Material Description	Substance	Substance Concentration within each component (% w/w)(maximum)	Max Weight of component (≤) in grams
Leaflets and Orifice	Graphite Substrate with Pyrolytic Carbon Coating	Graphite Substrate with Pyrolytic Carbon Coating	100	3.1
Cuff	Polyester Fabric	Polyester Blend	100	
Sutures	PTFE-impregnated braided polyester	PTFE	3-7	1.1
		PET	93-97	
Ring Retainers	Cobalt Chromium Alloy	Nickel	35	5.0
		Cobalt	35	
		Chromium	20	
		Molybdenum	10	
Spring	Cobalt Chromium Alloy	Cobalt	40	0.6
		Chromium	20	
		Nickel	15	
		Molybdenum	7	
		Manganese	2	
		Silicon	1.2 MAX	
		Carbon	0.15 MAX	
		Beryllium	0.10 MAX	
		Phosphorous	0.015 MAX	
		Iron	Balance	

CAUTION: One or more components of the Masters Series™ Mechanical Heart Valve contains the following substance defined as CMR 1B in a concentration above 0.1% weight by weight:

Cobalt; Chemical Abstracts Service CAS No. 7440-48-4; EC No. 231-158-0

Current scientific evidence supports that medical devices manufactured from metal alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

Masters Series™ Valve with PTFE Cuff Components and Material Breakdown

Component	Material Description	Substance	Substance Concentration within each component (% w/w)(maximum)	Max Weight of component (≤) in grams
Leaflets and Orifice	Graphite Substrate with Pyrolytic Carbon Coating	Graphite Substrate with Pyrolytic Carbon Coating	100	3.1
Cuff	PTFE Fabric	PTFE	100	
Sutures	PTFE-impregnated braided polyester	PTFE	3-7	2.3
		PET	93-97	
Ring Retainers	Cobalt Chromium Alloy	Nickel	35	0.4
		Cobalt	35	
		Chromium	20	
		Molybdenum	10	
Spring	Cobalt Chromium Alloy	Cobalt	40	0.6
		Chromium	20	
		Nickel	15	
		Molybdenum	7	
		Manganese	2	
		Silicon	1.2 MAX	
		Carbon	0.15 MAX	
		Beryllium	0.10 MAX	
		Phosphorous	0.015 MAX	
		Iron	Balance	

CAUTION: One or more components of the Masters Series™ Mechanical Heart Valve contains the following substance defined as CMR 1B in a concentration above 0.1% weight by weight:

Cobalt; Chemical Abstracts Service CAS No. 7440-48-4; EC No. 231-158-0

Current scientific evidence supports that medical devices manufactured from metal alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

Incident Reporting

If, in the course of use of this device, you have reason to believe that a serious incident occurred, please report it to the manufacturer. For customers in the European Union, report the serious incident to your national authority as well as to the manufacturer.

Disposal

This instructions for use is recyclable. Dispose of all packaging materials as appropriate. Dispose of valves and accessories per standard solid biohazard waste procedures.

Disclaimer of Warranty and Limitation of Remedy

Abbott Medical LLC (Abbott) has exercised reasonable care in the manufacture of this device; however, numerous factors beyond Abbott's control can affect this device and the results obtained from its use after the device leaves Abbott's possession, including the conditions under which the device is used, diagnosis of the patient, methods of administration or handling. ABBOTT, THEREFORE, MAKES NO WARRANTIES REGARDING THIS DEVICE AND EXPRESSLY DISCLAIMS ALL WARRANTIES, WHETHER EXPRESS OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. ABBOTT SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM THE USE OF THIS DEVICE. NO PERSON HAS ANY AUTHORITY TO ASSUME ANY LIABILITY ON BEHALF OF ABBOTT WITH RESPECT TO THE DEVICE OR TO MAKE ANY GUARANTEE, REPRESENTATION, OR WARRANTY WITH RESPECT TO THE DEVICE. The exclusions and limitations set forth above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected. Descriptions of specifications in any Abbott printed material are intended solely to generally describe the device at the time of manufacture and do not constitute any express warranties. ABBOTT WILL NOT BE RESPONSIBLE FOR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE RESULTING FROM REUSE, REPROCESSING, OR RE-STERILIZATION OF THIS DEVICE.

Symbols

The symbols below and harmonized symbols may be found on the product or product label. For harmonized symbols, refer to the Universal Symbols Glossary at medical.abbott/manuals.

Symbol	Description
Poly	Standard Polyester Cuff
POLY HP	Standard Polyester HP Cuff
PTFE	Standard PTFE Cuff
EXP POLY	Expanded Polyester Cuff
EXP POLY HP	Expanded Polyester HP Cuff
MITRAL	Mitral

Symbol	Description
AORTIC	Aortic
	Rotatable
 medical.abbott/manuals	Follow instructions for use on this website
	Date
	Healthcare center or physician
	Patient identification
	Physician telephone
 UA.TR.101	Ukrainian mark of conformity to the technical regulations, where 101 is the assessment body.
MD	Medical Device
STERILE 	Contents sterile: steam
CE 2797	European conformity, affixed in accordance with the relevant provisions of European Council Regulation 2017/745 (NB 2797). Hereby, Abbott Medical declares that this device is in compliance with the relevant provisions of this regulation.
 medical.abbott/manuals	Patient Information Website
UK CA 0086	United Kingdom Conformity Assessed
UKRP	Responsible person in the United Kingdom
Mechanical Heart Valve	Mechanical heart valve
	Importer
UDI	Unique Device Identifier
ID-EFFECTIVE	Effective Orifice Diameter
ID-INFLOW	Inflow Orifice Diameter
OD-CUFF	Cuff Outer Diameter
 (7440-48-4)	Contains hazardous substances (Cobalt CAS 7440-48-4)
Double-ended handle, sizer rings	Double-ended handle, sizer rings

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St. Jude Medical Puerto Rico LLC
Abbott Medical UK Limited