Peel-Away Introducer

405104, 405108, 405112, 405116, 405118, 405119, 405120,
405122, 405124, 405128, 405129, 405136, 405144, 405145,
405146, 405147, 405149, 405153, 405154, 405254, 405269,
405270, 405404, 405408, 405412, 405416, 405418, 405420,
405422, 405424, 405428

Instructions for Use

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

Pat. <http://www.abbott.com/patents>

© 2023 Abbott. All Rights Reserved.


Peel-Away Introducer

Instructions for Use

Description

The Peel-Away Introducers provide cardiovascular leads and catheters access to the venous vascular system for cardiac surgery procedures. Peel-Away Introducer sheaths have a straight shaft. The sheath is the outer component and the dilator fits inside the sheath with the distal end protruding out. A needle, which is connected to the syringe, is used to gain access to the vessel. After blood flow is established, the syringe is removed and a guidewire is inserted through the needle and advanced into the vessel. The assembled sheath/dilator combination is advanced over the guidewire into the vessel. At the proximal end of the sheath is a t-handle hub, which is designed to split and facilitate the peeling action of the scored sheath tubing. The sheath tubing is designed to peel continuously and symmetrically for safe removal of the introducer while maintaining venous access. Refer to the package label for device specifications such as size, length, and diameter.

Figure 1. Peel-Away Introducer

- 1. Sheath
- 2. Di-Lock Dilator
- 3. Guidewire
- 4. Needle
- 5. Syringe

Table 1. Peel-Away Introducer Specifications

Device Model Number	Sheath Size (French)	Dilator Size (French)	Sheath Usable Length (cm)	Dilator Usable Length (cm)
Sheath, Di-Lock Dilator				
405404	6	6	14	19
405408	7	7	14	19
405412	8	8	14	19
405416	9	9	14	19
405418	9.5	9.5	14	19
405420	10	10	14	19
405422	10.5	10.5	14	19
405424	11	11	14	19
405428	12	12	14	19

Table 1. Peel-Away Introducer Specifications

Device Model Number	Sheath Size (French)	Dilator Size (French)	Sheath Usable Length (cm)	Dilator Usable Length (cm)
Sheath, Di-Lock Dilator, 80 cm Guidewire				
405269	7	7	23	28
405270	8	8	23	28
405254	9	9	23	28
Sheath, Di-Lock Dilator, 50 cm Guidewire, Needle, Syringe				
405104	6	6	14	19
405108	7	7	14	19
405112	8	8	14	19
405116	9	9	14	19
405118	9.5	9.5	14	19
405119	8.5	8.5	14	19
405120	10	10	14	19
405122	10.5	10.5	14	19
405124	11	11	14	19
405128	12	12	14	19
405129	8	8	14	19
405136	14	14	14	19
405144	16	16	14	19
405145	8	8	14	19
405146	8.5	8.5	14	19
405147	9	9	14	19
405149	10	10	14	19
405153	7	7	14	19
405154	7	7	14	20

Intended Purpose

The Peel-Away Introducers are intended to provide a transvenous conduit for the introduction of cardiac leads and catheters into the venous vascular system.

Indications for Use

The Peel-Away Introducers are used during cardiac surgery procedures in patients indicated for cardiac rhythm management devices for the treatment of heart failure, arrhythmias, or cardiac conduction abnormalities.

Patient Target Group

The Peel-Away Introducers are intended for adult patients who require cardiac device procedures involving the introduction of cardiac leads

and catheters into the venous vasculature. The Peel-Away Introducer devices have not been studied in pediatric patients.

Intended User

The Peel-Away Introducers should only be used by physicians trained in venous access procedures.

Clinical Benefits

The Peel-Away Introducers provide a transvenous conduit to physicians accessing the venous vascular system for the introduction of cardiac leads and catheters during cardiac surgery procedures. The adult patients in these procedures are indicated for cardiac rhythm management devices for the treatment of heart failure, arrhythmias, or cardiac conduction abnormalities. The peeling function of the introducer allows for safe removal of the introducer while maintaining venous access with the cardiac lead and catheter.

Contraindications

There are no known contraindications.

Warnings

This device is not compatible with magnetic resonance imaging (MRI) systems.

Precautions

- The French size specified represents the inner diameter of the introducer sheath. Do not attempt to insert a device having a distal tip or body size larger than the introducer size indicated.
- Do not alter this device.
- Federal law (U.S.) restricts this device to sale by or on the order of a physician.
- Do not use forceps to split the t-handle or peel the sheath, it may cause the device to withdraw from the body.
- The Peel-Away Introducer sheath should withdraw freely from the body. If resistance is encountered, do not force the peeling action of the sheath or its withdrawal. Make procedural corrections, then continue.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- This device should only be used by physicians trained in venous access procedures.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into an artery may cause excessive bleeding and other complications.
- Damage to the guidewire may result if withdrawn through the needle.
- The need for symmetrical peeling of the sheath is crucial. If it is difficult to withdraw the device when the sheath is peeled, withdraw the sheath completely from the introduction site while holding the device in position. Then peel the sheath symmetrically from the device by grasping both t-handle halves.

- This device is intended for one time use only. Do not reprocess or reuse the device. Reuse can cause patient injury and the communication of infectious disease(s) from one patient to another.
- Do not allow the guidewire to inadvertently advance totally into the patient.

Potential Adverse Events

- Pneumothorax
- Hematoma formation
- Brachial plexus injury
- Sepsis/infection
- Excessive bleeding
- Vessel damage
- Air embolus
- Thrombus

Directions

CAUTION: This device should only be used by physicians trained in venous access.

1. Open the package and place the contents in a sterile field.
2. Prepare the skin and drape in the area of anticipated venous puncture.
3. Flush the needle, sheath, and dilator before use.
4. Access the vessel with an 18 gauge thin wall needle in standard fashion and insert a 0.038 in. (0.9652 mm) diameter guidewire into the vessel.

NOTE: The syringe may be attached to the needle for added support during insertion and as an aid in observing venous return.

NOTE: Leave an appropriate amount of guidewire exposed. The exposed portion of guidewire must be longer than the length of the percutaneous dilator sheath assembly.

5. Advance the dilator sheath introducer assembly over the exposed guidewire into the vessel using a twisting motion.

CAUTION: Do not allow the guidewire to inadvertently advance totally into the patient.

6. Withdraw the vessel dilator and guidewire, leaving the introducer sheath in position. Place your thumb over the exposed orifice of the sheath to reduce air aspiration and blood loss.
7. Insert and advance the device.
8. With the device advanced, split the t-handle of the introducer sheath and peel the sheath apart while withdrawing it from the vessel. Be careful not to withdraw the device during this process.

Packaging and Shelf-Life

It is recommended that the product remain in the unopened package until time of use. Contents are sterile if the package is unopened. Do not resterilize. The expiration date is marked on the outside of the package. The product must be stored in a cool, dark, and dry location. The instructions for use is recyclable. Dispose of used product and packaging following standard solid biohazard waste procedures.

Incident Reporting

If, in the course of use of this device, you have reason to believe that a serious incident occurred, please report it to the manufacturer. For customers in the European Union, report the serious incident to your national authority as well as to the manufacturer.

Clinical Data Summary

A summary of the safety and clinical performance for this device is available at <https://ec.europa.eu/tools/eudamed>. Search for the device using the UDI-DI provided on the outer packaging of the device.


Warranty

Abbott Medical warrants that its products shall be free from defects in materials and workmanship under normal use. This warranty does not exceed the “Expiration” date stated on any product labeling. The authorized uses and approved methods of use of each of our products is set forth in the related “Instructions for Use” that accompany each product. Abbott Medical disclaims any responsibility and liability for the use of its products in a manner that has not been authorized or approved. Abbott Medical’s liability under this warranty is limited to replacing its products. The foregoing warranty excludes and is in lieu of all other warranties whether expressed or implied including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Abbott Medical disclaims any liability for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this product, other than as expressly provided by specific law. Abbott Medical neither assumes nor authorizes any other person to assume for it any other or additional liability for loss, damage, or expense in connection with this product. For more details please review complete Abbott Medical warranty policy available from Abbott Medical or on the back of an Abbott Medical invoice.

Symbols

The symbols below and harmonized symbols may be found on the product or product label. For harmonized symbols, refer to the Universal Symbols Glossary at medical.abbott/manuals.

Symbol	Definition
	Medical Device
	Unique Device Identifier
	Single sterile barrier system
	Importer
	Conformité Européenne (European Conformity). Affixed in accordance with European Council Regulation 2017/745. Hereby, Abbott Medical declares that this device is in compliance with the essential requirements and other relevant provisions of this regulation.

Symbol	Definition
 UA.TR.101	Ukrainian mark of conformity to the technical regulations, where 101 is the assessment body.
 medical.abott/manuals	Follow instructions for use on this website
 (7440-48-4)	One or more components of this device may contain the following substance defined as CMR 1B in a concentration above 0.1% weight by weight: Cobalt: Chemical Abstracts Service (CAS) No. 7440-48-4; EC No. 231-158-0. Current scientific evidence supports that medical devices manufactured from cobalt alloys and stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.
	Straight Dilator
	Guidewire
	Maximum outer diameter of recommended guidewire

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