CPS Direct™ Universal

Slittable Outer Guide Catheter

Instructions for Use

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

Pat. <http://www.abbott.com/patents>

© 2022 Abbott. All Rights Reserved.

CPS Direct™ Universal

Slittable Outer Guide Catheter

Instructions for Use

Description

The CPS Direct™ Universal slittable outer guide catheter facilitates left heart lead delivery during cardiac resynchronization therapy (CRT) procedures. The CPS Direct™ Universal slittable outer guide catheter provides access to the coronary venous system and acts as a conduit for contrast medium, implantable coronary leads, or other devices. The shaft of the catheter has a silicon-based coating.

The system includes:

- CPS Direct™ Universal slittable outer guide catheter with hemostasis valve, for use in cannulating the coronary sinus vein.
- Dilator, for dilating the vein and facilitating access to the vasculature.
- CPS Direct™ Valve Bypass Tool, for facilitating the entry of devices through the hemostasis valve.

Table 1. CPS Direct™ Universal Slittable Outer Guide Catheter Models

Model	Description	Length	Inner Diameter	Outer Diameter
DS2C018	Straight	47 cm	8.01 F (2.67 mm)	9.90 F (3.30 mm)
DS2C019	115°	47 cm	8.01 F (2.67 mm)	9.90 F (3.30 mm)
DS2C020	135°	47 cm	8.01 F (2.67 mm)	9.90 F (3.30 mm)
DS2C021	Wide	47 cm	8.01 F (2.67 mm)	9.90 F (3.30 mm)
DS2C022	X-Wide	47 cm	8.01 F (2.67 mm)	9.90 F (3.30 mm)
DS2C023	Right-Sided	47 cm	8.01 F (2.67 mm)	9.90 F (3.30 mm)
DS2C024 ¹	115°-S	47 cm	8.01 F (2.67 mm)	9.90 F (3.30 mm)
DS2C025	Straight	54 cm	8.01 F (2.67 mm)	9.90 F (3.30 mm)
DS2C026	115°	54 cm	8.01 F (2.67 mm)	9.90 F (3.30 mm)
DS2C027	135°	54 cm	8.01 F (2.67 mm)	9.90 F (3.30 mm)
DS2C028	Wide	54 cm	8.01 F (2.67 mm)	9.90 F (3.30 mm)
DS2C029	3D	42 cm	8.01 F (2.67 mm)	9.90 F (3.30 mm)
DS2C030 ¹	115°-S	54 cm	8.01 F (2.67 mm)	9.90 F (3.30 mm)

¹ Model not CE marked

Intended Purpose

The CPS Direct™ Universal slittable outer guide catheter is intended to provide a transvenous conduit for navigation of the coronary venous system anatomy and delivery pathway during cardiac surgery procedures for contrast media, implantable left heart pacing leads, or other delivery tools.

Indications for Use

The CPS Direct™ Universal slittable outer guide catheter is a delivery tool used during cardiac implant procedures in patients indicated for cardiac resynchronization therapy (CRT) for treatment of heart failure.

Table 2. Accessories and Intended Purpose

Accessory (Model Number)	Intended Purpose
CPS™ Universal Slitter (DS2A003)	Slit the catheter after delivery of the left ventricular pacing lead.
CPS Direct™ Valve Bypass Tool (410194)	Enables the entry of devices through the hemostasis valve.

Contraindications

There are no known contraindications for CPS tools and accessories.

Patient Target Group

CPS tools and accessories are used in adult patients with heart failure, who are undergoing implantation of a cardiac resynchronization therapy system.

Intended User

This device should be used by or under the supervision of physicians trained in the technique of percutaneous venous access.

Warnings

- This is a single-use disposable medical device. Do not resterilize or reuse. Resterilizing or reuse of this device may negatively impact material integrity and structural stability as well as increase risk of infection.
- The contents of the package are sterile if package is unopened and undamaged. If the sterile packaging is damaged or unintentionally opened, do not use the catheter.
- The guide catheter dilator must be removed or retracted within the guide catheter tip to minimize the risk of perforation when in the right chamber of the heart.
- Temporary pacemaker should be readily available during implant as necessary.
- Do not use if stored outside the protective outer carton. Store in a cool, dark, and dry place.
- Use the catheter prior to the Use Before date on the package label.

Precautions

- Do not alter the device.
- Failure to follow the instructions may compromise catheter performance and result in complications.
- Refer to the instructions supplied with any interventional devices to be used in conjunction with the catheter for their intended uses, contraindications, and potential complications.
- Inspect all components prior to use.
- Always inspect the catheter for damage prior to insertion or re-insertion. Do not use if package is damaged. Examine all parts to ensure that its size and shape are suitable for the specific procedure for which it is to be used.
- Due to the relative stiffness of the catheter, damage to the walls of the vessels may include perforations or dissections. If resistance is encountered, do not force the catheter. Withdraw the catheter, correct the difficulty, and reinsert.
- This catheter should be maneuvered under fluoroscopic guidance. The long-term risk of protracted fluoroscopy (i.e. fluoroscopic radiation) has not been established. Therefore, careful consideration must be given for the use of the device.
- Precautions to prevent clotting should be taken whenever entering the vascular system. The patient's clotting system should be evaluated prior to catheterization.
- Use anticoagulant as training and experience dictate.
- Flush the lumen, sideport, and stopcock with normal saline prior to use in the patient. Once the catheter has been introduced into the patient, maintain the patency of the lumen by periodic aspiration and flushing of the lumen or by using a continuous infusion of normal saline through the sideport.
- Upon withdrawal of the guidewire or other inserted devices, aspiration followed by flushing through the stopcock is recommended to remove any fibrin accumulation.
- Not intended for use with high-pressure injection equipment.
- The sealing force of the hemostasis valve may alter or impair the function of some guidewires.
- Avoid twisting or excessive force on guidewires, which may cause core wire to penetrate the coil spring of the guidewire and damage the catheter. Exercise caution when using movable core guidewires as this can also cause damage to the vessels.
- Do not attempt to advance or withdraw the guidewire if unusual resistance is felt. Use fluoroscopy to determine cause.
- If excessive resistance is encountered while using the splitter, do not continue splitting the outer guide catheter.
- Stretching and/or kinking of the catheter may result in damage. Handle with care.
- Dispose of this device according to local laws and regulations.
- The device should only be used by experienced physicians trained in endovascular procedures.
- Individual patient anatomy and physician technique may require procedural variations.

- Guidewires or other inserted devices should be inserted or removed slowly through the hemostasis valve. Damage may occur due to rapid insertion or removal resulting in blood flow through the hemostasis valve. Do not advance or withdraw these devices when resistance is met. Determine cause by fluoroscopy and/or other methods and take remedial action.
- It is recommended that a guidewire be used to advance the guiding catheter into the venous system, right atrium or coronary sinus.
- Do not attempt to insert a catheter or other device with an outer diameter larger than the size specified into the catheter, hemostasis valve, or the valve bypass tool. Refer to the Catheter Models table and package labeling for size and diameter information.
- Proper precautions should be taken to prevent bleeding and air embolism. Infusion through the side port should be done only after all air is removed from the unit.
- This device should only be used with equipment that complies with international safety standards.
- Do not expose to solvents as this can compromise essential material and design characteristics leading to device failure.
- This catheter should be maneuvered under fluoroscopic guidance. The long-term risk of protracted fluoroscopy (i.e. fluoroscopic radiation) has not been established. Therefore, careful consideration must be given for the use of the device.

Potential Adverse Events

This device should be used by physicians familiar with percutaneous catheter introduction. Complications related to the use of the device and/or any catheterization procedure include, but are not limited to, the following:

Table 3. Potential Adverse Events

Adverse Event	Categories
Arrhythmia	▪ Induced atrial ectopy or arrhythmia ▪ Induced ventricular ectopy or arrhythmia ▪ Induced AV block ▪ Lead dysfunction (sensing/threshold Issue)
Cardiac perforation	▪ Cardiac tamponade ▪ Pericardial effusion ▪ Perforation of the myocardium ▪ Pericarditis
Death	
Dissection	▪ Coronary sinus dissection ▪ Coronary vein dissection
Embolism	▪ Air embolus ▪ Intravascular foreign body ▪ Dislodgement of intracardiac thrombus
Extra cardiac stimulation	▪ Diaphragmatic stimulation ▪ Intermittent phrenic nerve stimulation
Hypersensitivity	▪ Allergic reaction
Infection	▪ Endocarditis ▪ Pericarditis
Pleural perforation	▪ Hemothorax ▪ Pneumothorax
Prolonged surgery	▪ Lead dislodgement ▪ Failed or abandoned procedure
Tricuspid valve perforation	
Vascular perforation	▪ Arteriovenous fistula ▪ Arterial perforation ▪ Coronary sinus or coronary vein perforation ▪ Hematoma ▪ Venous perforation
Vascular thrombosis	▪ Thrombophlebitis ▪ Venous occlusion

How Supplied

The CPS Direct™ Universal slittable outer guide catheter is supplied sterile, for single use only, provided the package is unopened and undamaged.

Storage

Handle with care. The CPS Direct™ Universal slittable outer guide catheter should be stored in a well-ventilated area under conditions that protect it from extreme temperatures and humidity.

Store the catheter at 20°C (68°F). During storage, temperature excursions are permitted between 15°C and 25°C (59°F and 77°F). During transportation, the catheter can be exposed to temperature excursions between -5°C and 50°C (23°F and 122°F) that last up to 24 hours.

In addition, cartons containing these items should be protected from water and should not be crushed. Do not remove from outer box. Rotate stock regularly.

Instructions for Use

Inserting the CPS Direct™ Slittable Outer Guide Catheter

1. Review the label carefully for product identification and outer guide catheter configuration.
2. Flush all components of the outer guide catheter with normal saline solution.
3. Follow normal accepted practice for vessel puncture and guidewire insertion.
4. Once venous access is achieved, insert a guidewire through the needle and withdraw the needle.
5. Advance the dilator and the outer guide catheter over the guidewire and into the vein.
6. Advance the guidewire into the right atrium.
7. While retaining the guidewire in place, remove the dilator from the outer guide catheter.
8. Advance the outer guide catheter over the guidewire and into the right atrium.
9. Retract the guidewire, leaving the outer guide catheter in the right atrium.
10. Position the outer guide catheter into the coronary sinus using the preferred method. An inner catheter may be used to facilitate placement.
11. Insert and advance a catheter or lead into the desired position. A valve bypass tool, such as the CPS Direct™ Valve Bypass Tool, may be used to enhance introduction of a catheter or lead through the hemostasis valve as follows:
 - a. Insert the valve bypass tool shaft completely into the hemostasis valve until it rests on the cap of the hub.
 - b. Insert a catheter or lead through the lumen of the valve bypass tool and well into the sheath lumen.
 - c. While securing the catheter or lead, remove the valve bypass tool from the hemostasis valve by sliding it proximally up and off the catheter/lead.
 - d. Immediately aspirate for the presence of air and flush the outer guide catheter through the sideport.

Removing the CPS Direct™ Universal Slittable Outer Guide Catheter

NOTE:

- When removing catheters, take care to avoid dislodging the lead.
 - If an inner catheter is used, remove the inner catheter first, then the outer catheter.
1. Attach the distal portion of the CPS™ Universal Slitter to the lead by securing the lead in the channel on the top of the slitter. Hold the slitter and lead securely in one hand. Grasp the catheter to be slit with the other hand.
 2. Pull the catheter to be slit against the slitter blade and away from the slitter until the catheter is completely cut.

NOTE:

- If an inner catheter is used, hold the outer catheter steady while pulling back on the inner catheter. Repeat steps 1 and 2 to remove the outer catheter.
 - To decrease the risk of dislodging the lead, do not push the slit down the catheter. The slit must remain in a stable, fixed position while the catheter is pulled back against the blade.
3. Once the lead is visible and the catheters are removed, carefully remove the lead from the slit.
 4. Using fluoroscopy, verify the lead position.

Performance Characteristics

The CPS Direct™ Universal family of outer guide catheters is designed for transient use only during surgical procedures to facilitate left heart lead delivery. They are compatible with other products in the CPS Cardiac Positioning System family, which is an inter-compatible system of tools designed to efficiently and predictably deliver the left heart lead. The CPS catheter has not been evaluated for MR compatibility and its use in MR environment is not determined.

Incident Reporting

If, in the course of use of this device, you have reason to believe that a serious incident occurred, please report it to the manufacturer. For customers in the European Union, report the serious incident to your national authority as well as to the manufacturer.

Clinical Data Summary

A summary of the safety and clinical performance for this device is available at <https://ec.europa.eu/tools/eudamed>. Search for the device using the UDI-DI provided on the outer packaging of the device.

Clinical Benefits

The CPS Direct™ Universal slittable outer guide catheter enables the implanting physician to cannulate the coronary sinus and provides a conduit for contrast medium, implantable leads, and other delivery tools. It is an essential element in the success of implanting a left ventricular lead for successful resynchronization therapy.

Symbol Definitions

The symbols below and harmonized symbols may be found on the product or product label. For harmonized symbols, refer to the Universal Symbols Glossary at medical.abbott/manuals.

Symbol	Definition
	Caution
 medical.abbott/manuals	Follow instructions for use on this website
	CAUTION: Federal law (U.S.) restricts this device to sale by or on the order of a physician.
	Affixed in accordance with European Union Medical Device Regulation 2017/745. Hereby, Abbott Medical declares that this device is in compliance with the essential requirements and other relevant provisions of this Regulation.
	Temperature limit NOTE: The temperature range indicates excursions during transportation.
	Storage temperature
	Keep away from sunlight

Symbol	Definition
	Keep dry
	Non-pyrogenic
	Outer diameter
	Inner diameter
	Length
	Catalog number
	Do not re-use
	Do not re-sterilize
	Sterilized using ethylene oxide
	Use-by date
	Date of manufacture
	Slittable Outer Guide Catheter
	Do not use if package is damaged
	Packaging unit
	Medical Device
	Batch code
	Unique device identification number
	Single sterile barrier system
	Manufacturing facility
	Manufacturer
	Ukrainian mark of conformity to the technical regulations, where 101 is the assessment body.
	Importer

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