Amplatzer™ Septal Occluder

Instructions for Use

AR: سداة الحاجز تعليمات الاستخدام	HR: Septalni okluder Upute za uporabu	RO: Dispozitiv de închidere a septului Instrucțiuni de utilizare
BG: Оклудер за затваряне на септални дефекти Инструкции за употреба	HU: Szeptumtapasz Használati utasítás	RU: Окклюдер для закрытия дефектов межпредсердной перегородки Инструкция по эксплуатации
BS: Septalni okluder Upute za korištenje	IT: Occlusore settale Istruzioni per l'uso	SK: Septálny okluzor Návod na použitie
CS: Septální okludér Návod k použití	KO: 중격 폐색기 사용 설명서	SL: Septalno zapiralo Navodila za uporabo
DA: Septal okklusjonsanordning Brugsanvisning	LT: Pertvaros užvėrimo prietaisas Naudojimo instrukcija	SQ: Okluduesi septal Udhëzime për përdorimin
DE: Septum-okkluder Bedienungsanleitung	LV: Starpsienas nosprostotājs Lietošanas instrukcija	SR: Septalni okluder Uputstvo za upotrebu
EL: Συσκευή απόφραξης διαφράγματος Οδηγίες χρήσης	MK: Септален оклудер Упатство за употреба	SV: Okklusjonsanordning för septum Bruksanvisning
ES: Ocluser septal Instrucciones de uso	NL: Septale occluder Gebruiksaanwijzing	TR: Septal Oklüder Kullanım Talimatları
ET: Vaheseina sulguseade Kasutusjuhend	NO: septal lukkeanordning Bruksanvisning	UK: Септальний оклюдер Інструкція з використання
FI: Välliseinäaukon sulkulaite Käyttöohjeet	PL: Okluder do zamykania ubytków przegrody Instrukcja użycia	VI: Dù đóng Lỗ thông liên nhĩ Hướng dẫn Sử dụng
FR: Dispositif d'occlusion septal Mode d'emploi	PT: Dispositivo de Oclusão Septal Instruções de utilização	ZH: 间隔封堵器 使用说明

STERILE EO

CE
2797

Abbott

Symbols and Definitions / رموز وتعريفات / Символи и характеристики / Simboli i definicije / Symboly a definice / Symboler og definitioner / Symbole und Definitionen / Σύμβολα και ορισμοί / Símbolos y definiciones / Sümbold ja määratlused / Merkkien selitykset / Symboles et définitions / Simboli i definicije / A szimbólumok magyarázata / Simboli e definizioni / 기호 및 정의 / Simboliai ir apibrėžimai / Simboli un definīcijas / Символи и дефиниции / Symbolen en definities / Symboler og definisjoner / Symbole i definicje / Símbolos e definições / Simboluri și definiții / Условные обозначения и определения / Symboly a definície / Simboli in opredelitve / Simbole dhe përkufizimet / Simboli i definicije / Symboler och definitioner / Semboller ve Tanımları / Символи та визначення / Ký hiệu và Định nghĩa / 符号及其定义

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Amplatzer™ Septal Occluder

EN: Instructions for Use

CAUTION: Read all instructions carefully. Failure to follow these instructions, warnings, and precautions may lead to patient injury or occluder damage.

Device Description

The Amplatzer™ Septal Occluder is a self-expanding occluder made of braided nitinol wire designed for use in occluding atrial septal defects (ASD) and fenestrations following a Fontan procedure. The occluder has a short waist, a right atrial disc, and a left atrial disc. The diameter of the waist corresponds to the atrial defect size and is designed to fill the ASD. The occluder is filled with polyester patches that are securely sewn to the right atrial disc, the left atrial disc, and the waist by a polyester thread. The occluder consists of radiopaque marker bands at each end of the device, one attached to the left atrial disc and one connected to the end screw on the right atrial disc. The Amplatzer™ Septal Occluder is sterile and for single use only.

Refer to the figures and table in Appendix A: Supplemental Information for more information about the occluder. Table 2 shows the device and delivery sheath dimensions. Figure 1 and Figure 2 illustrate the following occluder and delivery system components:

Appendix A Figure 1

- | | |
|---|---------------------------|
| 1 | Right Atrial Disc |
| 2 | Left Atrial Disc |
| 3 | Waist |
| 4 | End Screw and Marker Band |
| 5 | Marker Band |

Appendix A Figure 2

- | | |
|----|---|
| 6 | Loader |
| 7 | Hemostasis Valve with Extension Tube and Stopcock |
| 8 | Sheath |
| 9 | Dilator |
| 10 | Delivery Cable |
| 11 | Plastic Vise |

Amplatzer™ Septal Occluder Performance Characteristics

The Amplatzer™ Septal Occluder is a self-expanding occluder made of braided nitinol wire. The occluder has a short waist designed to fill the ASD or fenestration and two retention discs that stabilize the occluder within the defect. Polyester patches are securely sewn within each disc and the waist, promoting occlusion. The braided nitinol wire and polyester patches promote occlusion and tissue ingrowth. Radiopaque marker bands provide visual guidance under fluoroscopy, aiding in precise occluder deployment. The end screw allows the attachment, recapture, and release of the occluder, enabling precise placement within the anatomy.

Intended Purpose

The Amplatzer™ Septal Occluder is a device intended for the nonsurgical occlusion of ostium secundum atrial septal defects or fenestrations.

Indications

The Amplatzer™ Septal Occluder is indicated for patients with ostium secundum atrial septal defect or fenestrations.

Target Patient Population

Patients with an ostium secundum atrial septal defect with clinical evidence of right ventricular (RV) volume overload (such as 1.5:1 degree of left-to-right shunt, or RV enlargement on echocardiography), or patients who have undergone a fenestrated Fontan procedure and require closure of the fenestration.

Intended User

Physicians trained in standard transcatheter techniques

Intended Clinical Benefit

The following benefits are expected to be achieved:

- Minimally invasive alternative to surgical closure associated with lower complication rates
- Reduced hospital stay compared to surgery
- Immediate atrial septal defect closure with subsequent alleviation of associated morbidities

Summary of Safety and Clinical Performance

A summary of the safety and clinical performance for this device is available at ec.europa.eu/tools/eudamed: This is the SSCP location after the launch of the European Database on Medical Devices/Eudamed. Search for the device using the Basic UDI-DI: 5415067AMP1100DV.

Limitations

No known limitations

Contraindications

The Amplatzer™ Septal Occluder is contraindicated for:

- patients known to have extensive congenital cardiac anomaly that can only be adequately repaired with cardiac surgery.
- patients with active endocarditis or other systemic infections requiring treatment.
- patients who are unable to tolerate intra-procedural anticoagulation or post-procedural anti-platelet therapy.
- patients whose vascular access site is inadequate to accommodate appropriate sheath size.
- patients known to have demonstrated intracardiac thrombus (especially left atrial or left atrial appendage thrombi).
- patients whose size or condition (for example, too small for transesophageal echocardiography [TEE] probe, catheter size, vasculature size, active infection) would cause the patient to be a poor candidate for cardiac catheterization.
- patients with defect margins less than 5 mm to the coronary sinus, inferior vena cava rim, an atrioventricular valve, or the right upper lobe pulmonary vein.

Warnings

- The Amplatzer™ Septal Occluder consists of a nickel-titanium alloy, which is generally considered safe. However, in vitro testing has demonstrated that nickel is released from this occluder for a minimum of 60 days. Patients who are allergic to nickel may have an allergic reaction to this occluder, especially those with a history of metal allergies. Certain allergic reactions can be serious; patients should be instructed to notify their physicians immediately if they suspect they are experiencing an allergic reaction such as difficulty breathing or inflammation of the face or throat. Some patients may also develop an allergy to nickel if this occluder is implanted. For more information, see Product Materials.
- DO NOT use the Amplatzer™ Septal Occluder after the Use-by date stated on the package label.
- Do not use the occluder if the sterile packaging is open or damaged.

- This occluder was sterilized using ethylene oxide and is for single use only. Do not reuse or resterilize. Attempts to resterilize the occluder can cause a malfunction, insufficient sterilization, or harm to the patient.
- Physicians must be prepared to deal with urgent situations such as pericardial effusion or occluder embolization that can require removal of the occluder.
- Remove embolized devices. Do not remove an embolized occluder through intracardiac structures unless the occluder is fully recaptured inside a sheath.
- Do not use a power injection syringe to inject contrast solution through the sheath.
- Implantation of this occluder may not replace the need for anticoagulation therapy in patients with ASD and paradoxical embolism.
- The use of transthoracic, transesophageal, or intracardiac echocardiographic imaging (TTE, TEE, or ICE) is required before, during, and after implantation. Imaging is mandatory as an aid in positioning the sheath through the most suitable defect to achieve complete closure of the defect. If TEE is used, the patient's esophageal anatomy must be adequate for the placement and manipulation of the TEE probe.
- When using a stop-flow technique with a compliant balloon catheter, do not inflate the balloon beyond the cessation of the shunt. DO NOT OVERINFLATE.
- Physicians should be aware of the risk of erosion. Erosion may be the result of a complex set of interactions between factors including, but not limited to the underlying anatomic substrate, retro-aortic rim of less than 5 mm in any echocardiographic plane, and the dynamic hemodynamic counter-motion between the atrium and the aorta. However, there is insufficient data including inadequate echocardiographic information about the cases already reported to determine etiology of erosion.
- Do not select an occluder size greater than 1.5 times the atrial septal defect diameter as determined by echocardiography, since it may increase the risk of erosion.
- Do not release the occluder from the delivery cable if the occluder does not conform to its original configuration, if the occluder position is unstable, or if the occluder interferes with any adjacent cardiac structure such as superior vena cava (SVC), pulmonary vein (PV), mitral valve (MV), coronary sinus (CS), or aorta (AO). Recapture the occluder and redeploy. If still unsatisfactory, recapture the occluder and replace with a new occluder.
- Placement of the Amplatzer™ Septal Occluder may impact future cardiac interventions (for example, transseptal puncture and mitral valve repair).

Precautions

- Use standard transcatheter techniques when using Amplatzer™ products.
- This occluder should only be used by physicians who have been trained in standard transcatheter techniques. The physician should determine which patients are candidates for procedures that use this occluder.
- The physician should exercise clinical judgment in situations that involve the use of anticoagulants or anti-platelet drugs before, during, and/or after use of this occluder.
- Prolonged procedures can cause increased exposure to anesthesia, contrast media, and/or radiation.
- Store in a dry place.

Patient Selection

Certain patients may be at higher risk for complications such as tissue erosion and occluder embolization. If higher-risk patients have occluders implanted, closer follow-up is warranted (see Post-Procedure Care). Higher-risk patients include:

- patients with deformation of the occluder at the aortic root.
- patients with high defects (minimal aortic and superior rims).
- patients with IVC rim deficiency (risk of occluder embolization).

Patients with multiple ASDs:

- Closure of multiple ASDs should only be attempted by physicians who have gained sufficient experience (more than 10–15 cases) to undertake more technically challenging procedures.
- If two large ASDs are separated by more than a 7 mm rim of tissue, implantation of two occluders may be justified.
- If multiple fenestrations are close to each other, one occluder may be used to cover all defects when placed in the largest defect.

Special Populations

- **Pregnancy:** The safety and effectiveness of this occluder has not been established during pregnancy. Fluoroscopic X-ray guidance is used during placement of the device. The risk of increased X-ray exposure for patients who are pregnant must be weighed against the potential benefits of this technique.
- **Nursing mother:** The safety and effectiveness of this occluder has not been established in lactating mothers. There has been no quantitative assessment for the presence of leachables in breast milk.

MR Conditional

Non-clinical testing has demonstrated that the Amplatzer™ Septal Occluder device is MR Conditional. A patient with the Amplatzer™ Septal Occluder can be safely scanned in an MR system under the following conditions:

- Static magnetic field of 1.5 Tesla (1.5T) or 3.0 Tesla (3.0T)
- Maximum spatial gradient field of 19 T/m (1900 G/cm)
- Maximum MR system reported, whole-body averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode).

Under the scan conditions defined above, the occluder is expected to produce a maximum temperature rise of less than or equal to 5°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the occluder extends approximately 16 mm from the Amplatzer™ Septal Occluder when imaged with a gradient echo pulse sequence in a 3.0T MR system.

Potential Adverse Events

Potential adverse events that may occur during or after a procedure placing this device include, but are not limited to:

- | | |
|--|---------------------------------------|
| • air embolus | • heart block |
| • allergic reaction/toxic effect due to contrast media, medication, or metal | • hematoma |
| • apnea | • hemolysis |
| • arrhythmia | • hypertension |
| • arteriovenous fistula | • hypotension |
| • bleeding | • infection |
| • brachial plexus injury | • myocardial infarction |
| • cardiac perforation | • perforation of vessel or myocardium |
| • cardiac tamponade | • pericardial effusion |
| • cardiovascular injury | • peripheral embolism |
| • death | • peripheral pulse loss |
| | • phrenic nerve injury |

- device embolization
- device malformation
- device migration
- dissection
- endocarditis
- erosion
- fever
- foreign material embolic event
- groin pain
- headache/migraine
- pleural effusion
- pseudoaneurysm
- residual shunt
- stroke/transient ischemic attack
- thromboembolism
- thrombus
- valve damage
- valvular insufficiency
- vascular access site injury

Compatibility

The Amplatzer™ Septal Occluder is compatible with 45-degree models of the Amplatzer™ TorqVue™ Delivery System, Amplatzer™ Trevisio™ Intravascular Delivery System, and the Amplatzer™ TorqVue™ Exchange System.

Device Lifetime

The Amplatzer™ Septal Occluder is expected to remain implanted permanently. The device lifetime for the Amplatzer™ Septal Occluder is at least 10 years, as demonstrated by bench testing of the occluder.

Patient Identification Card

A patient identification card is present within the packaging of the occluder. Hospital personnel should complete the front side of the Patient Identification Card by printing the patient's name, implant date, physician's name, healthcare facility information, and phone number in the designated areas. A peel-away device label is provided in the sterile packaging, which includes occluder information (model number, lot number, UDI number, and device name). Apply the peel-away device label to the back of the card in the designated area. If a peel-away label is not available, write the information in the designated area on the back of the card.

NOTE: Write the occluder lot number in the SN/LOT field.

After the card is completed, separate the top portion from the Patient Identification Card by tearing along the perforated line. Provide the completed Patient Identification Card to the patient. Instruct the patient to carry the Patient Identification Card at all times and to present the card to hospital staff in emergency situations. To obtain a replacement card if a patient loses or damages their card, contact Abbott Medical Technical Support. Once the replacement card is received, handwrite the required information in the respective areas and provide the new card to the patient.

Product Materials

The following materials are intended to come into contact with tissue:

Table 1. Product Composition

Component	Material Description	Substance	Maximum Substance Concentration in Component (% w/w)	Maximum Substance Mass Based on Largest Occluder (milligrams)
Braided Nitinol Wires	Nitinol	Nickel	<57	1100.5
		Titanium	<43	830.2
		Carbon	<0.1	2.0
		Cobalt	<0.1	2.0
		Copper	<0.1	2.0
		Chromium	<0.1	2.0
		Hydrogen	<0.1	2.0
		Iron	<0.1	2.0
		Niobium	<0.1	2.0
		Nitrogen	<0.1	2.0
		Oxygen	<0.1	2.0
End Screw	316L/316 LVM Stainless Steel	Iron	<69	42.2
		Chromium	<18	11.0
		Nickel	<14	8.6
		Molybdenum	<3	1.9
		Manganese	<2	1.3
		Silicon	<1	0.7
		Nitrogen	<0.1	0.1
		Phosphorus	<0.1	0.1
		Carbon	<0.1	0.1
		Sulfur	<0.1	0.1
Marker Band	Platinum / Iridium Alloy	Platinum	80	30.7
		Iridium	20	7.7
Polyester Thread	Polyester	Polyester Blend	100	18.8
Polyester Patch	Polyester	Polyethylene Terephthalate	100	93.2
CAUTION: One or more components of this occluder may contain the following substance defined as CMR 1B in a concentration above 0.1% weight by weight (w/w): Cobalt: Chemical Abstracts Service (CAS) No. 7440-48-4; EC No. 231-158-0. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.				

Directions for Use

All figures and tables referenced in this section are located in Appendix A: Supplemental Information.

Materials recommended for use with an Amplatzer™ Septal Occluder:

- 30–60 cc luer-lock syringe
- 0.035 inch guidewire
- Compliant balloon catheter (lower pressure inflation)

Pre-Procedure Care

It is recommended that at least 24 hours prior to the procedure that the patient receives aspirin (81 mg or 325 mg) or an alternative anti-platelet/anticoagulant if the patient has an aspirin intolerance.

CAUTION: When placing the occluder using the compatible delivery systems, refer to the instructions for use provided with the delivery systems.

General instructions for the Amplatzer™ Septal Occluder are provided below:

Procedure

1. Prepare the patient for a standard transcatheter procedure. Maintain a recommended activated clotting time (ACT) of greater than 200 seconds. Effective anticoagulation therapy should be maintained during and after the procedure as determined by the physician.
2. Puncture the femoral vein and perform a standard right-heart catheterization.
3. Perform an angiogram to demonstrate the atrial defect. Catheterize the left atrium using a 45° LAO position and cranial angulation of 35°–45°. Inject contrast medium into the right upper lobe pulmonary vein.
NOTE: Occluder size and placement are based on the defect location.
4. Use the J-tip guidewire to gain access to the left atrium.
5. Echo guidance: The procedure should be performed under echo guidance to allow for comprehensive assessment of all rims (with specific emphasis on adequacy of the anterior-superior rim) and cardiac structures to enable appropriate placement and position of the occluder, and to assess whether acute ASD closure has been achieved without pathologic interference or impingement on important surrounding cardiac structures. Use echocardiography to determine the diameter of the defect once a stop-flow condition has been achieved with a compliant balloon catheter. Refer to the instructions for use included with the compliant balloon catheter for more information.
6. Determine the diameter of the defect:
 - a. Insert a compliant balloon catheter over the guidewire into the left atrium.
 - b. Under echocardiographic guidance, place the balloon across the defect and inflate the balloon with diluted contrast medium until the shunt across the defect ceases.
 - c. Deflate the balloon until flow is visible, and then reinflate the balloon until the shunting ceases.

WARNING: Do not inflate the balloon beyond the stop-flow point or beyond the balloon's maximum inflation volume. Inflation beyond the stop-flow point may cause distention of the defect and/or balloon damage.

- d. Measure the defect using echocardiographic imaging.

NOTE: A waist in the balloon could appear without the cessation of flow. This would occur if there is more than one atrial septal defect. Echocardiography measurement should occur based on stop-flow, not the appearance of a waist.

7. Once the diameter of the defect is determined, select an occluder size equal to or one size greater than the diameter of the defect. Refer to Table 2 in Appendix A: Supplemental Information.

WARNING: Do not select a device size greater than 1.5 times the atrial septal defect diameter determined by echocardiography, since it may increase the risk of erosion.

8. Remove the balloon catheter and leave the 0.035 inch guidewire in place.
9. Select the appropriate delivery system for the occluder size (see Table 2 in Appendix A: Supplemental Information).
10. Prepare the delivery system for use according to its instructions for use.
11. Insert the dilator into the delivery sheath and tighten the rotating luer. Advance the dilator and delivery sheath over the guidewire, through the defect, and into the left upper pulmonary vein. Verify the correct position of the sheath by a test hand injection of contrast medium or by echocardiography.
CAUTION: If strong resistance is met during manipulation, discontinue the procedure and determine the cause of resistance before proceeding. If the cause of the resistance cannot be determined, withdraw all components.
12. Attach the hemostasis valve to the proximal end of the loader and flush with sterile saline. Slowly remove the dilator to prevent ingress of air. Allow blood backflow to purge all air from the system.
13. Slowly remove the guidewire and flush the delivery sheath with sterile saline.
14. Prepare the occluder for use:
 - Inspect the outer sterile pouch and verify that it is unopened and undamaged. Do not use the occluder if the outer pouch is open or damaged.
 - Gently open the outer, then the inner, pouches.
15. Advance the delivery cable through the hemostasis valve and loader.
16. Attach the occluder to the distal end of the delivery cable by rotating the occluder clockwise approximately five turns until the occluder is secure. Then rotate the occluder counterclockwise 1/8 of one turn to facilitate release of the occluder.
17. Immerse the occluder and loader in sterile saline and pull back on the delivery cable to retract the occluder inside the loader.
18. Attach the loader to the delivery sheath and tighten the rotating luer to lock the components together.

CAUTION: Attach the loader firmly to the sheath to ensure that there are no gaps between the inner surfaces of the two components.

19. Advance the delivery cable and occluder through the delivery sheath until the occluder reaches the tip of the delivery sheath. Do not rotate the delivery cable.
CAUTION: Do not advance the delivery cable and occluder if resistance is felt.
20. Use angiography and echocardiography for guidance. Hold the delivery cable in place while retracting the delivery sheath to deploy the left atrial disc and part of the connecting waist. Pull the occluder gently against the atrial septum. This can be felt and observed by echocardiography.
21. Maintain a slight tension on the delivery cable while retracting the delivery sheath approximately 5–10 cm to deploy the right atrial disc.
22. Use angiography and echocardiography to confirm that the occluder is in place and to evaluate for residual shunt or valve insufficiency.

WARNING: Do not release the occluder from the delivery cable if the occluder does not conform to its original configuration or if occluder position is unstable or interferes with any adjacent cardiac structure (such as SVC, PV, MV, CS, AO). Advance the delivery sheath to recapture the occluder and redeploy. If the occluder position is still unsatisfactory, recapture the occluder and replace it with a new occluder or abort the procedure.

NOTE: The occluder can be deployed a maximum of two times.

23. Use imaging to verify correct position of the occluder.

24. When the occluder placement is confirmed the occluder may be released.
- Detach the occluder by turning the plastic vise counterclockwise.

NOTE: In the unlikely event that occluder detachment is not possible, advance the delivery sheath against the right atrial disc to secure the occluder and to facilitate detachment.

- Remove the delivery cable and delivery sheath from the patient.

Post-Procedure Care

- All patients should be kept overnight for observation. A transthoracic echocardiogram should be performed before discharge.
- Patients with any observed pericardial effusion following the occluder implantation should be closely monitored with serial echocardiogram performed until the pericardial effusion resolves.
- Clinical follow-up with a cardiologist and echocardiograms are recommended at implant, one day post-implant, pre-discharge, and again at one week, one month, six months, and 12 months post-implant. Clinical follow-up with a cardiologist annually thereafter is also recommended.
- Patients should be educated to seek immediate medical attention that includes an echocardiogram, if they develop signs or symptoms of hemodynamic instability such as chest pain, arrhythmia, fainting, or shortness of breath.
- Patients should be instructed to avoid strenuous activity for a minimum of one month post-device implant or as directed by their physician. Strenuous activities may lead to the increased risk of adverse events including erosion. Patients should be reminded that if they experience any symptoms of shortness of breath or chest pain at any time, and especially after strenuous activity, they should seek medical care immediately.
- Patients should take appropriate endocarditis prophylaxis for six months following occluder implantation. The decision to continue endocarditis prophylaxis beyond six months is at the discretion of the physician.
- Patients should be treated with anti-platelet/anticoagulation therapy (such as aspirin) for six months post-implant. The decision to continue anti-platelet/ anticoagulation therapy beyond six months is at the discretion of the physician.

Storage and Handling Information

There are no special storage requirements beyond what is typical for medical devices.

Reporting Device Incidents

If, in the course of using this device, you have reason to believe that a serious incident occurred, please report it to the manufacturer. For customers in the European Union, report the serious incident to your national authority as well as to the manufacturer.

Disposal

Amplatzer™ Septal Occluder components should be appropriately classified as biohazards and disposed of in compliance with applicable facility procedures and local and country laws and regulations.

Warranty

Abbott Medical warrants to buyer that, for a period equal to the validated shelf life of the product, this product shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Abbott Medical's obligation under this warranty is limited to replacing or repairing at its option, at its factory, this product if returned within the warranty period to Abbott Medical and after confirmed to be defective by the manufacturer.

EXCEPT AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY AS TO MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

See the Terms and Conditions of Sale for further information.

Australian Warranty

This section is only provided in English.

This warranty is given by Abbott Medical. To make inquiries regarding this warranty, use the contact information for either the manufacturer or Australian Sponsor.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

The benefits to you given by this warranty are in addition to your other rights and remedies under the Australian Consumer Law.

Abbott Medical warrants to the buyer that the goods, for a period equal to the validated shelf life of the goods (the "Warranty Period"), shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Unless otherwise obligated by law or mandatory consumer guarantees, under this warranty Abbott Medical will, at its option, replace or repair the goods at its factory, if the buyer, at the buyer's expense, returns the goods within the Warranty Period to Abbott Medical at the address provided with the product and after Abbott Medical confirms that the goods are defective.

To the maximum extent permitted by law, if a mandatory term is implied by law or a mandatory consumer guarantee applies to the goods, and the goods are not of a kind ordinarily acquired for personal, domestic or household use or consumption, Abbott Medical's liability for the breach of the term or guarantee is limited to, at Abbott Medical's option, either replacement or repair of the goods or payment of the costs of replacing or repairing the goods.

EXCEPT FOR MANDATORY TERMS IMPLIED BY LAW, MANDATORY CONSUMER GUARANTEES, OR AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED. See the Terms and Conditions of Sale for further information.

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™ Indicates a trademark of the Abbott group of companies.
‡ Indicates a third party trademark, which is property of its respective owner.
© 2023 Abbott. All Rights Reserved.
Abbott Medical UK Limited