# Absorb

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English

# Absorb

## Bioresorbable Vascular Scaffold System

INFORMATION FOR PRESCRIBERS  
Table of Contents

|      |                                                                                                                               |
|------|-------------------------------------------------------------------------------------------------------------------------------|
| 1.0  | DEVICE DESCRIPTION                                                                                                            |
|      | Table 1: Drug Content in Absorb BVS                                                                                           |
|      | Table 2: <i>In vitro</i> Device Specifications                                                                                |
|      | Table 3: TagAlert Window Indicator Options, Instructions for Use and Disposal                                                 |
| 2.0  | HOW SUPPLIED                                                                                                                  |
| 3.0  | INDICATIONS                                                                                                                   |
| 4.0  | CONTRAINDICATIONS                                                                                                             |
| 5.0  | WARNINGS                                                                                                                      |
| 6.0  | PRECAUTIONS                                                                                                                   |
|      | 6.1 Scaffold Handling – Precautions                                                                                           |
|      | 6.2 Scaffold Placement – Precautions                                                                                          |
|      | 6.3 Use in Conjunction with Other Procedures                                                                                  |
|      | 6.4 Scaffold / System Removal – Precautions                                                                                   |
|      | 6.5 Post Implant – Precautions                                                                                                |
|      | 6.6 MRI Statement                                                                                                             |
|      | 6.7 Drug Interactions                                                                                                         |
|      | 6.8 Pregnancy                                                                                                                 |
| 7.0  | ADVERSE EVENTS                                                                                                                |
|      | 7.1 Observed Adverse Events                                                                                                   |
|      | 7.2 Potential Adverse Events                                                                                                  |
| 8.0  | CLINICAL INVESTIGATIONS OF ABSORB BVS SYSTEM                                                                                  |
|      | 8.1 ABSORB Clinical Trial                                                                                                     |
|      | 8.2 ABSORB Cohort A                                                                                                           |
|      | 8.2.1 Methodology and Current Status                                                                                          |
|      | 8.2.2 Clinical Outcome at 5 Years                                                                                             |
|      | Table 4: Clinical Outcome at 5 Years                                                                                          |
|      | 8.2.3 Angiographic, IVUS, and OCT Outcomes at 180 Days and 2 Years                                                            |
|      | 8.2.4 Vasomotor Function Results at 2 Years                                                                                   |
|      | Figure 1: Results of Acetylcholine and Methergine testing in proximal, scaffolded, and distal segments                        |
|      | 8.2.5 Discussion                                                                                                              |
|      | 8.3 ABSORB Cohort B                                                                                                           |
|      | 8.3.1 Methodology and Current Status                                                                                          |
|      | 8.3.2 Clinical Outcomes for Full Cohort B                                                                                     |
|      | Table 5: Hierarchical Clinical Outcomes for Cohort B (ITT Population)                                                         |
|      | 8.3.3 Angiographic, IVUS, and OCT Outcomes at 180 Days, 1, 2, and 3 Years                                                     |
|      | 8.3.4 Discussion                                                                                                              |
|      | 8.4 ABSORB EXTEND                                                                                                             |
|      | 8.4.1 Methodology and Current Status                                                                                          |
|      | 8.4.2 Summary of Clinical Results out to 1 year (Interim Data Cutoff)                                                         |
|      | Table 6: ABSORB EXTEND Subject Counts of Ischemia-driven Adverse Events through 393 Days (All Subjects Registered Population) |
|      | 8.4.3 Summary of Clinical Results out to 2 years (Interim Data Cutoff)                                                        |
|      | Table 7: ABSORB EXTEND Subject Counts of Ischemia-driven Adverse Events through 758 Days (250 Subjects Registered Population) |
|      | 8.5 Propensity Score Adjusted Analysis of Absorb BVS System and XIENCE V                                                      |
|      | 8.5.1 6-month Propensity Analysis                                                                                             |
|      | Table 8: MACE at 6 Months from the IPSW-Adjusted Data                                                                         |
|      | Table 9: Angiographic Outcomes at 6 Months from the IPSW-Adjusted Data                                                        |
|      | 8.5.2 1-year Propensity Analysis                                                                                              |
|      | Table 10: Clinical outcomes at 1 year from the IPSW-Adjusted Cohorts                                                          |
|      | 8.5.3 2-year Propensity Analysis                                                                                              |
|      | Table 11: Clinical outcomes at 2 years from the IPSW-Adjusted Cohorts                                                         |
|      | 8.6 Benefit of the Absorb Bioresorbable Scaffold Technology                                                                   |
| 9.0  | PATIENT SELECTION AND TREATMENT                                                                                               |
|      | 9.1 Individualization of Treatment                                                                                            |
| 10.0 | CLINICIAN USE INFORMATION                                                                                                     |
|      | 10.1 Inspection Prior to Use                                                                                                  |
|      | 10.2 Materials Required                                                                                                       |
|      | 10.3 Vessel and Lesion Selection                                                                                              |
|      | Table 12: Target Vessel Diameter Ranges and Absorb BVS Diameter to be Used                                                    |
|      | 10.4 Preparation                                                                                                              |
|      | 10.4.1 Packaging Removal                                                                                                      |
|      | 10.4.2 Dual Layer Sheath Removal                                                                                              |
|      | 10.4.3 Guide Wire Lumen Flush                                                                                                 |
|      | 10.4.4 Delivery System Preparation                                                                                            |
|      | 10.5 Delivery Procedure                                                                                                       |
|      | 10.6 Deployment Procedure                                                                                                     |
|      | 10.7 Further Expansion of the Deployed Scaffold                                                                               |
|      | 10.8 Removal Procedure                                                                                                        |
| 11.0 | TRADEMARKS                                                                                                                    |

1.0 DEVICE DESCRIPTION

The Absorb Bioresorbable Vascular Scaffold (BVS) System includes:

- A pre-mounted polymer poly (L-lactide) (PLLA) scaffold coated with a blend of the antiproliferative drug everolimus and polymer poly (D,L-lactide) (PDLLA) in a 1:1 ratio. The available dose of everolimus on the scaffold is shown in Table 1.

Table 1: Drug Content in Absorb BVS

| Scaffold Diameter (mm) | Scaffold Length (mm) | Drug Dose (µg) |
|------------------------|----------------------|----------------|
| 2.5, 3.0               | 8                    | 76             |
| 2.5, 3.0               | 12                   | 114            |
| 2.5, 3.0               | 18                   | 181            |
| 2.5, 3.0               | 23                   | 228            |
| 2.5, 3.0               | 28                   | 276            |
| 3.5                    | 12                   | 135            |
| 3.5                    | 18                   | 197            |

| Scaffold Diameter (mm) | Scaffold Length (mm) | Drug Dose (µg) |
|------------------------|----------------------|----------------|
| 3.5                    | 23                   | 246            |
| 3.5                    | 28                   | 308            |

- Four radiopaque markers located on the end rings of the scaffold mark the scaffold length prior to deployment and after expansion in the artery because the Absorb BVS is not visible under fluoroscopy.
- Two radiopaque markers, located underneath the balloon, fluoroscopically mark the working length of the balloon and the location of the undeployed scaffold of the scaffold delivery system.
- Absorb BVS System has a rapid exchange (RX) scaffold delivery system.
- Two proximal delivery system shaft markers (95 cm and 105 cm proximal to the distal tip) indicate the relative position of the delivery system to the end of brachial or femoral guiding catheter. Working catheter length is 143 cm.
- A shaft color change denotes the guide wire exit notch.

Table 2: *In vitro* Device Specifications

| Scaffold Diameter (mm) | Scaffold Length (mm) | * Minimum Guiding Catheter Compatibility (ID) | ** <i>In vitro</i> Scaffold Nominal Pressure |     | Rated Burst Pressure – RBP |      | Scaffold Free Area (%) |
|------------------------|----------------------|-----------------------------------------------|----------------------------------------------|-----|----------------------------|------|------------------------|
|                        |                      |                                               | atm                                          | kPa | atm                        | kPa  |                        |
| 2.5                    | 8                    | 6F (0.070"/1.8 mm)                            | 6                                            | 608 | 16                         | 1621 | 68                     |
| 2.5                    | 12                   | 6F (0.070"/1.8 mm)                            | 6                                            | 608 | 16                         | 1621 | 68                     |
| 2.5                    | 18                   | 6F (0.070"/1.8 mm)                            | 6                                            | 608 | 16                         | 1621 | 68                     |
| 2.5                    | 23                   | 6F (0.070"/1.8 mm)                            | 6                                            | 608 | 16                         | 1621 | 68                     |
| 2.5                    | 28                   | 6F (0.070"/1.8 mm)                            | 6                                            | 608 | 16                         | 1621 | 68                     |
| 3.0                    | 8                    | 6F (0.070"/1.8 mm)                            | 7                                            | 709 | 16                         | 1621 | 72                     |
| 3.0                    | 12                   | 6F (0.070"/1.8 mm)                            | 7                                            | 709 | 16                         | 1621 | 73                     |
| 3.0                    | 18                   | 6F (0.070"/1.8 mm)                            | 7                                            | 709 | 16                         | 1621 | 73                     |
| 3.0                    | 23                   | 6F (0.070"/1.8 mm)                            | 7                                            | 709 | 16                         | 1621 | 73                     |
| 3.0                    | 28                   | 6F (0.070"/1.8 mm)                            | 7                                            | 709 | 16                         | 1621 | 73                     |
| 3.5                    | 12                   | 6F (0.070"/1.8 mm)                            | 6                                            | 608 | 16                         | 1621 | 73                     |
| 3.5                    | 18                   | 6F (0.070"/1.8 mm)                            | 6                                            | 608 | 16                         | 1621 | 73                     |
| 3.5                    | 23                   | 6F (0.070"/1.8 mm)                            | 6                                            | 608 | 16                         | 1621 | 73                     |
| 3.5                    | 28                   | 6F (0.070"/1.8 mm)                            | 6                                            | 608 | 16                         | 1621 | 74                     |

\* See individual manufacturer specifications for (F) equivalent.

\*\* Assure full deployment of the scaffold (see section 10.6 **Deployment Procedure**). Deployment pressures should be based on lesion characteristics.

- TagAlert is a non-sterile temperature monitor for the shipping and storage of the Absorb BVS System. Before use of this product, check the TagAlert indicator located through the window in the back of the product box. The indicator should only show an OK as indicated in the digital display (Table 3A). If any other screen is present (Table 3B), do not use product.

Table 3: TagAlert Window Indicator Options, Instructions for Use and Disposal

|   | Window Indicator | Instructions for Use | Disposal Instructions                                                                                                   |
|---|------------------|----------------------|-------------------------------------------------------------------------------------------------------------------------|
| A |                  | Use                  | <br>TagAlert is a battery-powered device. Remove TagAlert from carton and dispose in accordance with local regulations. |
| B |                  | Do not use           |                                                                                                                         |

2.0 HOW SUPPLIED

**Sterile** – This device is E-beam radiation-sterilized. Non-pyrogenic. Do not use if the package is open or damaged.

This single-use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.

**Contents** – One (1) Absorb Bioresorbable Vascular Scaffold System; one (1) TagAlert temperature monitor

**Storage** – Store at or below 25°C (77°F); excursions permitted to 30°C (86°F).

3.0 INDICATIONS

The Absorb BVS is a temporary scaffold indicated for improving coronary luminal diameter that will eventually resorb and potentially facilitate normalization of vessel function in patients with ischemic heart disease due to *de novo* native coronary artery lesions. The treated lesion length should be less than the nominal scaffolding length (8 mm, 12 mm, 18 mm, 23 mm, 28 mm) with reference vessel diameters ≥ 2.0 mm and ≤ 3.8 mm.

4.0 CONTRAINDICATIONS

The Absorb BVS System is contraindicated for use in:

- Patients in whom antiplatelet and / or anticoagulant therapy is contraindicated
- Patients with a known hypersensitivity or contraindication to aspirin, both heparin and bivalirudin, clopidogrel, ticlopidine, prasugrel, and ticagrelor, everolimus, poly (L-lactide), poly (D,L-lactide), or platinum, or with contrast sensitivity, who cannot be adequately premedicated

5.0 WARNINGS

- The potential long-term benefits of the Absorb coronary device in comparison to the available therapies are to be confirmed in additional clinical studies.
- This product should not be used in patients who are not likely to comply with the recommended antiplatelet therapy.
- Judicious selection of patients is necessary, since the use of this device carries the associated risk of scaffold thrombosis, vascular complications, and / or bleeding events.
- Oral administration of everolimus in combination with cyclosporine has been associated with increased serum cholesterol and triglycerides. Therefore, patients should be monitored for changes in lipid profiles.
- Persons allergic to poly (L-lactide), poly (D,L-lactide), everolimus, or platinum may suffer an allergic reaction to this implant.
- It is not recommended to treat patients having a lesion with excessive tortuosity proximal to or within the lesion.
- Devices (i.e., guide sheaths) that decrease the inner diameter of the guide catheter outside of the Absorb BVS System minimum guide catheter compatibility (Table 2) must not be used with the Absorb BVS System. Do not insert a 5-in-6, or a 6-in-7 guide sheath into a 6F or 7F guiding catheter, as doing so will result in an inner diameter that is too small for use with the Absorb BVS System.
- Balloon dilatation of any cells of a deployed Absorb BVS may cause scaffold damage.
- Careful selection of the target lesion reference vessel diameter to the scaffold diameter, and adequate lesion preparation prior to scaffold implantation is recommended to minimize tissue potential damage to the scaffold during placement. It is not recommended to treat patients having a lesion that prevents complete inflation of an angioplasty balloon (example: a severely calcified lesion that has not had adequate lesion preparation), or a lesion with greater than 40% residual stenosis after predilatation by visual estimation.

6.0 PRECAUTIONS

6.1 Scaffold Handling – Precautions

- For single use, single insertion only.** Do not resterilize or reuse. Note the product "Use by" date on the package.
- Do not remove the scaffold from the delivery system**, as removal may damage the scaffold and / or lead to scaffold embolization. The scaffold system is intended to perform together as a system.
- This delivery system should not be used in conjunction with other stents.
- Special care must be taken not to handle or in any way disrupt the scaffold on the balloon. This is most important during catheter removal from packaging, placement over the guide wire, and advancement through the rotating hemostatic valve adapter and guiding catheter hub.
- Do not manipulate, touch, or handle the scaffold** with your fingers, as this may cause coating damage, contamination, or dislodgement of the scaffold from the delivery balloon.
- Use only the appropriate balloon inflation media. Do not use air or any gaseous medium to inflate the balloon, as this may cause uneven expansion and difficulty in deployment of the scaffold.
- Implantation of the scaffold should be performed **only** by physicians who have received appropriate training.
- Scaffold placement should only be performed at hospitals where access to emergency coronary artery bypass graft surgery (CABG) is available.
- Subsequent restenosis may require repeat dilatation of the arterial segment containing the scaffold. The long-term outcome following repeat dilatation of endothelialized scaffolds is unknown at present.

6.2 Scaffold Placement – Precautions

- Do not prepare or preinflate the delivery system prior to scaffold deployment, other than as directed. Use balloon pumping technique described in section 10.4.4 Delivery System Preparation.
- Size the reference target lesion diameter appropriately to avoid overexpanding the scaffold to ensure adequate scaffold apposition. This will reduce the risk of causing scaffold damage.
- When introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the scaffold from the balloon.
- Do not torque the catheter more than one (1) full turn.
- Use caution when advancing the Absorb BVS across the lesion. Multiple attempts to cross a lesion may lead to scaffold damage or dislodgement.
- Implanting a scaffold may lead to dissection of the vessel distal and / or proximal to the scaffold and may cause acute closure of the vessel, requiring additional intervention (CABG, further dilatation, placement of additional scaffolds, or other).
- Do not expand the scaffold if it is not properly positioned in the vessel. (See section 6.4 Scaffold / System Removal – Precautions.)
- When treating multiple lesions within the same vessel, scaffold / stent the distal lesion prior to scaffolding the proximal lesion. Scaffolding / stenting in this order obviates the need to cross the proximal scaffold in placement of the distal scaffold / stent, and reduces the chance of damaging or dislodging the proximal scaffold.
- Avoid scaffolding across any side branches  $\geq 2.0$  mm in diameter:
  - The inflated balloon diameter of the system used to deploy the scaffold should approximate the diameter of the vessel. To ensure full expansion of the scaffold, the balloon should be inflated to a minimum of nominal pressure.
  - Do not exceed the Rated Burst Pressure (RBP) as indicated on the product label. Monitor balloon pressures during inflation. Use of pressures higher than specified on the product label may result in a ruptured balloon, with possible intimal damage and dissection.
- Scaffold retrieval methods (use of additional wires, snares, and / or forceps) may result in additional trauma to the coronary vasculature and / or the vascular access site. Complications may include bleeding, hematoma, or pseudoaneurysm.
- When multiple Absorb Bioresorbable Vascular Scaffolds and drug-eluting stents are required, only Absorb Bioresorbable Vascular Scaffolds or everolimus-eluting stents must be used. Potential interaction with other drug-eluting stents or coated stents has not been evaluated and should be avoided.
- The extent of the patient's exposure to drug and polymer is directly related to the number of scaffolds implanted. A patient can receive up to a total scaffolding length of 94 mm.
- The safety and effectiveness of the Absorb BVS in patients with prior brachytherapy of the target lesion or the use of brachytherapy for treated site restenosis in an Absorb BVS have not been established. Both vascular brachytherapy and the Absorb BVS alter arterial remodeling. The combination between these two treatments has not been established.
- Postdilatation with a noncompliant balloon is recommended following instructions in section 10.7, as long as the postdilated segment is within the allowable expansion limits of the scaffold.
- An unexpanded scaffold may be retracted into the guiding catheter **one time only**. An unexpanded scaffold should not be reintroduced into the artery once it has been pulled back into the guiding catheter. Subsequent movement in and out through the distal end of the guiding catheter should not be performed, as the scaffold may be damaged or dislodged during retraction back into the guiding catheter.
- Should resistance be felt at any time during removal of the undeployed Absorb BVS System, please refer to the steps provided in section 6.4 Scaffold / System Removal – Precautions.
- In the rare event of acute occlusion following scaffold placement, a bailout implant may be inserted and deployed within the scaffold such that the Absorb BVS is completely covered by the bailout implant. All abrupt closures must be treated as an emergency per the hospital standard of care.  
**Note:** It is recommended that bailouts be done with a metallic everolimus-eluting stent of appropriate size.

6.3 Use in Conjunction with Other Procedures

- While vessel preparation in complex lesions may include the use of various mechanical atherectomy devices, the safety and effectiveness have not been formally established in clinical trials with use of either mechanical atherectomy devices (directional atherectomy catheters, rotational atherectomy catheters) or laser angioplasty catheters in conjunction with Absorb BVS implantation.

6.4 Scaffold / System Removal – Precautions

- Scaffold delivery system removal prior to scaffold deployment:
  - If removal of a scaffold system is required prior to deployment, ensure that the guide catheter is coaxially positioned relative to the scaffold delivery system and cautiously withdraw the scaffold delivery system into the guiding catheter. Should **unusual resistance** be felt at any time when withdrawing the scaffold into the guide catheter, the scaffold delivery system and the guide catheter should be **removed as a single unit**. This should be done under direct visualization with fluoroscopy.
- Withdrawal of the scaffold delivery system / postdilatation balloon from the deployed scaffold:
  1. Deflate the balloon by pulling negative on the inflation device. Larger and longer balloons will take more time (up to 30 seconds) to deflate than smaller and shorter balloons. Confirm balloon deflation through fluoroscopy.
  2. Position the inflation device to "negative" or "neutral" pressure.
  3. Stabilize guide catheter position just outside coronary ostium and anchor in place. Maintain guide wire placement across scaffold segment.
  4. Gently remove the scaffold delivery system with slow and steady pressure.
  5. Tighten the rotating hemostatic valve.**Notes:**
  - 1) If, during withdrawal of the catheter from the deployed scaffold, resistance is encountered, use the following steps to improve balloon rewrap:
    - Re-inflate the balloon up to nominal pressure, deflate and change pressure to neutral.
    - Repeat steps 1 through 5 above.
  - 2) After successful withdrawal of the balloon from the deployed scaffold, should any resistance be felt at any time when withdrawing the scaffold delivery system or postdilatation balloon into the guide catheter, **remove the entire system as a single unit**.

- Failure to follow these steps and / or applying excessive force to the delivery system can potentially result in loss or damage to the scaffold and / or delivery system components.
- If it is necessary to retain guide wire position for subsequent artery / lesion access, leave the guide wire in place and remove all other system components.

6.5 Post Implant – Precautions

- If necessary to cross a newly deployed scaffold with a guide wire, balloon, delivery system, or imaging catheters, exercise care to avoid disrupting the scaffold geometry.
- Antiplatelet therapy should be administered post-procedure (see section 9.1 Individualization of Treatment).

6.6 MRI Statement

Comparison to the XIENCE PRIME metallic stents\* indicates that the Absorb BVS is MR Conditional and can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla or less
- Spatial gradient of 2500 Gauss/cm or less
- Maximum whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode) for up to 15 minutes of scanning for each sequence

The Absorb BVS should not migrate in this MRI environment nor should there be safety concerns related to scaffold or marker heating. No concerns related to the effect of the above MRI environment on the drug content or drug release rate are anticipated. Due to its material properties, the same MRI condition restrictions as are appropriate for the XIENCE PRIME stent apply more conservatively to the Absorb BVS. It is unknown if MRI image quality will be compromised if the area of interest is in the exact same area or relatively close to the position of the Absorb BVS.

\*Nonclinical testing was completed on XIENCE PRIME units in single and in overlapped configurations up to 71 mm in length and is on file with Abbott Vascular.

6.7 Drug Interactions

Everolimus is extensively metabolized by the cytochrome P4503A4 (CYP3A4) in the gut wall and liver, and is a substrate for the countertransporter P-glycoprotein. Everolimus has also been shown to reduce the clearance of some prescription medications when it was administered orally along with cyclosporine (CsA). Hence, everolimus, when prescribed as an oral medication, may interact with other medications that include (but are not restricted to) inhibitors and inducers of CYP3A4 isozymes; absorption and subsequent elimination of everolimus may be influenced by drugs that affect these pathways. Formal drug interaction studies have not been performed with the Absorb BVS System. Therefore, due consideration should be given to the potential for both systemic and local drug interactions in the vessel wall, when deciding to place the Absorb BVS in a subject taking a drug with known interaction with everolimus.

Everolimus, when prescribed as an oral medication may interact with the following drugs or foods<sup>1,2</sup>:

- CYP3A4 / P-gp inhibitors that may increase everolimus blood concentrations:
  - Immunosuppressant cyclosporine
  - Antifungal agents (e.g., fluconazole, ketoconazole, itraconazole ketoconazole)
  - Macrolide antibiotics (e.g., clarithromycin, erythromycin)
  - Calcium channel blockers (e.g., verapamil, nifedipine, diltiazem)
  - Protease inhibitors (e.g., nelfinavir, indinavir, amprenavir)
  - Other substances (e.g., cisapride, metoclopramide, bromocriptine, cimetidine, danazol)
- CYP3A4 inducers that may decrease everolimus drug concentrations:
  - Antibiotics (e.g., rifampin, rifabutin)
  - Anticonvulsants (e.g., carbamazepine, phenobarbital, phenytoin)
  - Non-nucleoside reverse transcriptase inhibitors (e.g., efavirenz, nevirapine)
  - Glucocorticoids (e.g., dexamethasone, prednisone, prednisolone)
  - HMGCoA reductase inhibitors (e.g., simvastatin, lovastatin)
  - Other (e.g., St. John's Wort)

6.8 Pregnancy

This product (the Absorb BVS System) and everolimus have not been tested in pregnant or nursing women or in men intending to father children. Effects on the developing fetus have not been studied. While there is no contraindication, the risks and reproductive effects are unknown at this time.

7.0 ADVERSE EVENTS

7.1 Observed Adverse Events

Adverse events observed in the ABSORB clinical trials that are related to the key clinical outcomes of death, cardiac death, myocardial infarction (MI)(Q-wave and non-Q-wave), target lesion revascularization (TLR) (percutaneous coronary intervention [PCI] or coronary artery bypass graft), scaffold thrombosis, and ischemia-driven Major Adverse Cardiac Event (MACE) (composite of cardiac death, MI, ischemia-driven TLR [ID-TLR]) are presented in Table 4. All other adverse events are included in section 7.2.

7.2 Potential Adverse Events

Adverse events that may be associated with PCI, treatment procedures and the use of a coronary scaffold in native coronary arteries include the following, but are not limited to:

- Abrupt closure
- Access site complications
- Myocardial infarction
- Allergic reaction or hypersensitivity to contrast agent or polymer poly (L-lactide) (PLLA), polymer poly (D,L-lactide) (PDLLA), and drug reactions to antiplatelet drugs or contrast agent
- Aneurysm
- Arterial perforation
- Arterial rupture
- Arteriovenous fistula
- Arrhythmias, including atrial and ventricular
- Bleeding complications, which may require transfusion
- Cardiac arrest
- Cardiac, pulmonary, or renal failure
- Cardiac tamponade
- Coronary artery spasm
- Coronary or scaffold embolism
- Coronary or scaffold thrombosis
- Death
- Dissection of the coronary artery
- Distal emboli (air, tissue or thrombotic)
- Emergent or non-emergent surgery
- Fever
- Hypotension / hypertension
- Infection and pain
- Injury to the coronary artery

- Ischemia, myocardial
- Nausea and vomiting
- Palpitations
- Pericardial effusion
- Peripheral ischemia (due to vascular injury)
- Peripheral nerve injury
- Pulmonary edema
- Pseudoaneurysm
- Renal insufficiency / failure
- Restenosis of scaffolded segment
- Shock
- Stroke / cerebrovascular accident and TIA
- Total occlusion of coronary artery
- Unstable or stable angina pectoris
- Vascular complications, including entry site, which may require vessel repair
- Ventricular arrhythmias, including ventricular fibrillation and ventricular tachycardia
- Vessel dissection

Adverse events associated with daily oral administration of everolimus can be found in the informational brochure for the drug<sup>3,4</sup> and include the following, but are not limited to:

- Abdominal pain
- Acne
- Anemia
- Angioneurotic edema
- Coagulopathy
- Diarrhea
- Edema
- Embryo-fetal toxicity
- Hemolysis
- Hemolytic uremic syndrome
- Hepatic disorders
- Hepatitis
- Hypercholesterolemia
- Hyperlipidemia
- Hypertension
- Hypertriglyceridemia
- Hypogonadism male
- Infection
- Interstitial lung disease
- Jaundice
- Leukopenia
- Laboratory test alterations (elevations of serum creatinine, blood glucose, and lipids; decreases in hemoglobin, neutrophils, and platelets)
- Liver function test abnormal
- Lymphocele
- Myalgia
- Nausea
- Oral ulcerations
- Pain
- Pancreatitis
- Pericardial effusion
- Pleural effusion
- Pneumonia
- Non-infectious pneumonitis
- Pulmonary alveolar proteinosis
- Pyelonephritis
- Rash
- Renal failure
- Renal tubular necrosis
- Sepsis
- Surgical wound complication
- Thrombocytopenia
- Thrombotic thrombocytopenic purpura
- Urinary tract infection
- Venous thromboembolism
- Viral, bacterial, and fungal infections
- Vomiting
- Wound infection

8.0 CLINICAL INVESTIGATIONS OF ABSORB BVS SYSTEM

8.1 ABSORB Clinical Trial

The ABSORB clinical investigation is a prospective, single-arm, open-label, multicenter, international clinical study to evaluate safety and performance of the Absorb BVS in coronary arteries. The study was designed to enroll up to 30 patients (Cohort A) followed by approximately 80 patients (Cohort B) in approximately 10 clinical sites in Europe and Asia Pacific region.

**Eligible Criteria:** Patients at least 18 years of age with evidence of angina, myocardial ischemia, or a positive functional test; female patients with childbearing potential with a negative pregnancy test within 7 days prior to the index procedure and effective contraception during participation in this clinical investigation; and patients who signed an informed consent prior to enrollment. Angiographic inclusion criteria included: Nominal vessel diameter of 3.0 mm; lesion length  $\leq 8$  mm by visual estimation extended to  $\leq 14$  mm for the 3.0 x 18 mm scaffold; % diameter stenosis (%DS) of  $\geq 50\%$  and  $< 100\%$ ; TIMI flow of  $\geq 1$ . Key angiographic exclusion criteria included: aorto-ostial location; left main location within 2 mm of the origin of the LAD or LCX; excessive tortuosity proximal to or within the lesion; extreme angulation ( $\geq 90^\circ$ ); heavy calcification; restenosis from previous intervention; target vessel containing thrombus; other clinically significant lesions in the target vessel or side branch.

**Treatment Strategy:** Predilatation of the target lesion was mandatory. Planned overlapping of Absorb BVS was disallowed. Any bailout must be done with overlapping a XIENCE V stent of appropriate length, and if not available, a CYPHER® sirolimus-eluting stent; bailout with Absorb BVS was not permitted. Postdilatation was performed at operator discretion, but only using balloons sized to fit within the boundaries of the scaffold.

**Antiplatelet Regimen:** Subjects not on chronic antiplatelet medications had to receive a loading dose of clopidogrel bisulfate  $\geq 300$  mg and aspirin  $\geq 300$  mg to 72 hours prior to the index procedure if possible, but no later than 1 hour after the procedure. All patients were required to receive anticoagulation and other therapy during scaffold implantation according to the standard of care at the clinical site. All patients were to be maintained on 75 mg clopidogrel bisulfate daily for a minimum of 6 months and  $\geq 75$  mg of aspirin daily for the length of the clinical investigation (5 years). Patients who developed sensitivity to clopidogrel bisulfate were to be switched to ticlopidine hydrochloride at a dose in accordance with standard hospital practice.

8.2 ABSORB Cohort A

8.2.1 Methodology and Current Status

A total of 30 patients in Cohort A were enrolled between March 7, 2006 and July 18, 2006 at four clinical sites in Europe and New Zealand.

The ABSORB Cohort A Trial implanted the Absorb Cohort A BVS in patients with a single *de novo* native coronary artery lesion. Enrollment started with the 3.0 x 12 mm trial. The 3.0 x 18 mm size became available later during enrollment and was used in only 2 patients.

Clinical follow-up through 30, 180, 270 days, 1 year, 18 months, 2 years, 3 years, 4 years, and 5 years, and angiographic, IVUS, IVUS-VH, palpography, OCT data at 180 days and 2 years are available. MACE data at 18 months and coronary vasomotor test data at 2 years are also available. No further clinical observations will be performed for ABSORB Cohort A since 5-year follow-up was the last visit time point.

8.2.2 Clinical Outcome at 5 Years

At 5 years, the ischemia-driven MACE (defined as the composite endpoint of cardiac death, MI or ID-TLR) rate was 3.4% (Table 4). There was only one non-Q-wave MI (peak troponin 2.21  $\mu\text{g/L}$ ) related to the treatment of a non-flow-limiting stenosis (QCA diameter stenosis 42%) of an Absorb BVS implanted 46 days earlier, in a patient who had one episode of angina at rest without electrographic evidence of ischemia. For a perceived safety reason,

<sup>3</sup> Certican® drug prescribing sheet, date of most recent amendment May 21, 2013. Novartis.

<sup>4</sup> Aflinor® summary of product characteristics, date of authorization November 2013. Novartis.

the polymeric scaffold was covered by a drug-eluting metallic stent. No new MACE was recorded between 6 months and 5 years. No instances of scaffold thrombosis arose according to the protocol or the Academic Research Consortium (ARC) definitions.

| Table 4: Clinical Outcome at 5 Years                                                                                                                                                                         |                 |                    |                  |                  |                  |                  |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|--------------------|------------------|------------------|------------------|------------------|
| Hierarchical Rates                                                                                                                                                                                           | 6 Months (N=30) | 12 Months (N=29) † | 2 Years (N=29) † | 3 Years (N=29) † | 4 Years (N=29) † | 5 Years (N=29) † |
| Cardiac Death, %                                                                                                                                                                                             | 0               | 0                  | 0                | 0                | 0                | 0                |
| MI, % (n)                                                                                                                                                                                                    | 3.3 (1)*        | 3.4 (1)*           | 3.4 (1)*         | 3.4 (1)*         | 3.4 (1)*         | 3.4 (1)*         |
| Q-Wave MI                                                                                                                                                                                                    | 0               | 0                  | 0                | 0                | 0                | 0                |
| Non Q-Wave MI                                                                                                                                                                                                | 3.3 (1)*        | 3.4 (1)*           | 3.4 (1)*         | 3.4 (1)*         | 3.4 (1)*         | 3.4 (1)*         |
| Ischemia-Driven TLR, %                                                                                                                                                                                       | 0               | 0                  | 0                | 0                | 0                | 0                |
| by PCI                                                                                                                                                                                                       | 0               | 0                  | 0                | 0                | 0                | 0                |
| by CABG                                                                                                                                                                                                      | 0               | 0                  | 0                | 0                | 0                | 0                |
| Ischemia-Driven MACE (cardiac death, MI or Ischemia-driven TLR), % (n)                                                                                                                                       | 3.3 (1)*        | 3.4 (1)*           | 3.4 (1)*         | 3.4 (1)*         | 3.4 (1)*         | 3.4 (1)*         |
| Scaffold thrombosis, %                                                                                                                                                                                       | 0               | 0                  | 0                | 0                | 0                | 0                |
| <b>Notes:</b>                                                                                                                                                                                                |                 |                    |                  |                  |                  |                  |
| • Denominator excludes subjects without follow-up data and who didn't experience any of the following events: Death, MI, any revascularization (TLR, TVR, NTVR).                                             |                 |                    |                  |                  |                  |                  |
| † One subject officially withdrew from the study after 6 months.                                                                                                                                             |                 |                    |                  |                  |                  |                  |
| * This subject also underwent a TLR, not qualified as ID-TLR (DS = 42%) followed by post-procedural troponin elevation qualified as non-Q MI and died from his Hodgkin's disease at 888 days post-procedure. |                 |                    |                  |                  |                  |                  |

8.2.3 Angiographic, IVUS, and OCT Outcomes at 180 Days and 2 Years

QCA results were collected from baseline, post procedure, 180 days and 2 years. The 180-day data demonstrated an acceptable in-scaffold late loss of  $0.43 \pm 0.37$  mm<sup>2</sup> possibly driven by bioactive remodeling or mechanical late recoil. From 180 days to 2 years, the change in in-scaffold late loss was very low ( $0.48 \pm 0.28$  mm<sup>2</sup>). The 180-day grey-scale intravascular ultrasound (IVUS) analysis showed a significant reduction in the average lumen area ( $6.04 \pm 1.12$  mm<sup>2</sup> after procedure vs.  $5.19 \pm 1.33$  mm<sup>2</sup> at 180-day;  $p < 0.001$ ). The vessel area remained constant between baseline and 180-day follow-up ( $13.49 \pm 3.74$  mm<sup>2</sup> vs.  $13.79 \pm 3.84$  mm<sup>2</sup>), demonstrating the absence of significant expansive or constrictive remodeling. At 2 years, the main observation provided by grey-scale IVUS was the increase in minimal luminal area and average luminal area / volume together with a significant decrease in plaque area / volume between 6 months and 2 years. The 180-day OCT group (n = 13) provided a total of 671 apparent struts for visual evaluation, which showed that 93% of struts evaluated were well apposed to the vessel wall, and 99% of the struts were covered by tissue. The serial OCT group (n = 7) had serial data postprocedure, at 180 days and at 2 years (intention-to-treat population). The number of apparent struts decreased from 403 at baseline to 368 at 180 days and to 264 at 2-years follow-up (34.5% reduction over two years), all well covered and apposed to the vessel wall. Precinical analysis directly comparing pre-explanation OCT to post-explanation histology indicates there is no correlation between the presence or absence of OCT-visible features (apparent struts) and the presence or absence of polylactide polymer in the vessel. Hence, visual OCT features (apparent struts) are not always indicative of lack of resorption, but absence of the apparent struts by OCT confirms resorption. The lumen shape was regular with smooth, well delineated borders in all cases, and no intraluminal tissue was observed. Importantly, minimal and mean luminal area decreased significantly between post procedure and 180 days, but numerically enlarged between 180 days and 2 years.


8.2.4 Vasoconstrictor Function Results at 2 Years

Vasoconstrictor function proximal, within, and distal to the treated (scaffolded) segments at 2 years, was evaluated with either the endothelium-independent vasoconstrictor methylethylgonovine maleate (Methergine, Novartis, Basel, Switzerland), or the endothelium-dependent vasoactive agent acetylcholine (Ovisot, Daiichi-Sankyo, Tokyo, Japan), depending on local practice.

In the Methergine group (n = 7), significant vasoconstriction was observed in proximal and scaffolded segments (Figure 1). In the acetylcholine group (n = 9), although the overall change in vessel dimension was not statistically significant, vasodilatation was observed in five patients. Nitrates induced a significant vasodilatation in the scaffolded and distal segments, relative to baseline (pre-acetylcholine) at 2 years. These results suggest the potential restoration of vasoconstrictor function in the treated segments when the Absorb Cohort A BVS had resorbed. Further data on vascular responses will be collected from ongoing clinical trials.

Paired comparisons between the different time points were done by a Wilcoxon's signed rank test for continuous variables.

Figure 1: Results of Acetylcholine and Methergine testing in proximal, scaffolded, and distal segments

8.2.5 Discussion

In this prospective, single-group, open-label study, the Absorb BVS demonstrated acute success and safety following the first 30 days after implantation: 100% procedure success; 94% device success; and successful revascularization of the target lesion evident from post procedure angiographic and IVUS analyses. The excellent clinical safety persists through 5 years with no cardiac deaths, no ID-TLRs or scaffold thrombosis recorded, and only one non-Q-wave MI by 5 years. Despite the discontinuation of thienopyridine drugs, no scaffold thrombosis occurred, and no additional ischemia-driven target lesion restenosis was evident by 5 years.

The use of multiple imaging modalities yielded several important findings. At 180 days, an acceptable angiographic in-scaffold late loss of 0.43 mm was seen; this was higher than previously observed in XIENCE V at 6 months (SPIRIT FIRST), and possibly driven by bioactive remodeling or mechanical late recoil. The IVUS results showed low intrascaffold neointimal hyperplasia: 5.32% in-scaffold volume obstruction, 4.09 mm<sup>2</sup> in-scaffold neointimal volume. This confirmed a positive drug effect in inhibiting restenosis.

At 2 years, angiographic results showed an acceptable 2-year in-scaffold late loss of 0.48 mm with minimal changes from 6 months (0.43 mm) to 2 years. By contrast, both the IVUS and OCT data showed late luminal enlargement from 6 months to 2 years. The contrasting findings in late lumen dimensional changes between angiography and intracoronary imaging await confirmation in future clinical studies.

At 2 years, reduction in molecular weight and mass had occurred to such an extent that echogenicity was lost and that struts were no longer recognizable by intravascular ultrasound, leaving behind few IVUS-visible features. OCT-visible features related to vessel healing were present in some patients. In these patients, OCT showed an optically homogeneous vessel wall structure that, taken together with the potential restoration of vessel movement, suggests healing of the artery.

8.3 ABSORB Cohort B

8.3.1 Methodology and Current Status

Based upon the strong safety profile observed for the Absorb Cohort A BVS System in Cohort A, Cohort B of the ABSORB trial was initiated on March 19, 2009 to evaluate the Absorb BVS System in a prospective, open-label, multicenter registry. Subjects with up to two *de novo* native coronary artery lesions in separate epicardial vessels with visually estimated nominal vessel diameters of 3.0 mm and lesion(s) length  $\leq 14$  mm were enrolled, and received a single 3.0 x 18 mm Absorb BVS per lesion treated. Twelve clinical sites located in Europe, Australia, and New Zealand participated in this study. Enrollment of 101 patients (45 patients in Group 1 with imaging follow-up at 180 days and 2 years, and 56 patients in Group 2 with imaging follow-up at 1 year and 3 years) in Cohort B was completed on November 6, 2009. Currently, clinical data are available for post-procedure, 6 months, 9 months, 1 year, 2 years and 3 years. Additionally, 2-year angiographic and IVUS results are available for Group 1 subjects, and similar data at 1 year and 3 years are available for Group 2 subjects. Follow-up is ongoing. Patients will be followed through 5 years.

8.3.2 Clinical Outcomes for Full Cohort B

Table 5 shows clinical outcomes through 3 years for all subjects in the ABSORB Cohort B Trial (101 subjects).

| Table 5: Hierarchical Clinical Outcomes for Cohort B (ITT Population)   |                    |                     |                     |                   |                     |                     |
|-------------------------------------------------------------------------|--------------------|---------------------|---------------------|-------------------|---------------------|---------------------|
|                                                                         | 30 Days<br>N = 101 | 6 Months<br>N = 101 | 9 Months<br>N = 101 | 1 year<br>N = 101 | 2 years<br>N = 100* | 3 years<br>N = 100* |
| Cardiac Death, %                                                        | 0                  | 0                   | 0                   | 0                 | 0                   | 0                   |
| MI, % (n)                                                               | 2.0 (2)            | 3.0 (3)             | 3.0 (3)             | 3.0 (3)           | 3.0 (3)             | 3.0 (3)             |
| Q-Wave MI                                                               | 0                  | 0                   | 0                   | 0                 | 0                   | 0                   |
| Non Q-Wave MI                                                           | 2.0 (2)            | 3.0 (3)             | 3.0 (3)             | 3.0 (3)           | 3.0 (3)             | 3.0 (3)             |
| Ischemia-driven TLR, % (n)                                              | 0                  | 2.0 (2)             | 2.0 (2)             | 4.0 (4)           | 6.0 (6)             | 7.0 (7)             |
| by PCI                                                                  | 0                  | 2.0 (2)             | 2.0 (2)             | 4.0 (4)           | 6.0 (6)             | 7.0 (7)             |
| by CABG                                                                 | 0                  | 0                   | 0                   | 0                 | 0                   | 0                   |
| Ischemia-driven MACE, (cardiac death, MI or ischemia-driven TLR), % (n) | 2.0 (2)            | 5.0 (5)             | 5.0 (5)             | 6.9 (7)           | 9.0 (9)             | 10.0 (10)           |
| Scaffold thrombosis (%)                                                 | 0                  | 0                   | 0                   | 0                 | 0                   | 0                   |

Data are % (number of patients); MACE = (cardiac death, MI or Ischemia-driven TLR)  
\* One subject was lost to follow-up after the 2-year visit

For the full Cohort B (101 patients), there were three non-Q-wave MIs that occurred, two during hospitalization, and one at 43 days post-procedure. Seven ID-TLR by PCI have been reported through 3 years. The overall MACE rate at 3 years was 10.0% and there have been no cases of cardiac death or scaffold thrombosis per protocol or per ARC definitions.

The 3-year full Cohort B (101 patients) outcomes are numerically higher than the Cohort A (30 patients) outcomes at 3 years. The clinical outcomes at 36 months reveal a MACE rate of 3.4% and 10.0% in the ABSORB Cohort A and full Cohort B groups, respectively. Furthermore, there were no cases of cardiac death or scaffold thrombosis in either group. The Cohort B clinical results continue to support the performance and safety established in the Cohort A investigation.

8.3.3 Angiographic, IVUS, and OCT Outcomes at 180 Days, 1, 2 and 3 Years

The 180-day angiographic results from Cohort B Group 1 (Cohort B1) demonstrated a late loss of 0.19 mm (N = 42) which compares well to the 0.10 mm late loss (N = 23) of the 3.0 x 18 mm XIENCE V Everolimus Eluting Coronary Stent System (XIENCE V EECSS) observed in the SPIRIT FIRST-1n-Min trial, and show favorable improvement from the 0.43 mm late loss (N = 26) from Cohort A.

At 180 days, the IVUS results from Cohort B1 showed limited intrascaffold neointimal hyperplasia, the volume obstruction (VO) was 1.2% (N = 40) and the neointimal hyperplasia area was 0.08 mm<sup>2</sup> (N = 40). These results compare favorably with XIENCE V EECSS in SPIRIT FIRST (VO of 8.0% [N = 21] and NIH area of 0.56 mm<sup>2</sup> [N = 21]) and with Cohort A (VO of 5.3% [N = 25] and NIH area of 0.29 mm<sup>2</sup> [N = 25]). As observed in Cohort A, the 6-month IVUS results showed a significant reduction in the average lumen area ( $6.60 \pm 1.22$  mm<sup>2</sup> [N = 37] after procedure vs.  $6.37 \pm 1.12$  mm<sup>2</sup> [N = 37] at 180-day;  $p = 0.0049$ ). The vessel area remained comparable between baseline and 180-day follow-up ( $14.22 \pm 3.75$  mm<sup>2</sup> [N = 37] vs.  $14.49 \pm 3.67$  mm<sup>2</sup> [N = 37]) demonstrating the absence of significant expansive or constrictive remodeling.

The 1-year angiographic results from Cohort B Group 2 (Cohort B2) demonstrated a late loss of 0.27 mm (N = 56) which compares well to the 0.23 mm late loss (N = 22) of the 3.0 x 18 mm XIENCE V observed in the SPIRIT FIRST trial. At 1 year, the IVUS results from Cohort B2 (N = 54) showed limited intrascaffold neointimal hyperplasia, the percentage of in-scaffold area obstruction and the neointimal hyperplasia area were minimal (1.43% and 0.09 mm<sup>2</sup>, respectively). Over time there was a statistically significant but not clinically relevant increase in total plaque area with a concomitant increase in vessel area. Furthermore, average lumen area remained stable from post-procedure ( $6.31 \pm 0.95$  mm<sup>2</sup>) to 1 year ( $6.33 \pm 1.17$  mm<sup>2</sup>). By OCT, there was no significant change in the mean and minimum scaffold area between post procedure and 1 year (N = 21). There was a reduction in minimal lumen area and mean luminal area decreased significantly as a result of mean neointimal area growth ( $1.34 \pm 0.67$  mm<sup>2</sup>) but was not clinically significant. The lumen area stenosis increase was statistically significant from 20.2% post-procedure to 26.9% at 1 year but not clinically relevant. These data further supported that there was no reduction in scaffold area.

Imaging results at 2 years from Cohort B Group 1 (Cohort B1) demonstrated in-scaffold angiographic late loss of 0.27 mm (N = 38), which matches the late loss observed at 1 year for Cohort B Group B2 but is lower than the 0.33 mm (N = 83) reported from XIENCE V in SPIRIT II. The IVUS data at 2 years (N = 33) revealed a finding unique to Absorb BVS, which is mean scaffold area enlargement between 6 months and 2 years (2.18 vs. 7.08 mm<sup>2</sup>;  $p < 0.0001$ ). This enlargement was greater than the observed increase in neointimal hyperplasia and was accompanied by late increased mean lumen area (6.36 versus 6.85 mm<sup>2</sup>;  $p = 0.0105$ ). IVUS findings were reported with comparative paired (serial) analysis at post-procedure, 6 months, and 2 years. Similarly, OCT data (N = 23) were reported with paired (serial) analysis, and showed that from 6 months to 2 years the mean strut core area decreased from 0.20 to 0.15 mm<sup>2</sup> ( $p < 0.0001$ ), the mean neointimal area increased but remained minimal at 2 years ( $1.43$  versus  $2.11$  mm<sup>2</sup>;  $p < 0.0001$ ) and mean lumen area was unchanged. Moreover, the mean scaffold area increased from 7.47 to 8.24 mm<sup>2</sup> ( $p = 0.0155$ ) from post-procedure to 2 years indicating loss of mechanical integrity of the scaffold and potential expansion of the vessel. Scaffold strut coverage by neointima was near completed at 6 months (98%) and 2 years (99%).

At 3 years post-procedure, IVUS results from Cohort B Group 2 (Cohort B2, N=44) revealed scaffold enlargement (average scaffold area of  $6.35 \pm 0.99$  mm<sup>2</sup> at 1 year vs.  $7.08 \pm 1.55$  mm<sup>2</sup> at 3 years;  $p < 0.0001$ ) and increased average lumen area ( $6.35 \pm 1.17$  mm<sup>2</sup> at 1 year vs.  $6.81 \pm 1.62$  mm<sup>2</sup> at 3 years;  $p = 0.0006$ ). The late lumen enlargement is of particular interest because the vessel area was preserved (average vessel area of  $14.43 \pm 2.64$  mm<sup>2</sup> at 1 year vs.  $14.58 \pm 2.67$  mm<sup>2</sup> at 3 years;  $p = 0.407$ ) despite an increase in neointimal area ( $0.08 \pm 0.13$  mm<sup>2</sup> at 1 year vs.  $0.28 \pm 0.41$  mm<sup>2</sup> at 3 years;  $p < 0.0001$ ). OCT results at 3 years (N=18) also revealed scaffold enlargement between 1 and 3 years (mean scaffold area of  $7.51 \pm 0.95$  mm<sup>2</sup> at 1 year vs.  $8.84 \pm 2.15$  mm<sup>2</sup> at 3 years;  $p = 0.0008$ ). Lumen area remained relatively constant ( $6.01 \pm 1.29$  mm<sup>2</sup> at 1 year vs.  $6.09 \pm 1.67$  mm<sup>2</sup> at 3 years;  $p = 0.679$ ). The vessel area was not measured because OCT cannot demarcate the external elastic lamina. Overall, the neointimal hyperplasia remained small although there was an increase in neointimal hyperplasia between 1 year and 3 years (mean neointimal area of  $1.41 \pm 0.68$  mm<sup>2</sup> at 1 year vs.  $2.35 \pm 0.68$  mm<sup>2</sup> at 3 years;  $p < 0.0001$ ).

8.3.4 Discussion

A total of 101 patients were enrolled with clinical data available out to 3 years for the full Cohort B. In addition, imaging data at 180 days and 2 years are available for the 45 patients in Cohort B1, and at 1 year and 3 years for the 56 patients in Cohort B2. Overall, the Absorb BVS System has demonstrated acute success (98% procedure success and 100% device success) and safety up to 3 years after implantation. At 3 years, the MACE rate remained low at 3.4% and 10.0% in the ABSORB Cohort A and Cohort B groups, respectively. Furthermore, there were no cases of cardiac death or scaffold thrombosis in either cohort. Based on these results, the performance and safety continues to be established in the Cohort A and Cohort B studies.

Imaging of Cohort B indicated that angiographic late loss appears to stabilize between 1 and 3 years. In Group B1 and Group B2, IVUS results indicate a late expansion of scaffold and lumen dimensions between 6 months and 2 years, and between 1 and 3 years, respectively, which was consistent with the scaffold expansion observed by OCT. These results suggest scaffold bioresorption accompanied by a loss of structural integrity. This phenomenon is unique to bioresorbable vascular scaffolds and is in contrast with serial IVUS imaging of metallic DES, where late lumen expansion has not been observed. Follow-up is ongoing. Patients will be followed through 5 years.

8.4 ABSORB EXTEND

The primary objective of the ABSORB EXTEND study is to continue the assessment of the safety and performance of the Absorb BVS System in a larger patient population across more global geographies. This trial collects clinical data without the effects of extensive imaging. Treated lesions are slightly more complex than those included in prior ABSORB trials. Longer lesions are permitted by using planned overlap of the Absorb BVS, or by using longer Absorb BVS sizes as they become available. The following section summarizes the study design and interim clinical outcomes from the ABSORB EXTEND study.

<sup>5</sup> Claessen BE, Beijk MA, Legrand V, Ruzyllo W, Manari A, Varenne O, Suttorp MJ, Tijssen JG, Miquel-Hebert K, Veldhof S, Henriques JP, Serruys PW, Piek JJ. Two-year clinical, angiographic, and intravascular ultrasound follow-up of the XIENCE V everolimus-eluting stent in the treatment of patients with *de novo* native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009 Aug;2(4):339-47.

8.4.1 Methodology and Current Status

ABSORB EXTEND is a prospective, single-arm, open-label clinical study that is planned to register up to 1,000 subjects at up to 10 global sites. Each subject is to receive treatment in up to a maximum of two *de novo* native coronary lesions with each lesion in a different epicardial vessel. Subjects having target lesion length  $\leq$  28 mm and reference vessel sizes that are suitable to be treated with an Absorb BVS can be registered in ABSORB EXTEND. The 3.0 x 18 mm Absorb BVS has been available from the start of the trial registration. The 3.0 x 28 mm, 2.5 x 18 mm, 2.5 x 28 mm, 3.5 x 12 mm, 3.5 x 18 mm and 3.5 x 28 mm have been incorporated into the trial as they became available. Enrollment was completed in ABSORB EXTEND with 807 patients in October 2013. Patients will be followed for up to 3 years.

8.4.2 Summary of Clinical Results out to 1 year (Interim Data Cutoff)

Of the total of 605 patients registered in ABSORB EXTEND as of December 03, 2012, 37-day (30 days  $\pm$  7-day window), 194-day (180 days  $\pm$  14-day window) and 392-day (365 days  $\pm$  28-day window) clinical follow-up data were available for 591, 537 and 450 subjects, respectively (Table 6). The composite endpoints are presented using the protocol definitions for MI. The component endpoints are based on non-hierarchical counts unless otherwise noted.

Table 6: ABSORB EXTEND Subject Counts of Ischemia-driven Adverse Events through 393 Days (All Subjects Registered Population)

| Events                  | ABSORB EXTEND<br>0-37 Days (N=591) | ABSORB EXTEND<br>0-194 Days (N=537) | ABSORB EXTEND<br>0-393 Days (N=450) |
|-------------------------|------------------------------------|-------------------------------------|-------------------------------------|
| Hierarchical MACE       | 2.4% (14/591)                      | 3.0% (16/537)                       | 4.2% (19/450)                       |
| Hierarchical TVF        | 2.4% (14/591)                      | 3.4% (18/537)                       | 4.7% (21/450)                       |
| Non-Hierarchical Counts |                                    |                                     |                                     |
| Cardiac Death           | 0.2% (1/591)                       | 0.2% (1/537)                        | 0.2% (1/450)                        |
| MI                      | 2.4% (14/591)                      | 2.8% (15/537)                       | 2.9% (13/450)                       |
| QMI                     | 0.7% (4/591)                       | 0.6% (3/537)                        | 0.9% (4/450)                        |
| NQMI                    | 1.7% (10/591)                      | 2.2% (12/537)                       | 2.0% (9/450)                        |
| ID-TLR                  | 0.3% (2/591)                       | 0.6% (3/537)                        | 1.8% (8/450)                        |
| ID-non-TL TVR           | 0.0% (0/591)                       | 0.6% (3/537)                        | 0.9% (4/450)                        |
| Scaffold Thrombosis     | 0-30 days                          | 0-194 days                          | 0-393 days                          |
| Definite                | 0.3% (2/591)                       | 0.4% (2/536)                        | 0.7% (3/449)                        |
| Probable                | 0.2% (1/591)                       | 0.2% (1/536)                        | 0.2% (1/449)                        |
| Definite/Probable       | 0.5% (3/591)                       | 0.6% (3/536)                        | 0.9% (4/449)                        |

Note: Subjects are only counted once for each type of event in each time period. Note: All table entries were calculated based on the data from data cut-off date of December 3, 2012.

The 37-day MACE and TVF rates of 2.4% were driven by MI events (10 NQMI and 3 QMI hierarchically) that were all related to the target vessel. There was one cardiac death per ARC definition. The overall non-hierarchical MI rate at 37 days was 2.4% (14/591) (10 NQMI and 4 QMI). The MACE rate increased to 3.0% at 6 months, attributed mostly to the 15 hierarchical MI and 1 cardiac death. The TVF rate at 6 months was 3.4% (18/537 [1 cardiac death, 15 MI, and 2 ID-Non-TL TVR hierarchically]). One patient who died at day 108 was adjudicated as cardiac death per the ARC definition for cardiac death (0.2% at 194 days, 1/537). No Absorb BVS was implanted in the target lesion in this patient because there was difficulty in crossing the lesion with Absorb BVS. A metallic everolimus-eluting stent was implanted instead. Over the same period, the ID-TLR rate was 0.6% (3/537) due to the occurrence of 2 ID-TLR events by PCI and 1 ID-TLR event by CABG. The ID-non-TL TVR rate was 0.6% (3/537) due to the occurrence of 2 ID-non-TL TVR events by CABG and 1 ID-non-TL TVR event by PCI. At 1 year, the MACE rate was 4.2% (19/450 [1 cardiac death, 13 MI (per protocol), and 5 ID-TLR hierarchically]). The TVF rate was 4.7% (21/450 [1 cardiac death, 13 MI (per protocol), 5 ID-TLR, and 2 ID-non-TL TVR hierarchically]). The overall non-hierarchical MI rate (per protocol) at 393 days was 2.9% (13/450) (9 NQMI and 4 QMI).

8.4.3 Summary of Clinical Results out to 2 years (Interim Data Cutoff)

In the subset of first 250 patients that reached 2-year follow-up (758 days) (19 July 2013 data cut), there was one case of cardiac death and the 758-day MACE and rate was low (7.3%). TLF and TVF rates were also low at 6.9% and 8.1%, respectively. Further, the ARC defined definite + probable scaffold thrombosis rate through 758-days for this population was 0.8%.

Table 7: ABSORB EXTEND Subject Counts of Ischemia-driven Adverse Events through 758 Days (250 Subjects Registered Population)

| Events                             | 30 days,<br>n=250<br>(30 $\pm$ 7 days) | 6 months,<br>n=250<br>(180 $\pm$ 14 days) | 1 year,<br>n=250<br>(365 $\pm$ 28 days) | 2 years,<br>n=250<br>(730 $\pm$ 28 days) |
|------------------------------------|----------------------------------------|-------------------------------------------|-----------------------------------------|------------------------------------------|
| Hierarchical, Per protocol         |                                        |                                           |                                         |                                          |
| MACE                               | 2.0% (5/250)                           | 2.8% (7/250)                              | 4.4% (11/250)                           | 7.3% (18/248)                            |
| TVF                                | 2.0% (5/250)                           | 3.2% (8/250)                              | 4.8% (12/250)                           | 8.1% (20/248)                            |
| TLF                                | 2.0% (5/250)                           | 2.8% (7/250)                              | 4.4% (11/250)                           | 6.9% (17/248)                            |
| Non-hierarchical counts            |                                        |                                           |                                         |                                          |
| MI (per protocol)                  | 2.0% (5/250)                           | 2.4% (6/250)                              | 2.8% (7/250)                            | 4.0% (10/248)                            |
| Cardiac death                      | 0.0% (0/250)                           | 0.4% (1/250)                              | 0.4% (1/250)                            | 0.4% (1/248)                             |
| Ischemia-driven TLR                | 0.4% (1/250)                           | 0.4% (1/250)                              | 2.0% (5/250)                            | 4.0% (10/248)                            |
| Ischemia-driven non-TL TVR         | 0.0% (0/250)                           | 0.4% (1/250)                              | 0.8% (2/250)                            | 2.0% (5/248)                             |
| Scaffold Thrombosis (ARC Def/Prob) | 0.4% (1/250)                           | 0.4% (1/249)                              | 0.8% (2/249)                            | 0.8% (2/246)                             |

1. Data are % (number of subjects with event/number of subjects)
2. TLR=target lesion revascularization TV=target vessel
3. MACE=major adverse cardiac events (composite endpoint of Cardiac Death, MI, ischemia-driven TLR)

4. TLF=target lesion failure (composite endpoint of Cardiac Death, TVMI, ischemia-driven TVR)
5. TVF=target vessel failure (composite endpoint of Cardiac Death, MI, ischemia-driven TLR, ischemia-driven Non-TLR TVR)

Conclusion: The interim outcomes from ABSORB EXTEND further demonstrate the 12 and 24 month safety and performance of Absorb BVS in a larger population of subjects where longer lesion can be treated. Patients will be followed for up to 3 years.

8.5 Propensity Score Adjusted Analysis of Absorb BVS System and XIENCE V

Rationale and Method: Abbott Vascular is using its large body of historical clinical data on XIENCE V stents as a control group for comparison against the currently available Absorb BVS clinical data utilizing the propensity score method. Propensity score analysis is a valid statistical method that is well accepted and commonly used in the clinical trial, and in scientific and statistical communities for the comparison of data from non-randomized studies. This approach adjusts for the baseline imbalance between the test and control groups, creates a better match between both groups, and hence, reduces real and potential bias in the comparative analysis. After propensity score adjustment, baseline characteristics between the two groups are more balanced, and therefore the comparisons of the clinical outcomes become more reliable. Key results for a 6-month and a 1-year propensity analysis comparing Absorb BVS to XIENCE V data are presented using the inverse propensity score weighted (IPSW) method. Propensity analysis at later time points will be conducted as additional data from ongoing trials become available. The results of this comparison are being confirmed through ongoing clinical trials.

8.5.1 6-month Propensity Analysis

Thirty-eight patients from ABSORB Cohort B Group 1 and 73 XIENCE V patients treated with 3.0 x 18 mm stents from SPIRIT FIRST and SPIRIT II with 6-month OCA data are used for this 6-month analysis. Their unadjusted 6-month MACE rates are 2.63% (1/38) and 4.11% (3/73), respectively.

Results: The propensity score adjustment yielded adjusted sample sizes of 34 (Absorb BVS) and 66 (XIENCE V) from the 38 and 73, respectively, in the unadjusted arms. The adjusted data were more balanced in baseline patient demographics and angiographic characteristics between both arms. At 6 months, their clinical outcomes and angiographic results were comparable. The adjusted 6-month MACE rates were 2.01% (0.7/34) and 5.27% (3.5/66), respectively, for the Absorb BVS and XIENCE V arms (Table 8) with a lower MACE observed with Absorb BVS, a trend that is consistent with the results from the unadjusted groups.

Table 8: MACE at 6 Months from the IPSW-Adjusted Data

| Clinical Outcome at 6 Months | Absorb BVS System<br>(N = 34) | XIENCE V<br>(N = 66)           | P-value |
|------------------------------|-------------------------------|--------------------------------|---------|
| MACE<br>[95% CI]             | 2.01% (0.7)<br>[0.00%, 6.73%] | 5.27% (3.5)<br>[0.00%, 10.66%] | 0.4398  |

The 6-month in-device late loss in the IPSW-adjusted arms was 0.14  $\pm$  0.18 mm for the Absorb BVS System vs. 0.12  $\pm$  0.26 mm for the XIENCE V ECSSS, representing a non-significant difference in outcomes between the Absorb BVS System and the XIENCE V ECSSS (P = 0.68). Similarly, differences in the in-device % diameter stenosis, in-segment late loss, and in-segment % diameter stenosis were not significant between the two devices in this propensity score adjusted analysis (Table 9).

Table 9: Angiographic Outcomes at 6 Months from the IPSW-Adjusted Data

| QCA Results at 6 Months        | Absorb BVS System<br>(N = 34) | XIENCE V<br>(N = 66) | P-value |
|--------------------------------|-------------------------------|----------------------|---------|
| In-device % Diameter Stenosis  | 17.59 $\pm$ 8.61              | 15.27 $\pm$ 9.50     | 0.2347  |
| In-segment % Diameter Stenosis | 23.87 $\pm$ 9.48              | 23.83 $\pm$ 12.23    | 0.987   |
| In-device Late Loss (mm)       | 0.14 $\pm$ 0.18               | 0.12 $\pm$ 0.26      | 0.6849  |
| In-segment Late Loss (mm)      | 0.11 $\pm$ 0.29               | 0.08 $\pm$ 0.28      | 0.7186  |

8.5.2 1-year Propensity Analysis

Prior to propensity adjustment, a total of 314 Absorb BVS pooled from Cohort B and an interim dataset from ABSORB EXTEND, and 905 XIENCE V patients pooled from SPIRIT FIRST, SPIRIT II, and SPIRIT III with 1-year follow-up were available, irrespective of the device sizes used for their treatment. Their unadjusted 1-year MACE rates are 5.41% (17/314) and 5.64% (51/905), respectively.

Results: After the propensity adjustment, 282 Absorb BVS and 583 XIENCE V patients remained for the comparative analysis. An improved balance was achieved between both cohorts with the most striking changes observed in the lesion characteristics, which were no longer statistically different between Absorb BVS and XIENCE V cohorts. The MI, ID-TLR, MACE, and TLF (Target Lesion Failure defined as cardiac death, target vessel related MI and ID-TLR) rates were comparable between both adjusted cohorts at 37, 193, and 393 days. A summary of key 1-year outcomes is shown in Table 10. After adjustment, the MACE rates remained comparable, with a numerically lower rate obtained for Absorb BVS (4.15%, 11.7/282) versus XIENCE V (5.63%, 32.8/583). The respective scaffold / stent thrombosis rates were not statistically different (0.40% vs. 0.52%, P=0.80).

Table 10: Clinical outcomes at 1 year from the IPSW-Adjusted Cohorts

| Clinical Outcomes at 1 year                             | Absorb BVS<br>(N=282)   | XIENCE V<br>(N=583)     | P-value |
|---------------------------------------------------------|-------------------------|-------------------------|---------|
| MACE<br>[95% CI]                                        | 4.15%<br>[1.82%, 6.47%] | 5.63%<br>[3.76%, 7.50%] | 0.35    |
| MI<br>[95% CI]                                          | 2.64%<br>[0.77%, 4.51%] | 2.52%<br>[1.25%, 3.80%] | 0.92    |
| ID-TLR<br>[95% CI]                                      | 1.95%<br>[0.33%, 3.56%] | 3.10%<br>[1.70%, 4.51%] | 0.33    |
| TLF<br>[95% CI]                                         | 4.15%<br>[1.82%, 6.47%] | 5.09%<br>[3.30%, 6.87%] | 0.54    |
| Definite/Probable Scaffold/Stent Thrombosis<br>[95% CI] | 0.40%<br>[0.0%, 1.14%]  | 0.52%<br>[0.0%, 1.11%]  | 0.80    |

Discussion: The propensity score adjusted method achieved the intent of balancing the baseline characteristics between both XIENCE V and Absorb BVS groups from different ABSORB and SPIRIT trials. The adjusted results showed that Absorb BVS is comparable to

the XIENCE V stent for in-stent late loss at 6 months and for clinical outcomes at 6 months and 1 year after implantation, and further confirms the comparability between Absorb BVS and XIENCE V in the above parameters found in the unadjusted data.

8.5.3 2-year Propensity Analysis

As of the data cut on 19 July 2013, prior to propensity score matched analysis, a total of 250 BVS and 887 XV patients were available through 2 years. Following matching, 178 BVS and 293 XV patients remained for the comparative analysis.

Results: After propensity score matched adjustment, 2-year MACE and TLF rates remained similar between the Absorb BVS cohort and the XIENCE V cohort (6.7% vs. 8.9% for MACE, and 6.2% vs. 8.2% for TLF, respectively). The BVS cohort exhibited a trend of lower ID-TLR and cardiac death rate as compared to the XIENCE V cohort, although the differences were not significant. The 2-year TVF rate was also numerically lower in the BVS group (7.3% vs. 12.3%, p = 0.09). The adjusted definite / probable scaffold / stent thrombosis rates were similar between Absorb BVS and XIENCE V at 2 years (0.6%, 95% CI [0.01%, 3.09%] vs. 1.4%, 95% CI [0.37%, 3.46%], p=0.65). A summary of the key data is shown below in Table 11.

Table 11: Clinical outcomes at 2 years from the IPSW-Adjusted Cohorts

| Clinical Outcomes at 2 years                      | Absorb BVS<br>(N=178)   | XIENCE V<br>(N=293)      | P-value |
|---------------------------------------------------|-------------------------|--------------------------|---------|
| Cardiac Death<br>[95% CI]                         | 0.0%<br>[0.00%, 2.05%]  | 1.4%<br>[0.37%, 3.46%]   | 0.30    |
| MI<br>[95% CI]                                    | 4.5%<br>[1.96%, 8.66%]  | 4.4%<br>[2.38%, 7.47%]   | 1.00    |
| ID-TLR<br>[95% CI]                                | 3.4%<br>[1.25%, 7.19%]  | 3.8%<br>[1.89%, 6.62%]   | 1.00    |
| TLF<br>[95% CI]                                   | 6.2%<br>[3.12%, 10.79%] | 8.2%<br>[5.32%, 11.94%]  | 0.47    |
| MACE<br>[95% CI]                                  | 6.7%<br>[3.53%, 11.48%] | 8.9%<br>[5.88%, 12.73%]  | 0.49    |
| TVF<br>[95% CI]                                   | 7.3%<br>[3.95%, 12.17%] | 12.3%<br>[8.76%, 16.60%] | 0.09    |
| Definite/Probable Scaffold Thrombosis<br>[95% CI] | 0.6%<br>[0.01%, 3.09%]  | 1.4%<br>[0.37%, 3.46%]   | 0.65    |

Note: Components (Cardiac death, MI and ID-TLR) are presented as non-hierarchical

Discussion: The comparability of Absorb BVS to XIENCE V demonstrated with in-scaffold / stent late loss at 6 months (Cohort B only), and clinical outcomes at 6 months and 1 year were maintained through 2-years post implantation. The collective data confirms the comparability between Absorb BVS and XIENCE V based on propensity analyses.

8.6 Benefit of the Absorb Bioresorbable Scaffold Technology

In the ABSORB Cohort A clinical investigation, Absorb BVS showed excellent long-term clinical outcomes with low MACE rates out to 5 years with absence of any cardiac death, ID-TLR, Q-wave MI, and scaffold thrombosis. The ABSORB Cohort B clinical investigation confirms these findings out to 3 years, including patency comparable to XIENCE V. Similarly, ABSORB EXTEND clinical results out to 2 years support Cohort A and B outcomes. Absorb BVS thus performs all the functions of a drug-eluting stent while offering future potential benefits resulting from the absence of a permanent metallic implant. The gradual disappearance of the scaffold allows for the possibility of vascular healing and restoration of vessel function as discussed in the Vasomotor Function Results at 2 years (sections 8.2.4 – 8.2.5). The absence of a permanent metallic implant in the vascular tissue may facilitate any required re-interventions on the target vessel / lesion or side branches either by percutaneous or surgical means, thus enabling a broader range of treatment options after bioresorption of the scaffold. Finally, unlike permanent metal implants, polymeric implants do not cause imaging artifacts during non-invasive CT or MR evaluation. This provides the additional benefit that a polymeric bioresorbable scaffold may be more compatible with the growing usage of non-invasive follow-up imaging than is the case with metallic stents, potentially facilitating patient management and providing economic benefits.

9.0 PATIENT SELECTION AND TREATMENT

9.1 Individualization of Treatment

The risks and benefits described above should be considered for each patient before using the Absorb BVS System. Patient selection factors to be assessed should include judgment regarding risk of antiplatelet therapy. Special consideration should be given to those patients with recently active gastritis or peptic ulcer disease.

Antiplatelet drugs should be used in combination with the Absorb BVS. Physicians should use the information from the SPIRIT clinical trials and the ABSORB trials, coupled with the current literature on drug-eluting stents / scaffolds and the specific needs of individual patients, to determine the specific antiplatelet / anticoagulation dose and duration to be used for their patients in general practice. A minimum of six months' duration is recommended.

It is very important that the patient comply with the postprocedural antiplatelet recommendations. Premature discontinuation of prescribed antiplatelet medication could result in a higher risk of thrombosis, MI, or death. Prior to PCI, if a surgical or dental procedure is anticipated that requires early discontinuation or temporary interruption of antiplatelet therapy, the interventionalist and patient should carefully consider whether an everolimus-eluting scaffold and its associated recommended antiplatelet therapy are the appropriate PCI choice. Following PCI, should a surgical or dental procedure be recommended, the risks and benefits of the procedure should be weighed against the possible risk associated with premature or temporary discontinuation of antiplatelet therapy. Patients who require premature discontinuation of antiplatelet therapy secondary to significant active bleeding should be monitored carefully for cardiac events, and once stabilized, have their antiplatelet therapy restarted as soon as possible per the discretion of their treating physicians.

10.0 CLINICIAN USE INFORMATION

10.1 Inspection Prior to Use

Prior to using the Absorb BVS System, carefully remove the system from the package and inspect for bends, kinks and other damage. Verify that the scaffold does not extend beyond the radiopaque balloon markers and that it is still well-crimped onto the balloon catheter. Do not use if any defects are noted.

10.2 Materials Required

- 6F/0.070"/1.8 mm minimum inner diameter guiding catheter(s) of appropriate shape

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### 10.6 Deployment Procedure

- CAUTION: Refer to product label for *in vitro* scaffold inner diameter, nominal pressure and RBP.**
- Prior to deployment, reconfirm the correct position of the scaffold relative to the target lesion using the radiopaque balloon markers.
- Deploy the scaffold slowly, by pressurizing delivery system in **2-atm increments, every 5 seconds**, until scaffold is completely expanded. Maintain pressure for 30 seconds. Fully expand the scaffold by inflating to nominal pressure at a minimum. Accepted practice generally targets an initial deployment pressure that would achieve a scaffold inner ratio of about 1.1 times the reference vessel diameter. **CAUTION: Do not exceed the labeled rated burst pressure RBP of 16 atm (1621 kPa) or maximum deployment diameter of the scaffold.**
- Fluoroscopic visualization during scaffold expansion should be used in order to properly judge the optimum scaffold diameter as compared to the proximal and distal native coronary artery diameters (reference vessel diameters). Optimal scaffold expansion and proper apposition require that the scaffold be in full contact with the arterial wall.
- If necessary, the delivery system can be repressurized or further pressurized to ensure complete apposition of the scaffold to the artery wall. Fully cover the entire lesion and balloon-treated area (including dissections) with the Absorb scaffold, allowing for adequate scaffold coverage into healthy tissue proximal and distal to the lesion.
- Deflate the balloon by pulling negative on the inflation device for 30 seconds. Confirm complete balloon deflation before attempting to move the delivery system. If unusual resistance is felt during scaffold delivery system withdrawal, pay particular attention to the guiding catheter position. **Note:** See section **10.8 Removal Procedure** for instruction on withdrawal of scaffold delivery system.
- Confirm scaffold position and deployment using standard angiographic techniques. For optimal results, the entire stenosed artery segment should be covered by the scaffold. Fluoroscopic visualization during scaffold expansion should be used in order to properly judge the optimum expanded scaffold diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the scaffold be in full contact with the artery wall, which can be facilitated with the use of routine angiography, intravascular ultrasound (IVUS), or optical coherence tomography (OCT).
- Postdilatation with a noncompliant balloon is recommended** following instructions in section **10.7**, as long as the postdilated segment is within the allowable expansion limits of the scaffold.

### 10.7 FURTHER EXPANSION OF THE DEPLOYED SCAFFOLD

- DEPLOYED SCAFFOLDS SHOULD NOT BE LEFT UNDERDILATED.** Deployed scaffolds should be well-apposed to the vessel wall. To achieve optimal scaffold apposition, postdilatation is recommended with a low profile, high pressure, noncompliant balloon, as long as the postdilated segment is within the allowable expansion limits of the scaffold. The scaffolded segment should be carefully recrossed with a prolapsed guide wire to avoid disrupting the scaffold geometry. Postdilatation must only be done with balloons sized to fit within the boundaries of the scaffold. **The compliance chart of the noncompliant balloon selected must be carefully reviewed prior to dilatation and an appropriate maximum pressure used to ensure that the scaffold is not over-dilated.** The scaffold dilatation limit is 0.5 mm above the nominal diameter for all scaffold sizes. **CAUTION: Do not dilate the scaffold beyond the dilatation limit which is 0.5 mm above the nominal diameter. Over-dilatation may result in scaffold damage.**

| Nominal Scaffold Diameter | Dilatation Limit                           |
|---------------------------|--------------------------------------------|
| 2.5 mm                    | 3.00 mm<br>Maximum postdilatation diameter |
| 3.0 mm                    | 3.50 mm<br>Maximum postdilatation diameter |
| 3.5 mm                    | 4.00 mm<br>Maximum postdilatation diameter |
- If more than one Absorb BVS is needed to cover the lesion and balloon-treated area, it is suggested that, to avoid the potential for gap restenosis, the scaffolds be overlapped by a minimum of 1 mm and a maximum of 4 mm. To ensure that there are no gaps between scaffolds, the balloon marker bands of the second Absorb BVS should be positioned inside the deployed scaffold prior to expansion. It is recommended not to use more than two Absorb BVS to treat one lesion.
- Ensure the final scaffold diameter matches the reference vessel diameter to **ENSURE GOOD SCAFFOLD APPPOSITION**. Reconfirm scaffold position and angiographic results. Repeat inflations until achieved.

### 10.8 Removal Procedure

**Withdrawal of the scaffold delivery catheter / postdilatation balloon from the deployed scaffold:**

- Deflate the balloon by pulling negative pressure on the inflation device. Larger and longer balloons will take more time (up to 30 seconds) to deflate than smaller and shorter balloons. Confirm balloon deflation under fluoroscopy.
  - Position the inflation device to "negative" or "neutral" pressure.
  - Fully open the rotating hemostatic valve.
  - Stabilize guide catheter position just outside coronary ostium and anchor in place. Maintain guide wire placement across scaffolded segment.
  - Gently remove the scaffold delivery system with slow and steady pressure.
  - Tighten the rotating hemostatic valve.
- If, during withdrawal of the catheter from the deployed scaffold, resistance is encountered, use the following steps to improve balloon rewrap:**
- Re-inflate the balloon up to nominal pressure, deflate and change pressure to neutral.
  - Repeat steps 1 – 5 above.
  - Re-evaluate the scaffolded region once the balloon is removed to assure optimal apposition.
- Note:** After successful withdrawal of the balloon from the deployed scaffold, should any resistance be felt at any time when withdrawing the scaffold delivery system or post-dilatation balloon into the guide catheter, **remove the entire system as a single unit**. See section **6.4 Scaffold / System Removal – Precautions** for specific delivery system removal instructions.

### 11.0 TRADEMARKS

Absorb, XIENCE PRIME, and XIENCE V are trademarks of the Abbott Group of Companies. Certican and Aflinor are registered trademarks of Novartis Pharmaceuticals, Inc. CYPHER is a registered trademark of Cordis Corporation in the U.S.

- for the target vessel
- 2 – 3 syringes (10 – 20 cc)
  - 1,000 u/500 cc heparinized normal saline (HepNS)
  - 0.014 inch (0.36 mm) x 175 cm (minimum length) guide wire
  - Rotating hemostatic valve with 0.096 inch (2.44 mm) minimum inner diameter
  - Contrast diluted 1:1 with normal saline
  - Inflation device
  - Three-way stopcock
  - Torque device
  - Guide wire introducer

### 10.3 Vessel and Lesion Selection

- Quantitative imaging is recommended for the assessment of target vessel diameter at baseline for appropriate Absorb BVS size selection.
- The target vessel diameter ranges to be treated in the procedure are indicated in **Table 12**, along with the Absorb BVS diameter to be used.

**Table 12: Target Vessel Diameter Ranges and Absorb BVS Diameter to be Used**

| Target Vessel Diameter Distal and Proximal | Absorb BVS Diameter to be Used |
|--------------------------------------------|--------------------------------|
| ≥ 2.0 mm and ≤ 3.0 mm                      | 2.5 mm                         |
| ≥ 2.5 mm and ≤ 3.3 mm                      | 3.0 mm                         |
| ≥ 3.0 mm and ≤ 3.8 mm                      | 3.5 mm                         |

- For cases where the combination of target vessel diameter and target lesion length is appropriate to be treated by more than one scaffold size, the selection of scaffold size can be made per the judgment of the physician.

### 10.4 Preparation

#### 10.4.1 Packaging Removal

**Note:** The foil pouch is the sterile barrier. Sterile product is contained within this one pouch—there is no secondary pouch.

- Peel the pouch open from the top corner.
- Carefully remove the delivery system from its protective tubing for preparation of the delivery system.
- Do not bend or kink the hypotube during removal.

#### 10.4.2 Dual Layer Sheath Removal

- While holding the distal catheter shaft with one hand, grasp only the yellow outer sheath with the other hand and gently slide the sheath distally.
- A longitudinal split on the inner sheath will open up and be visible.
- The stylet and dual layer sheath are removed simultaneously from the delivery system by continuing to slide the yellow sheath distally until the inner and outer layers of the dual layer sheath as well as the stylet are free from the catheter system. See Section **6.1 Scaffold Handling – Precautions**. Do not use the device if the sheath cannot be removed as indicated.
- Verify that the scaffold does not extend beyond the radiopaque balloon markers and no scaffold struts are lifted. **Do not use if any defects are noted.**

#### 10.4.3 Guide Wire Lumen Flush

- Flush the guide wire lumen with HepNS until fluid exits the guide wire exit notch. **Note:** Avoid manipulation of the scaffold while flushing the guide wire lumen as this may disrupt the placement of the scaffold on the balloon.

#### 10.4.4 Delivery System Preparation

- Prepare an inflation device / syringe with diluted contrast medium.
- Attach an inflation device / syringe to stopcock; attach it to the inflation port of the product. Do not bend the product hypotube when connecting to the inflation device / syringe.
- With the tip down, orient the delivery system vertically.
- Open the stopcock to delivery system; pull negative for 30 seconds; release to neutral for contrast fill.
- Close the stopcock to the delivery system; purge the inflation device / syringe of all air.
- Repeat steps 3 through 5 until all air is expelled. If bubbles persist, do not use the product.
- If a syringe was used, attach a prepared inflation device to the stopcock.
- Open the stopcock to the delivery system.
- Leave on neutral.

**Note:** The labeled scaffold diameter refers to expanded scaffold inner diameter.

### 10.5 Delivery Procedure

- Prepare the vascular access site according to standard practice.
- Predilate the lesion with a percutaneous transluminal coronary angioplasty catheter.** Predilatation is strongly recommended and can also be utilized to properly size the vessel. **Note:** Limit the length of the predilatation by the PTCA balloon to avoid creating a region of vessel injury that is outside the boundaries of the Absorb scaffold.
- Administer a standard dose of intracoronary nitroglycerine prior to finalizing the RVD within the target zone.
- Maintain neutral pressure on the inflation device. Open the rotating hemostatic valve as widely as possible.
- Backload the delivery system onto the proximal portion of the guide wire, while maintaining guide wire position across the target lesion.
- Advance the delivery system over the guide wire to the target lesion. Utilize radiopaque balloon markers to position the scaffold across the lesion; perform angiography to confirm scaffold position. **Note:** If removal of a scaffold system is required prior to deployment, ensure that the guide catheter is coaxially positioned relative to the scaffold delivery system and cautiously withdraw the scaffold delivery system into the guiding catheter. Should **unusual resistance** be felt at any time when withdrawing the scaffold into the guide catheter, the scaffold delivery system and the guide catheter should be **removed as a single unit**. This should be done under direct visualization with fluoroscopy.
- Tighten the rotating hemostatic valve. The scaffold is now ready to be deployed.

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<sup>5</sup> Claessen BE, Beijk MA, Legrand V, Ruzlylo W, Manari A, Varenne O, Sutrop MJ, Tijssen JG, Miquel-Hebert K, Velthof S, Henriques JP, Serruys PW, Piek JJ. Two-year clinical, angiographic, and intravascular ultrasound follow-up of the XIENCE V everolimus-eluting stent in the treatment of patients with *de novo* native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009 Aug;2(4):339-47.

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<sup>5</sup> Claessen BE, Beijk MA, Legrand V, Ruziyllo W, Manari A, Varenne O, Suttrop MJ, Tijssen JG, Miquel-Hebert K, Veldhof S, Henriques JP, Serruys PW, Piek JJ. Two-year clinical, angiographic, and intravascular ultrasound follow-up of the XIENCE V everolimus-eluting stent in the treatment of patients with *de novo* native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009 Aug;2(4):339-47.

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<sup>5</sup> Claessen BE, Beijk MA, Legrand V, Ruzulyo V, Manari A, Varenne O, Suttorp MJ, Tijssen JG, Miquel-Hebert K, Velthof S, Henriques JP, Serruys PW, Piek JJ. Two-year clinical, angiographic, and intravascular ultrasound follow-up of the XIENCE V everolimus-eluting stent in the treatment of patients with *de novo* native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009; 24(4):339-47.

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<sup>5</sup> Claessen BE, Beijk MA, Legrand V, Ruzlylo V, Manari A, Varenne O, Sutrop MJ, Tijssen JG, Miquel-Hebert K, Veldhof S, Henriques JP, Serruys PW, Piek JC. Two-year clinical, angiographic, and intravascular ultrasound follow-up of the XIENCE V everolimus-eluting stent in the treatment of patients with *de novo* native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009 Aug;2(4):339-47.

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<sup>5</sup> Claessen BE, Beijk MA, Legrand V, Ruzlylo V, Manari A, Varenne O, Suttorp MJ, Tijssen JG, Miquel-Hebert K, Velthof S, Henriques JP, Serruys PW, Piek JX. Two-year clinical, angiographic, and intravascular ultrasound follow-up of the XIENCE V everolimus-eluting stent in the treatment of patients with de novo native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009 Aug;2(4):339-47.

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<sup>5</sup> Claessen BE, Beijk MA, Legrand V, Ruzulyo W, Manari A, Varenke O, Suttorp MJ, Tjissen JG, Michel-Hebert K, Velthof S, Henriques JP, Serruys PW, Piek JJ. Two-year clinical, angiographic, and intravascular ultrasound follow-up of the XIENCE V everolimus-eluting stent in the treatment of patients with de novo native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009 Aug;2(4):339-47.

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<sup>8</sup> Claessen BE, Beijk MA, Legrand V, Ruzyllo W, Manari A, Varenne O, Suttorp MJ, Tijssen JG, Miquel-Hebert K, Velthof S, Henriques JP, Serruys PW, Piek JJ. Two-year clinical, angiographic, and intravascular ultrasound follow-up of the XIENCE V everolimus-eluting stent in the treatment of patients with de novo native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009 Aug;2(4):339-47.

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<sup>5</sup> Claessen BE, Beijl MA, Legrand V, Ruzlylo W, Manari A, Varenne O, Suttrop MJ, Tijssen JG, Michel-Hebert K, Velthof S, Henriques JP, Serruys PW, Piek JJ. Two-year clinical, angiographic and intracatheter follow-up of the XIENCE V evolutionary-eluting stent in the treatment of patients with *de novo* native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009 Aug;2(4):339-47.

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<sup>5</sup> Claessen BE, Beijik MA, Legrand V, Ruzlylo V, Manari A, Varenne O, Suttorp MJ, Tijssen JG, Miquel-Hebert K, Veldhof S, Henriques JP Serruys PW, Piek JF. Two-year clinical, angiographic, and intravascular ultrasound follow-up of the XIENCE V everolimus-eluting stent in the treatment of patients with *de novo* native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009 Aug;2(4):339-47.

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\* Claessen BE, Beijk MA, Legrand V, Ruzlyo W, Manari A, Varenne O, Suttorp MJ, Tijssen JG, Miquel-Hebert K, Velthof S, Henriques JP, Serruys PW, Piek JJ. Two-year clinical, angiographic, and intravascular ultrasound follow-up of the XIENCE V everolimus-eluting stent in the treatment of patients with *de novo* native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009 Aug;2(4):339-47.

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<sup>5</sup> Claessen BE, Beijl MA, Legrand V, Ruzlylo W, Manari A, Varenne O, Suttrop MJ, Tijssen JG, Miquel-Hebert K, Veldhof S, Henriques JP, Serruys PW, Piek JJ. Two-year clinical, angiographic, and intravascular ultrasound follow-up of the XIENCE V everollimust-eluting stent in the treatment of patients with de novo native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009 Aug;2(4):339-47.

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<sup>5</sup> Claessen BE, Beijk MA, Legrand V, Ruzlylo V, Manari A, Varenne O, Suttrop MJ, Tijssen JG, Miguel-Hebert K, Velthof S, Henriques JP, Serruys PW, Piek JH. Two-year clinical, angiographic, and intravascular ultrasound follow-up of the XIENCE V evelurom-eluting stent in the treatment of patients with *de novo* native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009 Aug;2(4):339-47.

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
Rapid Exchange; Greitisaitis keitimas; Hitra izmenjava

Non-pyrogenic; Nepirogeninis; Apirogeno

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### INFORMATION FOR PRESCRIBERS

#### Table of Contents

|        |                                                                                                                               |
|--------|-------------------------------------------------------------------------------------------------------------------------------|
| 1.0    | DEVICE DESCRIPTION                                                                                                            |
|        | Table 1: Drug Content in Absorb BVS                                                                                           |
|        | Table 2: <i>In vitro</i> Device Specifications                                                                                |
|        | Table 3: TagAlert Window Indicator Options, Instructions for Use and Disposal                                                 |
| 2.0    | HOW SUPPLIED                                                                                                                  |
| 3.0    | INDICATIONS                                                                                                                   |
| 4.0    | CONTRAINDICATIONS                                                                                                             |
| 5.0    | WARNINGS                                                                                                                      |
| 6.0    | PRECAUTIONS                                                                                                                   |
| 6.1    | Scaffold Handling – Precautions                                                                                               |
| 6.2    | Scaffold Placement – Precautions                                                                                              |
| 6.3    | Use in Conjunction with Other Procedures                                                                                      |
| 6.4    | Scaffold / System Removal – Precautions                                                                                       |
| 6.5    | Post Implant – Precautions                                                                                                    |
| 6.6    | MRI Statement                                                                                                                 |
| 6.7    | Drug Interactions                                                                                                             |
| 6.8    | Pregnancy                                                                                                                     |
| 7.0    | ADVERSE EVENTS                                                                                                                |
| 7.1    | Observed Adverse Events                                                                                                       |
| 7.2    | Potential Adverse Events                                                                                                      |
| 8.0    | CLINICAL INVESTIGATIONS OF ABSORB BVS SYSTEM                                                                                  |
| 8.1    | ABSORB Clinical Trial                                                                                                         |
| 8.2    | ABSORB Cohort A                                                                                                               |
|        | 8.2.1 Methodology and Current Status                                                                                          |
|        | 8.2.2 Clinical Outcome at 5 Years                                                                                             |
|        | Table 4: Clinical Outcome at 5 Years                                                                                          |
|        | 8.2.3 Angiographic, IVUS, and OCT Outcomes at 180 Days and 2 Years                                                            |
|        | 8.2.4 Vasomotor Function Results at 2 Years                                                                                   |
|        | Figure 1: Results of Acetylcholine and Methergine testing in proximal, scaffolded, and distal segments                        |
|        | 8.2.5 Discussion                                                                                                              |
| 8.3    | ABSORB Cohort B                                                                                                               |
|        | 8.3.1 Methodology and Current Status                                                                                          |
|        | 8.3.2 Clinical Outcomes for Full Cohort B                                                                                     |
|        | Table 5: Hierarchical Clinical Outcomes for Cohort B (ITT Population)                                                         |
|        | 8.3.3 Angiographic, IVUS, and OCT Outcomes at 180 Days, 1, 2, and 3 Years                                                     |
|        | 8.3.4 Discussion                                                                                                              |
| 8.4    | ABSORB EXTEND                                                                                                                 |
|        | 8.4.1 Methodology and Current Status                                                                                          |
|        | 8.4.2 Summary of Clinical Results out to 1 year (Interim Data Cutoff)                                                         |
|        | Table 6: ABSORB EXTEND Subject Counts of Ischemia-driven Adverse Events through 393 Days (All Subjects Registered Population) |
|        | 8.4.3 Summary of Clinical Results out to 2 years (Interim Data Cutoff)                                                        |
|        | Table 7: ABSORB EXTEND Subject Counts of Ischemia-driven Adverse Events through 758 Days (250 Subjects Registered Population) |
| 8.5    | Propensity Score Adjusted Analysis of Absorb BVS System and XIENCE V                                                          |
| 8.5.1  | 6-month Propensity Analysis                                                                                                   |
|        | Table 8: MACE at 6 Months from the IPSW-Adjusted Data                                                                         |
|        | Table 9: Angiographic Outcomes at 6 Months from the IPSW-Adjusted Data                                                        |
| 8.5.2  | 1-year Propensity Analysis                                                                                                    |
|        | Table 10: Clinical outcomes at 1 year from the IPSW-Adjusted Cohorts                                                          |
| 8.5.3  | 2-year Propensity Analysis                                                                                                    |
|        | Table 11: Clinical outcomes at 2 years from the IPSW-Adjusted Cohorts                                                         |
| 8.6    | Benefit of the Absorb Bioresorbable Scaffold Technology                                                                       |
| 9.0    | PATIENT SELECTION AND TREATMENT                                                                                               |
| 9.1    | Individualization of Treatment                                                                                                |
| 10.0   | CLINICIAN USE INFORMATION                                                                                                     |
| 10.1   | Inspection Prior to Use                                                                                                       |
| 10.2   | Materials Required                                                                                                            |
| 10.3   | Vessel and Lesion Selection                                                                                                   |
|        | Table 12: Target Vessel Diameter Ranges and Absorb BVS Diameter to be Used                                                    |
| 10.4   | Preparation                                                                                                                   |
| 10.4.1 | Packaging Removal                                                                                                             |
| 10.4.2 | Dual Layer Sheath Removal                                                                                                     |
| 10.4.3 | Guide Wire Lumen Flush                                                                                                        |
| 10.4.4 | Delivery System Preparation                                                                                                   |
| 10.5   | Delivery Procedure                                                                                                            |
| 10.6   | Deployment Procedure                                                                                                          |
| 10.7   | Further Expansion of the Deployed Scaffold                                                                                    |
| 10.8   | Removal Procedure                                                                                                             |
| 11.0   | TRADEMARKS                                                                                                                    |

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- A pre-mounted polymer poly (L-lactide) (PLLA) scaffold coated with a blend of the antiproliferative drug everolimus and polymer poly (D,L-lactide) (PDLLA) in a 1:1 ratio. The available dose of everolimus on the scaffold is shown in **Table 1**.

Table 1: Drug Content in Absorb BVS

| Scaffold Diameter (mm) | Scaffold Length (mm) | Drug Dose (µg) |
|------------------------|----------------------|----------------|
| 2.5, 3.0               | 8                    | 76             |
| 2.5, 3.0               | 12                   | 114            |
| 2.5, 3.0               | 18                   | 181            |
| 2.5, 3.0               | 23                   | 228            |
| 2.5, 3.0               | 28                   | 276            |
| 3.5                    | 12                   | 135            |
| 3.5                    | 18                   | 197            |

| Scaffold Diameter (mm) | Scaffold Length (mm) | Drug Dose (µg) |
|------------------------|----------------------|----------------|
| 3.5                    | 23                   | 246            |
| 3.5                    | 28                   | 308            |

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\*\* Assure full deployment of the scaffold (see section **10.6 Deployment Procedure**). Deployment pressures should be based on lesion characteristics.

- TagAlert is a non-sterile temperature monitor for the shipping and storage of the Absorb BVS System. Before use of this product, check the TagAlert indicator located through the window in the back of the product box. The indicator should only show an OK as indicated in the digital display (**Table 3A**). If any other screen is present (**Table 3B**), do not use product.

Table 3: TagAlert Window Indicator Options, Instructions for Use and Disposal

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### 5.0 WARNINGS

- The potential long-term benefits of the Absorb coronary device in comparison to the available therapies are to be confirmed in additional clinical studies.
- This product should not be used in patients who are not likely to comply with the recommended antiplatelet therapy.
- Judicious selection of patients is necessary, since the use of this device carries the associated risk of scaffold thrombosis, vascular complications, and / or bleeding events.
- Oral administration of everolimus in combination with cyclosporine has been associated with increased serum cholesterol and triglycerides. Therefore, patients should be monitored for changes in lipid profiles.
- Persons allergic to poly (L-lactide), poly (D,L-lactide), everolimus, or platinum may suffer an allergic reaction to this implant.
- It is not recommended to treat patients having a lesion with excessive tortuosity proximal to or within the lesion.
- Devices (i.e., guide sheaths) that decrease the inner diameter of the guide catheter outside of the Absorb BVS System minimum guide catheter compatibility (**Table 2**) must not be used with the Absorb BVS System. Do not insert a 5-in-6, or a 6-in-7 guide sheath into a 6F or 7F guiding catheter, as doing so will result in an inner diameter that is too small for use with the Absorb BVS System.
- Balloon dilatation of any cells of a deployed Absorb BVS may cause scaffold damage.
- Careful selection of the target lesion reference vessel diameter to the scaffold diameter, and adequate lesion preparation prior to scaffold implantation is recommended to minimize potential damage to the scaffold during placement. It is not recommended to treat patients having a lesion that prevents complete inflation of an angioplasty balloon (example: a severely calcified lesion that has not had adequate lesion preparation), or a lesion with greater than 40% residual stenosis after predilatation by visual estimation.

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6.2 Scaffold Placement – Precautions

- Do not prepare or preinflate the delivery system prior to scaffold deployment, other than as directed. Use balloon pumping technique described in section 10.4.4 Delivery System Preparation.
- Size the reference target lesion diameter appropriately to avoid overexpanding the scaffold to ensure adequate scaffold apposition. This will reduce the risk of causing scaffold damage.
- When introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the scaffold from the balloon.
- Do not torque the catheter more than one (1) full turn.
- Use caution when advancing the Absorb BVS across the lesion. Multiple attempts to cross a lesion may lead to scaffold damage or dislodgement.
- Implanting a scaffold may lead to dissection of the vessel distal and / or proximal to the scaffold and may cause acute closure of the vessel, requiring additional intervention (CABG, further dilatation, placement of additional scaffolds, or other).
- Do not expand the scaffold if it is not properly positioned in the vessel. (See section 6.4 Scaffold / System Removal – Precautions.)
- When treating multiple lesions within the same vessel, scaffold / stent the distal lesion prior to scaffolding the proximal lesion. Scaffolding / stenting in this order obviates the need to cross the proximal scaffold in placement of the distal scaffold / stent, and reduces the chance of damaging or dislodging the proximal scaffold.
- Avoid scaffolding across any side branches  $\geq 2.0$  mm in diameter:
- The inflated balloon diameter of the system used to deploy the scaffold should approximate the diameter of the vessel. To ensure full expansion of the scaffold, the balloon should be inflated to a minimum of nominal pressure.
- Do not exceed the Rated Burst Pressure (RBP) as indicated on the product label. Monitor balloon pressures during inflation. Use of pressures higher than specified on the product label may result in a ruptured balloon, with possible intimal damage and dissection.
- Scaffold retrieval methods (use of additional wires, snares, and / or forceps) may result in additional trauma to the coronary vasculature and / or the vascular access site. Complications may include bleeding, hematoma, or pseudoaneurysm.
- When multiple Absorb Bioresorbable Vascular Scaffolds and drug-eluting stents are required, only Absorb Bioresorbable Vascular Scaffolds or everolimus-eluting stents must be used. Potential interaction with other drug-eluting stents or coated stents has not been evaluated and should be avoided.
- The extent of the patient's exposure to drug and polymer is directly related to the number of scaffolds implanted. A patient can receive up to a total scaffolding length of 94 mm.
- The safety and effectiveness of the Absorb BVS in patients with prior brachytherapy of the target lesion or the use of brachytherapy for treated site restenosis in an Absorb BVS have not been established. Both vascular brachytherapy and the Absorb BVS alter arterial remodeling. The combination between these two treatments has not been established.
- Postdilatation with a noncompliant balloon is recommended following instructions in section 10.7, as long as the postdilated segment is within the allowable expansion limits of the scaffold.
- An unexpanded scaffold may be retracted into the guiding catheter **one time only**. An unexpanded scaffold should not be reintroduced into the artery once it has been pulled back into the guiding catheter. Subsequent movement in and out through the distal end of the guiding catheter should not be performed, as the scaffold may be damaged or dislodged during retraction back into the guiding catheter.
- Should resistance be felt at any time during removal of the undeployed Absorb BVS System, please refer to the steps provided in section 6.4 Scaffold / System Removal – Precautions.
- In the rare event of acute occlusion following scaffold placement, a bailout implant may be inserted and deployed within the scaffold such that the Absorb BVS is completely covered by the bailout implant. All abrupt closures must be treated as an emergency per the hospital standard of care.  
**Note:** It is recommended that bailouts be done with a metallic everolimus-eluting stent of appropriate size.

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- Scaffold delivery system removal prior to scaffold deployment:  
If removal of a scaffold system is required prior to deployment, ensure that the guide catheter is coaxially positioned relative to the scaffold delivery system and cautiously withdraw the scaffold delivery system into the guiding catheter. Should **unusual resistance** be felt at any time when withdrawing the scaffold into the guide catheter, the scaffold delivery system and the guide catheter should be **removed as a single unit**. This should be done under direct visualization with fluoroscopy.
- Withdrawal of the scaffold delivery system / postdilatation balloon from the deployed scaffold:
  1. Deflate the balloon by pulling negative on the inflation device. Larger and longer balloons will take more time (up to 30 seconds) to deflate than smaller and shorter balloons. Confirm balloon deflation under fluoroscopy.
  2. Position the inflation device to "negative" or "neutral" pressure.
  3. Stabilize guide catheter position just outside coronary ostium and anchor in place. Maintain guide wire placement across scaffold segment.
  4. Gently remove the scaffold delivery system with slow and steady pressure.
  5. Tighten the rotating hemostatic valve.**Notes:**
  - 1) If, during withdrawal of the catheter from the deployed scaffold, resistance is encountered, use the following steps to improve balloon rewrap:
    - Re-inflate the balloon up to nominal pressure, deflate and change pressure to neutral.
    - Repeat steps 1 through 5 above.
  - 2) After successful withdrawal of the balloon from the deployed scaffold, should any resistance be felt at any time when withdrawing the scaffold delivery system or postdilatation balloon into the guide catheter, **remove the entire system as a single unit**.

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**Treatment Strategy:** Predilatation of the target lesion was mandatory. Planned overlapping of Absorb BVS was disallowed. Any bailout must be done with overlapping a XIENCE V stent of appropriate length, and if not available, a CYPHER<sup>®</sup> sirolimus-eluting stent; bailout with Absorb BVS was not permitted. Postdilatation was performed at operator discretion, but only using balloons sized to fit within the boundaries of the scaffold.

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The ABSORB Cohort A Trial implanted the Absorb Cohort A BVS in patients with a single *de novo* native coronary artery lesion. Enrollment started with the 3.0 x 12 mm size. The 3.0 x 18 mm size became available later during enrollment and was used in only 2 patients.

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<sup>3</sup> Certican<sup>®</sup> drug prescribing sheet, date of most recent amendment May 21, 2013. Novartis.

<sup>4</sup> Aflinor<sup>®</sup> summary of product characteristics, date of authorization November 2013. Novartis.

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QCA results were collected from baseline, post procedure, 180 days and 2 years. The 180-day data demonstrated an acceptable in-scaffold late loss of 0.43 ± 0.37 mm<sup>2</sup> possibly driven by bioactive remodeling or mechanical late recoil. From 180 days to 2 years, the change in in-scaffold late loss was very low (0.48 ± 0.28 mm<sup>2</sup>).

The 180-day grey-scale intravascular ultrasound (IVUS) analysis showed a significant reduction in the average lumen area (6.04 ± 1.12 mm<sup>2</sup> after procedure vs. 5.19 ± 1.33 mm<sup>2</sup> at 180-day, *p* < 0.001). The vessel area remained constant between baseline and 180-day follow-up (13.49 ± 3.74 mm<sup>2</sup> vs. 13.79 ± 3.84 mm<sup>2</sup>), demonstrating the absence of significant expansive or constrictive remodeling. At 2 years, the main observation provided by grey-scale IVUS was the increase in minimal luminal area and average luminal area / volume together with a significant decrease in plaque area / volume between 6 months and 2 years.

The 180-day OCT group (n = 13) provided a total of 671 apparent struts for visual evaluation, which showed that 93% of struts evaluated were well apposed to the vessel wall, and 99% of the struts were covered by tissue. The serial OCT group (n = 7) had serial data postprocedure, at 180 days and at 2 years (intention-to-treat population). The number of apparent struts decreased from 403 at baseline to 368 at 180 days and to 264 at 2-years follow-up (34.5% reduction over two years), all well covered and apposed to the vessel wall. Precinical analysis directly comparing pre-explanation OCT to post-explanation histology indicates there is no correlation between the presence or absence of OCT-visible features (apparent struts) and the presence or absence of polylactide polymer in the vessel. Hence, visual OCT features (apparent struts) are not always indicative of lack of resorption, but absence of the apparent struts by OCT confirms resorption. The lumen shape was regular with smooth, well delineated borders in all cases, and no intraluminal tissue was observed. Importantly, minimal and mean luminal area decreased significantly between post procedure and 180 days, but numerically enlarged between 180 days and 2 years.

8.2.4 Vasmotor Function Results at 2 Years

Vasomotor function proximal, within, and distal to the treated (scaffolded) segments at 2 years, was evaluated with either the endothelium-independent vasoconstrictive methylergonovine maleate (Methergine, Novartis, Basel, Switzerland), or the endothelium-dependent vasoactive agent acetylcholine (Ovisot, Daiichi-Sankyo, Tokyo, Japan), depending on local practice.

In the Methergine group (n = 7), significant vasoconstriction was observed in proximal and scaffolded segments (Figure 1). In the acetylcholine group (n = 9), although the overall change in vessel dimension was not statistically significant, vasodilatation was observed in five patients. Nitrates induced a significant vasodilatation in the scaffolded and distal segments, relative to baseline (pre-acetylcholine) at 2 years. These results suggest the potential restoration of vasomotor function in the treated segments when the Absorb Cohort A BVS had resorbed. Further data on vascular responses will be collected from ongoing clinical trials.

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The 180-day angiographic results from Cohort B Group 1 (Cohort B1) demonstrated a late loss of 0.19 mm (N = 42) which compares well to the 0.10 mm late loss (N = 23) of the 3.0 x 18 mm XIENCE V Everolimus Eluting Coronary Stent System (XIENCE V EECSS) observed in the SPIRIT FIRST First-In-Man trial, and show favorable improvement from the 0.43 mm late loss (N = 26) from Cohort A.

At 180 days, the IVUS results from Cohort B1 showed limited intrascaffold neointimal hyperplasia, the volume obstruction (VO) was 1.2% (N = 40) and the neointimal hyperplasia area was 0.08 mm<sup>2</sup> (N = 40). These results compare favorably with XIENCE V EECSS in SPIRIT FIRST (VO of 8.0% [N = 21] and NIH area of 0.56 mm<sup>2</sup> [N = 21]) and with Cohort A (VO of 5.3% [N = 25] and NIH area of 0.29 mm<sup>2</sup> [N = 25]). As observed in Cohort A, the 6-month IVUS results showed a significant reduction in the average lumen area (6.60 ± 1.22 mm<sup>2</sup> [N = 37] after procedure vs. 6.37 ± 1.12 mm<sup>2</sup> [N = 37] at 180-day, *p* = 0.0049). The vessel area remained comparable between baseline and 180-day follow-up (14.22 ± 3.75 mm<sup>2</sup> [N = 37] vs. 14.49 ± 3.67 mm<sup>2</sup> [N = 37]) demonstrating the absence of significant expansive or constrictive remodeling.

The 1-year angiographic results from Cohort B Group 2 (Cohort B2) demonstrated a late loss of 0.27 mm (N = 56) which compares well to the 0.23 mm late loss (N = 22) of the 3.0 x 18 mm XIENCE V observed in the SPIRIT FIRST trial. At 1 year, the IVUS results from Cohort B2 (N = 54) showed limited intrascaffold neointimal hyperplasia, the percentage of in-scaffold area obstruction and the neointimal hyperplasia area were minimal (1.43% and 0.09 mm<sup>2</sup>, respectively). Over time there was a statistically significant but not clinically relevant increase in total plaque area with a concomitant increase in vessel area. Furthermore, average lumen area remained stable from post-procedure (6.31 ± 0.95 mm<sup>2</sup>) to 1 year (6.33 ± 1.17 mm<sup>2</sup>). By OCT, there was no significant change in the mean and minimum scaffold area between post procedure and 1 year (N = 21). There was a reduction in minimal lumen area and mean luminal area decreased significantly as a result of mean neointimal area growth (1.34 ± 0.67 mm<sup>2</sup>) but was not clinically significant. The lumen area stenosis increase was statistically significant from 20.2% post-procedure to 26.9% at 1 year but not clinically relevant. These data further supported that there was no reduction in scaffold area.

Imaging results at 2 years from Cohort B Group 1 (Cohort B1) demonstrated in-scaffold angiographic late loss of 0.27 mm (N = 38), which matches the late loss observed at 1 year for Cohort B Group B2 but is lower than the 0.33 mm (N = 83) reported from XIENCE V in SPIRIT II. The IVUS data at 2 years (N = 33) revealed a finding unique to Absorb BVS, which is mean scaffold area enlargement between 6 months and 2 years (6.42 versus 7.08 mm<sup>2</sup>, *p* < 0.0001). This enlargement was greater than the observed increase in neointimal hyperplasia and was accompanied by late increased mean lumen area (6.36 versus 6.85 mm<sup>2</sup>, *p* = 0.0105). IVUS findings were reported with comparative paired (serial) analysis at post-procedure, 6 months, and 2 years. Similarly, OCT data (N = 23) were reported with paired (serial) analysis, and showed that from 6 months to 2 years the mean strut core area decreased from 0.20 to 0.15 mm<sup>2</sup> (*p* < 0.0001), the mean neointimal area increased but remained minimal at 2 years (1.43 versus 2.11 mm<sup>2</sup>, *p* < 0.0001) and mean lumen area was unchanged. Moreover, the mean scaffold area increased from 7.47 to 8.24 mm<sup>2</sup> (*p* = 0.0155) from post-procedure to 2 years indicating loss of mechanical integrity of the scaffold and potential expansion of the vessel. Scaffold strut coverage by neointima was near completed at 6 months (98%) and 2 years (99%).

At 3 years post-procedure, IVUS results from Cohort B Group 2 (Cohort B2, N=44) revealed scaffold enlargement (average scaffold area of 6.35 ± 0.99 mm<sup>2</sup> at 1 year vs. 7.08 ± 1.55 mm<sup>2</sup> at 3 years, *p* < 0.0001) and increased average lumen area (6.35 ± 1.17 mm<sup>2</sup> at 1 year vs. 6.81 ± 1.62 mm<sup>2</sup> at 3 years, *p* = 0.0006). The late lumen enlargement is of particular interest because the vessel area was preserved (average vessel area of 14.43 ± 2.64 mm<sup>2</sup> at 1 year vs. 14.58 ± 2.67 mm<sup>2</sup> at 3 years, *p* = 0.407) despite an increase in neointimal area (0.08 ± 0.13 mm<sup>2</sup> at 1 year vs. 0.28 ± 0.41 mm<sup>2</sup> at 3 years, *p* < 0.0001). OCT results at 3 years (N=18) also revealed scaffold enlargement between 1 and 3 years (mean scaffold area of 7.51 ± 0.95 mm<sup>2</sup> at 1 year vs. 8.64 ± 2.15 mm<sup>2</sup> at 3 years, *p* = 0.0008). Lumen area remained relatively constant (6.01 ± 1.29 mm<sup>2</sup> at 1 year vs. 6.09 ± 1.67 mm<sup>2</sup> at 3 years, *p* = 0.679). The vessel area was not measured because OCT cannot demarcate the external elastic lamina. Overall, the neointimal hyperplasia remained small although there was an increase in neointimal hyperplasia between 1 year and 3 years (mean neointimal area of 1.41 ± 0.68 mm<sup>2</sup> at 1 year vs. 2.35 ± 0.68 mm<sup>2</sup> at 3 years, *p* < 0.0001).

8.3.4 Discussion

A total of 101 patients were enrolled with clinical data available out to 3 years for the full Cohort. In addition, imaging data at 180 days and 2 years are available for the 45 patients in Cohort B1, and at 1 year and 3 years for the 56 patients in Cohort B2. Overall, the Absorb BVS System has demonstrated acute success (98% procedure success and 100% device success) and safety up to 3 years after implantation. At 3 years, the MACE rate remained low at 3.4% and 10.0% in the ABSORB Cohort A and Cohort B groups, respectively. Furthermore, there were no cases of cardiac death or scaffold thrombosis in either cohort. Based on these results, the performance and safety continues to be established in the Cohort A and Cohort B studies.

Imaging of Cohort B indicated that angiographic late loss appears to stabilize between 1 and 3 years. In Group B1 and Group B2, IVUS results indicate a late expansion of scaffold and lumen dimensions between 6 months and 2 years, and between 1 and 3 years, respectively, which was consistent with the scaffold expansion observed by OCT. These results suggest scaffold bioresorption accompanied by a loss of structural integrity. This phenomenon is unique to bioresorbable vascular scaffolds and is in contrast with serial IVUS imaging of metallic DES, where late lumen expansion has not been observed<sup>5</sup>. Follow-up is ongoing. Patients will be followed through 5 years.

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4.4.2 Summary of Clinical Results out to 1 year (Interim Data Cutoff)

Of the total of 605 patients registered in ABSORB EXTEND as of December 30, 2012, 37-day (30 days  $\pm$  7-day window), 194-day (180 days  $\pm$  14-day window) and 393-day (365 days  $\pm$  28-day window) clinical follow-up data were available for 591, 537 and 450 subjects, respectively (Table 6). The composite endpoints are presented using the protocol definitions for MI. The component endpoints are based on non-hierarchical counts unless otherwise noted.

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Note: Subjects are only counted once for each type of event in each time period.  
Note: All table entries were calculated based on the data from data cut-off date of December 3, 2012.

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4.4.3 Summary of Clinical Results out to 2 years (Interim Data Cutoff)

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1. Data are % (number of subjects with event/number of subjects)
2. TLR=target lesion revascularization TV=target vessel
3. MACE=major adverse cardiac events (composite endpoint of Cardiac Death, MI,

- ischemia-driven TLR)
4. TLF=target lesion failure (composite endpoint of Cardiac Death, TVMI, ischemia-driven TVR)
5. TVF=target vessel failure (composite endpoint of Cardiac Death, MI, ischemia-driven TLR, ischemia-driven Non-TLR TVR)

**Conclusion:** The interim outcomes from ABSORB EXTEND further demonstrate the 12 and 24 month safety and performance of Absorb BVS in a larger population of subjects where longer lesion can be treated. Patients will be followed for up to 3 years.

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**Rationale and Method:** Abbott Vascular is using its large body of historical clinical data on XIENCE V stents as a control group for comparison against the currently available Absorb BVS clinical data utilizing the propensity score method. Propensity score analysis is a valid statistical method that is well accepted and commonly used in the clinical trial, and in scientific and statistical communities for the comparison of data from non-randomized studies. This approach adjusts for the baseline imbalance between the test and control groups, creates a better match between both groups, and hence, reduces real and potential bias in the comparative analysis. After propensity score adjustment, baseline characteristics between the two groups are more balanced, and therefore the comparisons of the clinical outcomes become more reliable. Key results for a 6-month and a 1-year propensity analysis comparing Absorb BVS to XIENCE V data are presented using the inverse propensity score weighted (IPSW) method. Propensity analysis at later time points will be conducted as additional data from ongoing trials become available. The results of this comparison are being confirmed through ongoing clinical trials.

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**Results:** The propensity score adjustment yielded adjusted sample sizes of 34 (Absorb BVS) and 66 (XIENCE V) from the 38 and 73, respectively, in the unadjusted arms. The adjusted data were more balanced in baseline patient demographics and angiographic characteristics between both arms. At 6 months, their clinical outcomes and angiographic results were comparable. The adjusted 6-month MACE rates were 2.01% (0.7/34) and 5.27% (3.5/66), respectively, for the Absorb BVS and XIENCE V arms (Table 8) with a lower MACE observed with Absorb BVS, a trend that is consistent with the results from the unadjusted groups.

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| Clinical Outcome at 6 Months | Absorb BVS System (N = 34) | XIENCE V (N = 66)           | P-value |
|------------------------------|----------------------------|-----------------------------|---------|
| MACE [95% CI]                | 2.01% (0.7) [0.00%, 6.73%] | 5.27% (3.5) [0.00%, 10.66%] | 0.4398  |

The 6-month in-device late loss in the IPSW-adjusted arms was 0.14  $\pm$  0.18 mm for the Absorb BVS System vs. 0.12  $\pm$  0.26 mm for the XIENCE V ECSS, representing a non-significant difference in outcomes between the Absorb BVS System and the XIENCE V ECSS (P = 0.68). Similarly, differences in the in-device % diameter stenosis, in-segment late loss, and in-segment % diameter stenosis were not significant between the two devices in this propensity score adjusted analysis (Table 9).

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**Results:** After the propensity adjustment, 282 Absorb BVS and 583 XIENCE V patients remained for the comparative analysis. An improved balance was achieved between both cohorts with the most striking changes observed in the lesion characteristics, which were no longer statistically different between Absorb BVS and XIENCE V cohorts. The MI, ID-TLR, MACE, and TLF (Target Lesion Failure defined as cardiac death, target vessel related MI and ID-TLR) rates were comparable between both adjusted cohorts at 37, 193, and 393 days. A summary of key 1-year outcomes is shown in Table 10. After adjustment, the MACE rates remained comparable, with a numerically lower rate obtained for Absorb BVS (4.15%, 11.7/282) versus XIENCE V (5.63%, 32.8/583). The respective scaffold / stent thrombosis rates were not statistically different (0.40% vs. 0.52%, P=0.80).

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| Clinical Outcomes at 1 year                          | Absorb BVS (N=282)   | XIENCE V (N=583)     | P-value |
|------------------------------------------------------|----------------------|----------------------|---------|
| MACE [95% CI]                                        | 4.15% [1.82%, 6.47%] | 5.63% [3.76%, 7.50%] | 0.35    |
| MI [95% CI]                                          | 2.64% [0.77%, 4.51%] | 2.52% [1.25%, 3.80%] | 0.92    |
| ID-TLR [95% CI]                                      | 1.95% [0.33%, 3.56%] | 3.10% [1.70%, 4.51%] | 0.33    |
| TLF [95% CI]                                         | 4.15% [1.82%, 6.47%] | 5.09% [3.30%, 6.87%] | 0.54    |
| Definite/Probable Scaffold/Stent Thrombosis [95% CI] | 0.40% [0.0%, 1.14%]  | 0.52% [0.0%, 1.11%]  | 0.80    |

**Discussion:** The propensity score adjusted method achieved the intent of balancing the baseline characteristics between both XIENCE V and Absorb BVS groups from different

ABSORB and SPIRIT trials. The adjusted results showed that Absorb BVS is comparable to the XIENCE V stent for in-stent late loss at 6 months and for clinical outcomes at 6 months and 1 year after implantation, and further confirms the comparability between Absorb BVS and XIENCE V in the above parameters found in the unadjusted data.

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**Results:** After propensity score matched adjustment, 2-year MACE and TLF rates remained similar between the Absorb BVS cohort and the XIENCE V cohort (6.7% vs. 8.9% for MACE, and 6.2% vs. 8.2% for TLF, respectively). The BVS cohort exhibited a trend of lower ID-TLR and cardiac death rate as compared to the XIENCE V cohort, although the differences were not significant. The 2-year TVF rate was also numerically lower in the BVS group (7.3% vs. 12.3%, p = 0.09). The adjusted definite / probable scaffold / stent thrombosis rates were similar between Absorb BVS and XIENCE V at 2 years (0.6%, 95% CI [0.01%, 3.09%] vs. 1.4%, 95% CI [0.37%, 3.46%], p=0.65). A summary of the key data is shown below in Table 11.

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| Clinical Outcomes at 2 years                   | Absorb BVS (N=178)   | XIENCE V (N=293)      | P-value |
|------------------------------------------------|----------------------|-----------------------|---------|
| Cardiac Death [95% CI]                         | 0.0% [0.00%, 2.05%]  | 1.4% [0.37%, 3.46%]   | 0.30    |
| MI [95% CI]                                    | 4.5% [1.96%, 8.66%]  | 4.4% [2.38%, 7.47%]   | 1.00    |
| ID-TLR [95% CI]                                | 3.4% [1.25%, 7.19%]  | 3.8% [1.89%, 6.62%]   | 1.00    |
| TLF [95% CI]                                   | 6.2% [3.12%, 10.79%] | 8.2% [5.32%, 11.94%]  | 0.47    |
| MACE [95% CI]                                  | 6.7% [3.53%, 11.48%] | 8.9% [5.88%, 12.73%]  | 0.49    |
| TVF [95% CI]                                   | 7.3% [3.95%, 12.17%] | 12.3% [8.76%, 16.60%] | 0.09    |
| Definite/Probable Scaffold Thrombosis [95% CI] | 0.6% [0.01%, 3.09%]  | 1.4% [0.37%, 3.46%]   | 0.65    |

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**Discussion:** The comparability of Absorb BVS to XIENCE V demonstrated with in-scaffold / stent late loss at 6 months (Cohort B only), and clinical outcomes at 6 months, and 1 year were maintained through 2-years post implantation. The collective data confirms the comparability between Absorb BVS and XIENCE V based on propensity analyses.

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In the ABSORB Cohort A clinical investigation, Absorb BVS showed excellent long-term clinical outcomes with low MACE rates out to 5 years with absence of any cardiac death, ID-TLR, Q-wave MI, and scaffold thrombosis. The ABSORB Cohort B clinical investigation confirms these findings out to 3 years, including patency comparable to XIENCE V. Similarly, ABSORB EXTEND clinical results out to 2 years support Cohort A and B outcomes. Absorb BVS thus performs all the functions of a drug-eluting stent while offering future potential benefits resulting from the absence of a permanent metallic implant. The gradual disappearance of the scaffold allows for the possibility of vascular healing and restoration of vessel function as discussed in the Vasomotor Function Results at 2 Years (sections 8.2.4 – 8.2.5). The absence of a permanent metallic implant in the vascular tissue may facilitate any required re-interventions on the target vessel / lesion or side branches either by percutaneous or surgical means, thus enabling a broader range of treatment options after bioresorption of the scaffold. Finally, unlike permanent metal implants, polymeric implants do not cause imaging artifacts during non-invasive CT or MR evaluation. This provides the additional benefit that a polymeric bioresorbable scaffold may be more compatible with the growing usage of non-invasive follow-up imaging than is the case with metallic stents, potentially facilitating patient management and providing economic benefits.

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It is very important that the patient comply with the postprocedural antiplatelet recommendations. Premature discontinuation of prescribed antiplatelet medication could result in a higher risk of thrombosis, MI, or death. Prior to PCI, if a surgical or dental procedure is anticipated that requires early discontinuation or temporary interruption of antiplatelet therapy, the interventionalist and patient should carefully consider whether an everolimus-eluting scaffold and its associated recommended antiplatelet therapy are the appropriate PCI choice. Following PCI, should a surgical or dental procedure be recommended, the risks and benefits of the procedure should be weighed against the possible risk associated with premature or temporary discontinuation of antiplatelet therapy.

Patients who require premature discontinuation of antiplatelet therapy secondary to significant active bleeding should be monitored carefully for cardiac events, and once stabilized, have their antiplatelet therapy restarted as soon as possible per the discretion of their treating physicians.

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10.2 Materials Required

- 6F/0.070"/1.8 mm minimum inner diameter guiding catheter(s) of appropriate shape for the target vessel
- 2 – 3 syringes (10 – 20 cc)
- 1,000 u/500 cc heparinized normal saline (HepNS)
- 0.014 inch (0.36 mm) x 175 cm (minimum length) guide wire
- Rotating hemostatic valve with 0.096 inch (2.44 mm) minimum inner diameter
- Contrast diluted 1:1 with normal saline
- Inflation device
- Three-way stopcock
- Torque device
- Guide wire introducer

10.3 Vessel and Lesion Selection

- Quantitative imaging is recommended for the assessment of target vessel diameter at baseline for appropriate Absorb BVS size selection.
- The target vessel diameter ranges to be treated in the procedure are indicated in Table 12, along with the Absorb BVS diameter to be used.

| Table 12: Target Vessel Diameter Ranges and Absorb BVS Diameter to be Used |                                |
|----------------------------------------------------------------------------|--------------------------------|
| Target Vessel Diameter Distal and Proximal                                 | Absorb BVS Diameter to be Used |
| ≥ 2.0 mm and ≤ 3.0 mm                                                      | 2.5 mm                         |
| ≥ 2.5 mm and ≤ 3.3 mm                                                      | 3.0 mm                         |
| ≥ 3.0 mm and ≤ 3.8 mm                                                      | 3.5 mm                         |

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1. Peel the pouch open from the top corner.
2. Carefully remove the delivery system from its protective tubing for preparation of the delivery system.
3. Do not bend or kink the hypotube during removal.

10.4.2 Dual Layer Sheath Removal

1. While holding the distal catheter shaft with one hand, grasp only the yellow outer sheath with the other hand and gently slide the sheath distally.
2. A longitudinal split on the inner sheath will open up and be visible.
3. The stylet and dual layer sheath are removed simultaneously from the delivery system by continuing to slide the yellow sheath distally until the inner and outer layers of the dual layer sheath as well as the stylet are free from the catheter system. See Section 6.1 Scaffold Handling – Precautions. Do not use the device if the sheath cannot be removed as indicated.
4. Verify that the scaffold does not extend beyond the radiopaque balloon markers and no scaffold struts are lifted. **Do not use if any defects are noted.**

10.4.3 Guide Wire Lumen Flush

1. Flush the guide wire lumen with HepNS until fluid exits the guide wire exit notch.  
**Note:** Avoid manipulation of the scaffold while flushing the guide wire lumen as this may disrupt the placement of the scaffold on the balloon.

10.4.4 Delivery System Preparation

1. Prepare an inflation device / syringe with diluted contrast medium.
2. Attach an inflation device / syringe to stopcock; attach it to the inflation port of the product. Do not bend the product hypotube when connecting to the inflation device / syringe.
3. With the tip down, orient the delivery system vertically.
4. Open the stopcock to delivery system; pull negative for 30 seconds; release to neutral for contrast fill.
5. Close the stopcock to the delivery system; purge the inflation device / syringe of all air.
6. Repeat steps 3 through 5 until all air is expelled. If bubbles persist, do not use the product.
7. If a syringe was used, attach a prepared inflation device to the stopcock.
8. Open the stopcock to the delivery system.
9. Leave on neutral.

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10.5 Delivery Procedure

1. Prepare the vascular access site according to standard practice.
2. **Predilate the lesion with a percutaneous transluminal coronary angioplasty catheter.** Predilatation is strongly recommended and can also be utilized to properly size the vessel.  
**Note:** Limit the length of the predilatation by the PTCA balloon to avoid creating a region of vessel injury that is outside the boundaries of the Absorb scaffold.
3. Administer a standard dose of intracoronary nitroglycerine prior to finalizing the RVD within the target zone.
4. Maintain neutral pressure on the inflation device. Open the rotating hemostatic valve as widely as possible.
5. Backload the delivery system onto the proximal portion of the guide wire, while maintaining guide wire position across the target lesion.
6. Advance the delivery system over the guide wire to the target lesion. Utilize radiopaque balloon markers to position the scaffold across the lesion; perform angiography to confirm scaffold position.  
**Note:** If removal of a scaffold system is required prior to deployment, ensure that the guide catheter is coaxially positioned relative to the scaffold delivery system and cautiously withdraw the scaffold delivery system into the guiding catheter. Should **unusual resistance be felt at any time** when withdrawing the scaffold into the guide catheter, the scaffold delivery system and the guide catheter should be **removed as a single unit**. This should be done under direct visualization with fluoroscopy.
7. Tighten the rotating hemostatic valve. The scaffold is now ready to be deployed.

10.6 Deployment Procedure

1. **CAUTION: Refer to product label for *in vitro* scaffold inner diameter, nominal pressure and RBP.**
2. Prior to deployment, reconfirm the correct position of the scaffold relative to the target lesion using the radiopaque balloon markers.
3. Deploy the scaffold slowly, by pressurizing delivery system in **2-atm increments, every 5 seconds**, until scaffold is completely expanded. Maintain pressure for 30 seconds. Fully expand the scaffold by inflating to nominal pressure at a minimum. Accepted practice generally targets an initial deployment pressure that would achieve a scaffold inner ratio of about 1.1 times the reference vessel diameter.  
**CAUTION: Do not exceed the labeled rated burst pressure RBP of 16 atm (1621 kPa) or maximum deployment diameter of the scaffold.**
4. Fluoroscopic visualization during scaffold expansion should be used in order to properly judge the optimum scaffold diameter as compared to the proximal and distal native coronary artery diameters (reference vessel diameters). Optimal scaffold expansion and proper apposition require that the scaffold be in full contact with the arterial wall.
5. If necessary, the delivery system can be repressurized or further pressurized to ensure complete apposition of the scaffold to the artery wall.  
Fully cover the entire lesion and balloon-treated area (including dissections) with the Absorb scaffold, allowing for adequate scaffold coverage into healthy tissue proximal and distal to the lesion.
6. Deflate the balloon by pulling negative on the inflation device for 30 seconds.  
Confirm complete balloon deflation before attempting to move the delivery system. If unusual resistance is felt during scaffold delivery system withdrawal, pay particular attention to the guiding catheter position.  
**Note:** See section 10.8 Removal Procedure for instruction on withdrawal of scaffold delivery system.
7. Confirm scaffold position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the scaffold. Fluoroscopic visualization during scaffold expansion should be used in order to properly judge the optimum expanded scaffold diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the scaffold be in full contact with the artery wall, which can be facilitated with the use of routine angiography, intravascular ultrasound (IVUS), or optical coherence tomography (OCT).
8. **Postdilatation with a noncompliant balloon is recommended** following instructions in section 10.7, as long as the postdilated segment is within the allowable expansion limits of the scaffold.

10.7 FURTHER EXPANSION OF THE DEPLOYED SCAFFOLD

1. **DEPLOYED SCAFFOLDS SHOULD NOT BE LEFT UNDERDILATED.**  
Deployed scaffolds should be well-apposed to the vessel wall. To achieve optimal scaffold apposition, postdilatation is recommended with a low profile, high pressure, noncompliant balloon, as long as the postdilated segment is within the allowable expansion limits of the scaffold. The scaffolded segment should be carefully recrossed with a prolapsed guide wire to avoid disrupting the scaffold geometry. Postdilatation must only be done with balloons sized to fit within the boundaries of the scaffold. **The compliance chart of the noncompliant balloon selected must be carefully reviewed prior to dilatation and an appropriate maximum pressure used to ensure that the scaffold is not over-dilated.** The scaffold dilatation limit is 0.5 mm above the nominal diameter for all scaffold sizes.  
**CAUTION: Do not dilate the scaffold beyond the dilatation limit which is 0.5 mm above the nominal diameter. Over-dilatation may result in scaffold damage.**
- | Nominal Scaffold Diameter | Dilatation Limit                           |
|---------------------------|--------------------------------------------|
| 2.5 mm                    | 3.00 mm<br>Maximum postdilatation diameter |
| 3.0 mm                    | 3.50 mm<br>Maximum postdilatation diameter |
| 3.5 mm                    | 4.00 mm<br>Maximum postdilatation diameter |
2. If more than one Absorb BVS is needed to cover the lesion and balloon-treated area, it is suggested that, to avoid the potential for gap stenosis, the scaffolds be overlapped by a minimum of 1 mm and a maximum of 4 mm. To ensure that there are no gaps between scaffolds, the balloon marker bands of the second Absorb BVS should be positioned inside the deployed scaffold prior to expansion. It is recommended not to use more than two Absorb BVS to treat one lesion.
  3. Ensure the final scaffold diameter matches the reference vessel diameter to **ENSURE GOOD SCAFFOLD APPPOSITION**. Reconfirm scaffold position and angiographic results. Repeat inflations until achieved.

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1. Deflate the balloon by pulling negative pressure on the inflation device. Larger and longer balloons will take more time (up to 30 seconds) to deflate than smaller and shorter balloons. Confirm balloon deflation under fluoroscopy.
  2. Position the inflation device to “negative” or “neutral” pressure.
  3. Fully open the rotating hemostatic valve.
  4. Stabilize guide catheter position just inside coronary ostium and anchor in place. Maintain guide wire placement across scaffolded segment.
  5. Gently remove the scaffold delivery system with slow and steady pressure.
  6. Tighten the rotating hemostatic valve.
- If, during withdrawal of the catheter from the deployed scaffold, resistance is encountered, use the following steps to improve balloon rewrap:**
- Re-inflate the balloon up to nominal pressure, deflate and change pressure to neutral.
  - Repeat steps 1 – 5 above.
  - Re-evaluate the scaffolded region once the balloon is removed to assure optimal apposition.
- Note:** After successful withdrawal of the balloon from the deployed scaffold, should **any resistance be felt at any time** when withdrawing the scaffold delivery system or post-dilatation balloon into the guide catheter, **remove the entire system as a single unit**. See section 6.4 Scaffold / System Removal – Precautions for specific delivery system removal instructions.

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<sup>5</sup> Claessen BE, Beijk MA, Legrand V, Ruzulyo W, Manari A, Varenne O, Suttrop MJ, Tijssen JG, Miquel-Hebert K, Veldhof S, Henriques JP, Serruys PW, Piek JJ. Two-year clinical, angiographic, and intravascular ultrasound follow-up of the XIENCE V everolimus-eluting stent in the treatment of patients with *de novo* native coronary artery lesions: the SPIRIT II trial. Circ Cardiovasc Interv. 2009 Aug;2(4): 339-47.

<sup>9</sup> Rosenbaum PR. Model-based direct adjustment. J Am Stat Assoc. 1987;82: 387-394

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